

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
Online Spares available at www.invacare.ie

JULY 2018

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® Dragon Standard Seating
Ireland
Max User Weight 140Kg (22 Stone)

Standard Features



- Standard seat unit with 50mm seat cushion and Tension adjustable back
- Seat Widths 440mm - 490mm (17" - 19")
- Seat Depths 410mm - 510mm (16" - 20")
- Detachable T-Armrests with Standard armpad
- Manually settable seat tilt and Manual settable recline (-10° - 30°)
- Standard 80° detachable legrest hangers with adapt footplate
- Shark 2 Remote with Standard Bracket
- Blue Metallic Frame colour
- 230mm puncture proof PU-foam sealed, Grey castors
- 318mm x 57mm (12.5" x 2.25") Puncture proof solid inlay Grey drive wheels
- Rear suspension
- 6 km/h (4mph) 2 Pole Motors and 50 Ah Batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock
- Wheel lock / parking brake
- 8A (200 - 250V) Charger



* Configurations may vary from the image shown

Standard Base Price €3,676

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
Customer weight =		Kg		
Code		Description		Ireland
AML0000		Dragon Basic Model		€3,676
AMZ0000		Dragon Special Basic Model - PLEASE COMPLETE DETAILS AT THE END OF THE FORM		POA
SEATING TYPE				
AML1150	<input checked="" type="checkbox"/>	Standard Seat Unit	◆	
SEAT CUSHION				
AML1336	<input checked="" type="checkbox"/>	Seat Cushion 50mm	◆	
BACK CUSHION				
AML1178	<input checked="" type="checkbox"/>	Tension adjustable Sling back	◆	
UPHOLSTERY TYPE				
AML0268	<input checked="" type="checkbox"/>	Black (Nylon)	◆	
SEAT HEIGHT (Without Seat Cushion)				
AML0362	<input checked="" type="checkbox"/>	450mm (17.75") (adjustable from 425mm - 500mm)	◆	
SEAT WIDTH				
AML0343	<input checked="" type="checkbox"/>	440mm - 490mm (430mm seat cushion) (17" - 19")	◆	
SEAT DEPTH				
AML0326		410mm (16")	○	Only available with Settable Recliner
AML0327		460mm (18")	○	
AML0330		510mm (20")	○	

BACKREST HEIGHT (without Seat Cushion)			
AML0448		480mm (19")	○
AML0454		540mm (21")	○
ARMREST (Pair)			
AML0606	x	Detachable T-Armrest, adjustable in height 290mm - 360mm	◆
ARMPAD (Pair)			
AML0647		Standard Armpad (325mm)	◆
AML0649		Extended Armpad (+50mm)	○
SEAT ACCESSORIES			
HEADREST - Adequate headrest option must be selected			
AML1353		Contour Headrest	€21
AML1427		Headrest with cheek support, sideways adjustable	€57
AML1495		Equivalent non-Invacare Headrest supporting the requirements of ISO7176-19:2001, provided independently from this order	○
SEAT TILT			
AML1002		Manually settable Tilt	◆
AML1014		Electrical Tilt (-4.5° to 25.5°)	€105
BACKREST RECLINER			
AML0958		Settable Recliner (-10° to 30°)	◆
AML0953		Pneumatic Recliner RH (5° to 30°)	Not available with RH control and Tray Bracket €315
AML0954		Pneumatic Recliner LH (5° to 30°)	Not available with LH control and Tray Bracket €315
AML0947		Electrical Recliner (4° - 20°)	€499
FOOTREST HANGER TYPE (Pair)			
AML0526		Standard Hanger 80°	◆
AML0529		Manual Height Adjustable (Elevating Legrest)	€42
FOOTPLATE TYPE (Pair)			
AML0513	x	Angle and depth adjustable Footplate (Adapt)	◆
REMOTE CONTROL			
STANDARD REMOTE SHARK 2 & BRACKET (Please also select bracketry option)			
AML1801	x	Shark 2 Remote (enables operation of up to 2 electrical adjustments)	◆
AML1207		Standard Bracket RH	○
AML1208		Standard Bracket LH	○
AML1707		Swing Away + height adjustable RH high end	○
AML1708		Swing Away + height adjustable LH high end	○
DUAL CONTROL (can only be selected as an addition to user remote)			
AML1891		Dual Control RH (price applies to dual control and bracket)	€78
AML1892		Dual Control LH (price applies to dual control and bracket)	€78
DRIVE UNIT OPTIONS			
FRAME COLOUR			
AML0259	x	Blue Metallic (Rear wheel rim Black and front fork Silver)	◆
FRONT WHEEL			
AML0811		230mm (9") Grey Pneumatic	○
AML0813		230mm (9") Grey Puncture Proof (Solid Inlay)	◆
REAR WHEEL			
AML0711		318mm x 57mm (12.5" x 2.25") Grey Pneumatic	○
AML0727		318mm x 57mm (12.5" x 2.25") Grey Puncture Proof (Solid Inlay)	◆
PERFORMANCE OPTIONS			
SUSPENSION			
AML1264	x	Rear Suspension	◆
MOTOR TYPE			
AML1282	x	6 km/h (4mph) Standard version	◆
BATTERY			
AML1248		50 Ah (AGM)	◆
AML1259		60 Ah (Gel)	€42
PERFORMANCE ACCESSORIES			
CHARGER			
AML1343		OptiCharge Charger 8A (200 - 250V)	◆
AML1313		OptiCharge Charger 10 (200 - 250V) Charging time will be 20% faster	€133
SAFETY OPTIONS			
ANTI TIPPERS			
AML1205	x	Anti tippers	◆
LIGHTS			
AML1254	x	Without Light system	◆

KERB RIDER			
AML1244		Central Kerb Rider with front fixation	€10
TRANSPORTATION KIT			
AML1494	x	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is recommended)	◆
FIXATION BELT			
AML1329	x	Lap Belt with Auto Lock	◆
ACCESSORIES			
LUGGAGE RACK			
AML1394	x	Rear Carrier	◆
CANE HOLDER			
AML1312		Cane Holder RH	€42
AML1363		Cane Holder LH	€42
TRAY - mounted opposite side to controller			
AML1404		Swing Away RH	€21
AML1405		Swing Away LH	€21
Total Price			



Special Instructions - Please supply the following information:

Name:

Email:

Phone:

Fax:

Description of special request (please include relevant option code or part number):