



Yes, you can:

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
Online Spares available at: www.invacare.ie

Powerchairs Ireland Prescription Form

January 2022
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

ESPRIT Action 4 NG Junior

Standard Features

- Max User Weight : 80kg (12.6 Stone)
- Alber E-FIX 12" x 2.1/4" propulsion
- Single crossbar
- Seat widths 305mm to 380mm 12" - 15"
- Seat depth 325mm (13") to 375mm (15") adjustable (set at 350mm)
- Seat height 460mm (18")
- Nylon seat upholstery
- Half folding 7° 430mm backrest
- Standard nylon sling type backrest upholstery
- 80° swing in / swing out legrest hangers
- Flip up angle and depth adjustable footplates
- Height adjustable, removable arm armrests
- Choice of 7 standard colours (with 8 Special colour options)
- Push to lock brakes
- Anti tippers with chassis jack and buckle belt

Standard Base price: €4,787



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
SFE0200	☐ ESPRIT ACTION 4NG JUNIOR	€4,787
MODEL*		
SFE0050	☐ Esprit Action 4NG - Junior Version	◆
CHAIR TYPE*		
SFE0100	☐ Esprit Action 4NG (NC 0200)	◆
SFZ0100	☐ Special Esprit Action 4NG	POA
CHASSIS*		
SFE0200	☐ Single Crossbar	◆

Customer Ref:

Order Number:

PROPULSION TYPE*			
SFE0300	☐	Alber E-FIX 12" x 2,1/4" (E26, 125 kg) Pneumatic tyres + solid insert + bumper	◆
USER REMOTE*			
SFE0400	☐	User remote control RH	◆
SFE0410	☐	User remote control LH	○
SFE0420	☐	User attendant remote (NC 0440, 0450, 0460, 0480, 0630, 5110, 5120)	€60
DUAL REMOTE			
SFE0440	☐	Intuitive Dual remote RH (NC 0480, 2000, 2010)	€916
SFE0450	☐	Intuitive Dual remote LH (NC 0480, 2000, 2010)	€916
SFE0460	☐	Intuitive single Dual remote (NC 2000, 2010)	€654
DUAL REMOTE BRACKET			
SFE0480	☐	Adjustable bracket for intuitive single DUAL remote	○
BATTERY*			
SFE0500	☐	Battery 12Ah	◆
POWER OPTIONS			
SFE0600	☐	Additional Battery 12Ah	€291
SFE0620	☐	Car additional charger	€360
SFE0630	☐	Remote protector	€66
SFE0640	☐	External charging plug	€80
SEAT WIDTH*			
SFE1000	☐	305mm (12")	○
SFE1010	☐	330mm (13")	○
SFE1020	☐	355mm (14")	○
SFE1030	☐	380mm (15")	○
SEAT DEPTH*			
SFE1200	☐	325mm to 375mm (13" - 15")	◆
SEAT HEIGHT*			
SFE1300	☐	460mm (18")	◆
SEAT UPHOLSTERY TYPE*			
SFE1400	☐	Nylon with hook and loop strips, width & depth adjustable	◆
BACKREST TYPE*			
SFE2000	☐	Folding backrest 430mm Height	◆
SFE2010	☐	Folding Backrest 510mm Height (NC 2120)	○
SFE2020	☐	Reclining (ratchet) 0° / 30°, 3 height 390mm / 410mm / 430mm	€133
BACKREST UPHOLSTERY TYPE*			
SFE2100	☐	Width & height adjustable nylon sling type	◆
SFE2120	☐	Contour adjustable (10 + 1 Straps) padded shaped cover	○
HANGERS RH*			
SFE3000	☐	80° short RH swing in / swing out (NC 3100)	○
SFE3010	☐	80° long RH swing in / swing out (NC 3100)	€41
HANGERS LH*			
SFE3100	☐	80° short LH swing in / swing out	◆
SFE3110	☐	80° long RH swing in / swing out	€41
FOOTPLATES RH*			
SFE3210	☐	Flip up angle and depth adjustable RH footplate	◆
SFE3220	☐	Flip up one piece footplate (NC 3310)	€81

FOOTPLATES LH*			
SFE3310	☐	Flip up angle and depth adjustable LH footplate	☐
ARMRESTS RH*			
SFE4000	☐	Removable RH desk	☐
SFE4010	☐	Removable RH full length armpad	☒
ARMRESTS LH*			
SFE4100	☐	Removable LH desk	☐
SFE4110	☐	Removable LH full length armpad	☒
SEATING OPTIONS			
SFE5021	☐	Essential cushion depth 375mm	☐
SFE5110	☐	Headrest (NC 5120)	€165
SFE5120	☐	Headrest with cheek support	€190
SFE5210	☐	Adjustable calf strap	☐
COLOURS*			
SFE6030	☐	Pearl White	☐
SFE6031	☐	Blue Pop	☐
SFE6032	☐	Hypnotic Purple	☐
SFE6033	☐	Red Matt	☐
SFE6034	☐	Matt Black	☐
SFE6035	☐	Polar Blue	☐
SFE6041	☐	Candy Pink	☐
SFE6051	☐	Special RAL Colour	€196
SWO0801	☐	Slate Grey	€196
SWO0806	☐	Electric Green	€196
SWO0820	☐	Anthracite Black	€196
SWO0607	☐	Dark Red	€196
SWO0515	☐	Pearl Grey	€196
SWO0647	☐	Intense Yellow	€196
SWO0648	☐	Vibrant Orange	€196
FRONT WHEELS*			
SFE7000	☐	175mm x 50mm (7" x 2") Pneumatic	☐
SFE7010	☐	175mm x 50mm (7" x 2") Solid	☒
MANUAL REAR WHEELS			
SFE7110	☐	560mm x 35mm (22" x 1.3/8") Pneumatic (NC 8100)	€239
SFE7140	☐	600mm x 35mm (24" x 1.3/8") Solid (NC 8100)	€239
HANDRIM TYPE FOR ADDITIONAL MANUAL REAR WHEEL			
SFE7200	☐	Aluminium pre-anodised handrim	☐
AXLE TYPE FOR ADDITIONAL MANUAL REAR WHEEL			
SFE7300	☐	Quick release axle	☐
SAFETY OPTIONS			
SFE8000	☐	Rear anti-tippers (NC 8010)	☐
SFE8010	☐	Rear anti-tippers with chassis jacks	☒
SFE8020	☐	Buckle belt	☒
SFE8025	☐	Transportation Kit	☒

OTHER OPTIONS			
SFE8100		Cane holder	€21
SFE8200		Acrylic flip up tray	€331

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:

SFE0200	ESPRIT ACTION 4NG JUNIOR
SFE0200	Single Crossbar
SFE0440	Intuitive Dual remote RH
SFE0450	Intuitive Dual remote LH
SFE0460	Intuitive single Dual remote
SFE0480	Adjustable bracket for intuitive single DUAL remote
SFE0630	Remote protector
SFE2000	Folding backrest 430mm Height
SFE2010	Folding Backrest 510mm Height
SFE2120	Contour adjustable (10 + 1 Straps) padded shaped cover
SFE3100	80° short LH swing in / swing out
SFE3310	Flip up angle and depth adjustable LH footplate
SFE5110	Headrest
SFE5120	Headrest with cheek support
SFE8010	Rear anti-tippers with chassis jacks
SFE8100	Cane holder



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Fax:

Description of special request (please include relevant option code or part number):