

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
Online Spares available at www.invacare.ie

JANUARY 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® FOX
Ireland
Max User Weight 127Kg (20 stone)

Standard Features



- Canvas sling seat with Black upholstery and backrest height 490mm
- Seat widths from 350mm - 500mm
- Seat depths from 370mm - 460mm
- Seat to floor heights 480mm or 510mm
- Detachable T-armrest with standard armpads
- Manual settable seat Tilt (Set at 4°, 0° - 8° achievable)
- Settable backrest Recline (85° - 115°)
- Standard 80° Legrest hangers with fixed flip-up footplates
- **LiNX** Control REM110 with standard bracket
- Choice of 5 colour options
- Solid Grey castors 200mm x 50mm (8" x 2"), Silver Rim
- Grey Drive Wheels, Puncture Proof Solid Inlay, (12 1/2" x 2 1/4"), Silver Rim (aluminium)
- 6 km/h (4mph) Motors with 34 Ah AGM Battery
- 8A (200 - 250V) Charger
- Bolted on anti tipper
- Wheel lock / parking brake
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Turning radius 1600mm - 1630mm depending on configuration



* Configurations may vary from the image shown

Standard Price : €3,045

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
Customer weight =		Kg		
Code		Description		Ireland
AXA0000		Fox Basic Model		€3,045
AXZ0000		Fox Special Basic Model		POA
SEATING TYPE				
AXA1140	X	Fox Seat Unit		◆
SEAT FRAME / SEAT WIDTH SETTING				
AXA0368		350mm (14") achieved by inrigged arms Seat frame 400mm	Not available with both RH & LH amputee supports or detachable T-armrests	○
AXA0361		400mm (16") (seat frame 1) + 50mm armrest adjustment	Not available with both RH & LH amputee supports	○
AXA0345		450mm (18") (seat frame 2) + 50mm armrest adjustment		○
AXA0350		500mm (20") (seat frame 3) + 50mm armrest adjustment only with seat board	Not available with sling seat or sling upholstery	○
SEAT DEPTH SETTING				
AXA0324		370mm to the front of the back posts Wheelbase setting: short / seat mounting position (COG): Rear		○
AXA0308		400mm to the front of the back posts Wheelbase setting: long / seat mounting position (COG): Front		○
AXA0327		430mm to the front of the back posts Wheelbase setting: long / seat mounting position (COG): Front		○
AXA0319		460mm to the front of the back posts Wheelbase setting: long / seat mounting position (COG): Front		○

SEAT TYPE				
AXA1187		Sling Seat	Not available with matrix Libra cushion	◆
AXA1188		Seat Board		○
SEAT CUSHION				
AXA1409		Without Seat Cushion		◆
AXA1336		Seat Cushion 50mm thick (Black Nylon)		€45
AXA2147		Matrix Libra Cushion, 90mm thick		€483
BACKREST UPHOLSTERY / CUSHION				
AXA1184		Sling Upholstery		◆
AXA1178		Tension Adjustable Backrest (Black Nylon Upholstery)		€53
AXA2174		Matrix Posture Back, height 400mm (Mesh fabric)	Not available with 500mm seat width	€519
AXA2175		Matrix Posture Back, height 500mm (Mesh fabric)	Not available with 500mm seat width	€519
AXA2176		Matrix Elite Back, height 500mm (Mesh fabric)	Not available with 500mm seat width	€519
AXA2177		Matrix Elite TR Back, height 500mm (Mesh fabric)		€519
AXA2178		Matrix Elite Deep Back, height 500mm (Mesh fabric)		€672
BACKREST HEIGHT for SLING UPHOLSTERY + TENSION ADJUSTABLE				
AXA0449		490mm		◆
AXA0457		570mm		€53
BACKREST SUPPORT FOR MATRIX BACKS				
AXA0413		Backrest frame short (470mm)		○
AXA0414		Backrest frame long (540mm)		○
ARMREST (Pair)				
AXA0606		T-armrest with depth (50mm) and width adjustable receiver frame width + 25mm each side	Not available with short armpad	◆
AXA0677		Inrigged Armrests (25mm each side)	Not available with long ergonomic armpad	€76
ARMPADS (Pair)				
AXA0650		Long Ergonomic Armpad (397mm)		○
AXA0647		Standard Armpad (325mm)		◆
AXA0649		Long Armpad (375mm)		○
AXA0645		Short Armpad (220mm)		○
SEAT ACCESSORIES				
HEADREST / NECKREST (Dartex)				
AXA1350		Neckrest, depth adjustable		€170
AXA1432		Headrest, sideways adjustable		€218
SEAT ADJUSTMENTS				
SEAT HEIGHT SETTING				
AXA0312		Low Seat To Floor height 480mm at 4° factory setting (-30mm / + 30mm depending on Tilt settings)		○
AXA0344		High Seat To Floor height 510mm at 4° factory setting (- 30mm / + 30mm depending on Tilt settings)		○
SEAT TILT				
AXA1002	X	Manually settable Seat Tilt (0° / 2° / 4° / 6° / 8°), Factory setting: 4°		◆
BACKREST RECLINER				
AXA0961	X	Recliner, settable in steps (85° / 95° / 105° / 115°)		◆
LEGREST / FOOTREST				
FOOTREST HANGERS (Pair)				
AXA0526		Standard 80° hangers		◆
AXA0529		Manual Elevating Legrests		€205
AMPUTEE SUPPORTS (Replaces legrest hanger)				
AXA0554		Amputee Support RH		€192
AXA0555		Amputee support LH		€192
FOOTPLATE TYPE (Pair)				
AXA0512		Fixed Angle, Flip-up Footplates		◆
AXA0513		Angle and Depth Adjustable Footplates (Adapt)		€87




LINX PRIMARY CONTROLS

Drive Only



REM110

Dual Control



REM050

Compact Dual Control



ACU200

LINX PRIMARY REMOTE & BRACKET			
AXA2201		REM110 Drive only	◆
AXA1207		Standard Bracket RH	○
AXA1208		Standard Bracket LH	○
AXA1707		Swing Away RH	€81
AXA1708		Swing Away LH	€81
AXA2227		Bracket RH mounted for Attendant	€65
AXA2228		Bracket LH mounted for Attendant	€65
LINX DUAL CONTROL & BRACKET: CAN ONLY BE SELECTED IN ADDITION TO LINX USER CONTROL			
AXA2220		Dual Control (REM050 as secondary remote)	€445
AXA2219		Compact Dual Control (ACU200)	€424
AXA1207		Standard Bracket RH	○
AXA1208		Standard Bracket LH	○
<div style="display: flex; justify-content: space-around; align-items: flex-start;"> <div style="text-align: center;"> <p>NEW 3" U-Shape</p> </div> <div style="text-align: center;"> <p>NEW 4" U-Shape</p> </div> <div style="text-align: center;"> <p>NEW Rubber dome</p> </div> <div style="text-align: center;"> <p>NEW Grooved Mushroom</p> </div> </div>			
BODYPOINT JOYSTICK KNOBS NEW!			
AXA2721		3" wide U-Shaped, Joystick Handle	€68
AXA2722		4" wide U-Shaped, Joystick Handle	€68
AXA2723		Wide Rubber Dome, Joystick Handle	€68
AXA2724		Grooved Mushroom Joystick Handle	€68
WHEELBASE SETTING			
AXA1226	X	Recommended factory setting, based on seat depth and seat to floor height <i>Wheel setting positions can be changed as desired after performing the appropriate risk assessment</i>	◆
SHROUD COLOUR			
AXA0201		Pearl Night Blue	○
AXA0233		Metallic Silver	○
AXA0228		Silky Matt Black	○
AXA0204		Transparent Cherry Red	○
AXA0222		White	○
CASTORS			
AXA0802		Castor Grey 200mm x 50mm (8" x 2") PU Tyre, Black Rim	○
AXA0803		Castor Grey 200mm x 50mm (8" x 2") Solid Inlay, Silver Rim	◆
DRIVE WHEELS			
AXA0709		Drive Wheel Grey, Puncture Proof Solid Inlay, (12 1/2" x 2 1/4"), Black Rim (plastic)	○
AXA0710		Drive Wheel Grey, Puncture Proof Solid Inlay, (12 1/2" x 2 1/4"), Silver Rim (aluminium)	◆
PERFORMANCE OPTIONS			
MOTOR TYPE			
AXA1294		3 km/h (2mph)	○
AXA1282		6 km/h (4mph)	◆
BATTERY			
AXA1486		34 Ah AGM	◆
AXA1255		40 Ah Gel	€117
PERFORMANCE ACCESSORIES			
CHARGER			
AXA1343		OptiCharge Charger 8A (200 - 250V)	◆
AXA1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€ 158
SAFETY ACCESSORIES			
WHEEL LOCK / PARKING BRAKE			
AXA1206	X	Manual Wheel Locks	◆
KERB RIDER			
AXA1244		Kerb Rider	€205
TRANSPORTATION KIT			
AXA1494	X	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is recommended)	◆
LAP BELT			
AXA1329	X	Lap Belt with Auto-Lock	◆
ACCESSORIES			
CANE HOLDER			
AXA1312		Cane holder RH	€25
AXA1363		Cane holder LH	€25

TRAY - mounted opposite side to controller		
AXA1404	Tray, Swing Away RH	€213
AXA1405	Tray, Swing Away LH	€213
OXYGEN BOTTLE HOLDER		
AXA1454	Oxygen bottle holder	€104
CASTOR LOCKS		
AXA1490	Castor Lock Levers	€32
		Total Price:

Seat height rear	Tilt (°)	Seat height front
45	0	45
45	2	47
45	4	48
45	6	50
45	8	51
48	0	48
48	2	50
48	4	51
48	6	53
48	8	54

* Seat Depth Measurement points according to ISO 7176-7:1998 are from the middle of the backrest upholstery to the front edge of the seat.

AXA0324	400 mm according to ISO 7176-7:1998
AXA0308	430 mm according to ISO 7176-7:1998
AXA0327	460 mm according to ISO 7176-7:1998
AXA0319	490 mm according to ISO 7176-7:1998

Note: Invacare Fox User and Service Manuals refer to seat depths lengths as listed in the above table



Special Instructions - Please supply the following information:

Name:

Email:

Phone:

Fax:

Description of special request (please include relevant option code or part number):