



Yes, you can:

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
www.invacare.ie

AUGUST 2022
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Küschall K-series Attract 2.0

Standard Features

- Max User Weight 130kg (20.5 stone)
- Crash Tested ISO 7176-Part 19
- Active Rigid Wheelchair with Height Adjustable Rigid Footrest
- Multiple adjustments, Front and Rear Seat Heights and Centre of Gravity
- Full Hydroformed Aluminium Front Frame 75° Front Angle
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims and Pneumatic Tyres
- Choice of 7 Seat Widths 360mm to 480mm
- Choice of 4 Seat Depths 375mm to 450mm
- Range of Front Seat Heights from 480mm to 500mm
- Range of Rear Seat Heights from 400mm to 480mm
- Lockdown Folding, Angle and Height Adjustable Backrest
- Black Nylon, Tension Adjustable Upholstery
- Without Push Handles, Composite Clothesguards
- Black Aluminium Axle Tube with 3° Camber
- Black Comfort Front Castor
- Standard Brake Lever Bent

Standard Base price: €1,151



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
DKB0070	<input type="checkbox"/> KÜSCHALL K-SERIES ATTRACT 2.0	€1,151
FRAME*		
DKB1050	<input type="checkbox"/> Tapered frame angle 75° to the ground aluminium	◆
DKB1080	<input type="checkbox"/> Tapered frame angle 90° to the ground aluminium (NC 5300, 5350)	○

SEAT WIDTH*			
DKB1320		360 mm	○
DKB1325		380 mm	○
DKB1330		400 mm	◆
DKB1335		420 mm	○
DKB1340		440 mm	○
DKB1345		460 mm	○
DKB1350		480 mm	○
SEAT DEPTH*			
DKB1420		375 mm	○
DKB1435		400 mm	◆
DKB1450		425 mm	○
DKB1465		450 mm	○
SEAT POSITION. Seat module can be moved forward and backward along the frame to achieve a more active or stable rear wheel position*			
DKB1500		Position 1 = 140 mm MOST ACTIVE POSITION	○
DKB1510		Position 2 = 125 mm	○
DKB1520		Position 3 = 110 mm	◆
DKB1530		Position 4 = 95 mm	○
BACK HEIGHT. Transportation kit is only possible with backrest heights of 405mm and above*			
DKB1620		315 mm (NC 7600)	○
DKB1625		330 mm (NC 7600)	○
DKB1630		345 mm (NC 7600)	○
DKB1635		360 mm (NC 7600)	○
DKB1640		375 mm (NC 7600)	○
DKB1645		390 mm (NC 7600)	○
DKB1650		405 mm	○
DKB1655		420 mm	○
DKB1660		435 mm	○
DKB1665		450 mm	○
DKB1670		465 mm	○
DKB1675		480 mm	○
FRONT SEAT HEIGHT. For 75° frames, the maximum lower leg must be front seat height minus 30mm and minus 40mm for 90° frames*			
DKB1755		480 mm	○
DKB1760		490 mm	○
DKB1765		500 mm	◆
REAR SEAT HEIGHT*			
DKB1820		400 mm (NC 3490, 3495, 6040)	○
DKB1825		410 mm (NC 3490, 3495)	○
DKB1830		420 mm	○
DKB1835		430 mm	○
DKB1840		440 mm	◆
DKB1845		450 mm	○
DKB1850		460 mm	○
DKB1855		470 mm	○
DKB1860		480 mm	○

LOWER LEG LENGTH Cushion height must be deducted*			
DKB1918		290 mm	○
DKB1920		300 mm	○
DKB1922		310 mm	○
DKB1924		320 mm	○
DKB1926		330 mm	○
DKB1928		340 mm	○
DKB1930		350 mm	○
DKB1932		360 mm	○
DKB1934		370 mm	○
DKB1936		380 mm	○
DKB1938		390 mm	○
DKB1940		400 mm	○
DKB1942		410 mm	○
DKB1944		420 mm	○
DKB1946		430 mm	◆
DKB1948		440 mm	○
DKB1950		450 mm	○
DKB1952		460 mm	○
DKB1954		470 mm	○
BACKREST MATERIAL*			
DKB2000		Aluminium	◆
BACKREST TYPE*			
DKB2100		Standard Backrest with Backbrace, lockdown folding for easy transportation	◆
BACKREST UPHOLSTERY*			
DKB2220		Tension Adjustable Backrest, Black Nylon Cover, back pocket/hook and loop fastener	◆
PUSH HANDLES*			
DKB2400		Without Push Handles	◆
DKB2410		Push Handles Fixed Short/Mini Hand Grip	€70
DKB2420		Push Handles Fixed Long Hand Grip	€70
DKB2430		Push Handles Height Adjustable, Rear Set	€150
DKB2450		Push Handles Foldable	€127
BACKREST ANGLE*			
DKB2520		78° angle from the seat	○
DKB2540		82° angle from the seat	○
DKB2560		86° angle from the seat	○
DKB2580		90° angle from the seat	◆
DKB2590		94° angle from the seat	○
SEAT UPHOLSTERY*			
DKB2600		Standard Nylon Seat Upholstery with one 5cm central band (Loop)	◆
SEAT CUSHION			
DKB2700		Standard Seat Cushion, 60 mm Double Foam with Visco Top Layer	€56

FOOTREST*			
DKB3200	<input type="checkbox"/>	Footrest Tube, Rigid - Aluminium Hydroformed Profile	◆
DKB3220	<input type="checkbox"/>	Carbon Footplate, Angle Adjustable (NC 3290)	€122
FOOTREST COVER			
DKB3290	<input type="checkbox"/>	Footrest Cover	€35
CALF STRAP*			
DKB3300	<input type="checkbox"/>	Calf strap mounted on the frame with hook and loop	◆
CLOTHESGUARDS/MUDGUARDS*			
DKB3405	<input type="checkbox"/>	Composite Clothesguards, Flip up, Fold in	◆
DKB3450	<input type="checkbox"/>	Composite Mudguards, Removable (NC 3490, 3495)	€91
ARMREST			
DKB3490	<input type="checkbox"/>	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included	€88
DKB3495	<input type="checkbox"/>	Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, no Clothesguard included	€88
FRAME COLOUR*			
DKB4000	<input type="checkbox"/>	Space Black (RAL 9005)	◆
FRAME DECAL TYPE*			
DKB4195	<input type="checkbox"/>	Tag (Small Kuschall brand transfer on front of frame) (NC 2222)	◆
DCB4197	<input type="checkbox"/>	Without Tag (Kuschall branding only) (NC 4300, 4310)	○
TAG COLOUR*			
DKB4300	<input type="checkbox"/>	Black	○
DKB4310	<input type="checkbox"/>	Red	○
DKB2222	<input type="checkbox"/>	No tag colour, select only with Option DCB4197 Without Tag (NC 4195)	○
FRONT CASTORS*			
DKB5050	<input type="checkbox"/>	100mm x 30mm (4" x 1.2") Everyday Black hub and tyre	○
DKB5100	<input type="checkbox"/>	100mm x 30mm (4" x 1.2") Comfort Black hub and tyre	○
DKB5300	<input type="checkbox"/>	125mm x 30mm (5" x 1.2") Everyday Black hub and tyre	○
DKB5350	<input type="checkbox"/>	125mm x 30mm (5" x 1.2") Comfort Black hub and tyre	◆
REAR WHEEL POSITION*			
DKB6020	<input type="checkbox"/>	Standard, Fixed Rear Wheel Position (NC 6100, 6110, 6120, 6130, 6140)	◆
DKB6040	<input type="checkbox"/>	Rear Wheel Extension (5 Positions) (NC 1111)	€166
CAMBER AXLE TYPE*			
DKB6060	<input type="checkbox"/>	Standard Aluminium Axle	◆
POSITION OF AXLE*			
DKB6100	<input type="checkbox"/>	1st Position from the front	○
DKB6110	<input type="checkbox"/>	2nd Position from the front	○
DKB6120	<input type="checkbox"/>	3rd Position from the front	○
DKB6130	<input type="checkbox"/>	4th Position from the front	○
DKB6140	<input type="checkbox"/>	5th Position from the front	○
DKB1111	<input type="checkbox"/>	Standard fixed position select with DKB6020	○
CAMBER ANGLE*			
DKB6200	<input type="checkbox"/>	0° Camber	○
DKB6210	<input type="checkbox"/>	3° Camber	◆
REAR WHEEL SIZE*			
DKB6410	<input type="checkbox"/>	24"	◆

REAR WHEEL TYPE*			
DKB6500	<input type="checkbox"/>	24" Active Wheel (Includes Light Aluminium handrims, light pneumatic grey tyres as standard) Alternative tyre selection available below	(NC 6700) ◆
DKB6510	<input type="checkbox"/>	24" Kuschall Starec Wheel, Black Cross Spokes, Rim & Hub	€45
REAR WHEEL TYRE*			
DKB6600	<input type="checkbox"/>	24" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection	◆
DKB6620	<input type="checkbox"/>	24" Primo Sentinel tyre. Pneumatic and puncture resistant	€56
DKB6660	<input type="checkbox"/>	24" Solid Profile Black Tyre	€78
HANDRIM*			
DKB6700	<input type="checkbox"/>	24" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance	€78
DKB6780	<input type="checkbox"/>	24" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance	◆
BRAKES*			
DKB7210	<input type="checkbox"/>	Standard Push to Lock Brakes, Lever Bent	◆
DKB7240	<input type="checkbox"/>	Active/Scissor Brakes	(NC 7300) €89
BRAKE LEVER EXTENSION			
DKB7300	<input type="checkbox"/>	Brake Lever Extension	€25
SAFETY OPTIONS			
DKB7500	<input type="checkbox"/>	Anti-Tipper Right	(NC 7520) €83
DKB7510	<input type="checkbox"/>	Anti-Tipper Left	(NC 7520) €83
DKB7520	<input type="checkbox"/>	Anti-Tipper Both sides	€166
DKB7600	<input type="checkbox"/>	Transportation Kit Includes Posture Belt and Labels	(NC 7800) ◆
DKB7800	<input type="checkbox"/>	Standard Pelvic belt	○
OTHER OPTIONS			
DKB8010	<input type="checkbox"/>	Front Frame Protectors; 2 Piece. Removable with hook and loop	€16
DKB8200	<input type="checkbox"/>	Cane Holder, Mounted Right Side	€39
DKB8210	<input type="checkbox"/>	Cane Holder, Mounted Left Side	€39

TOTAL PRICE	
--------------------	--

Non Conforming Options:

DKB3290	Footrest Cover
DKB3490	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included
DKB3495	Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, no Clothesguard included
DKB4195	Tag (Small Kuschall brand transfer on front of frame)
DKB4300	Black
DKB4310	Red
DKB2222	No tag colour, select only with Option DCB4197 Without Tag
DKB5300	125mm x 30mm (5" x 1.2") Everyday Black hub and tyre
DKB5350	125mm x 30mm (5" x 1.2") Comfort Black hub and tyre
DKB6040	Rear Wheel Extension (5 Positions)
DKB6100	1st Position from the front
DKB6110	2nd Position from the front
DKB6120	3rd Position from the front
DKB6130	4th Position from the front
DKB6140	5th Position from the front
DKB1111	Standard fixed position select with DKB6020
DKB6700	24" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance
DKB7300	Brake Lever Extension
DKB7520	Anti-Tipper Both sides
DKB7600	Transportation Kit Includes Posture Belt and Labels
DKB7800	Standard Pelvic belt



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email ordersireland@invacare.com or Fax: 01 810 7085

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):