

Quick guide  
to configurations



## Aquatec<sup>®</sup> Ocean VIP Configuration Guide

Invacare's modular, tilt-in-space  
shower chair commode

**AQUATEC**

[www.safepatienthandling.co.uk](http://www.safepatienthandling.co.uk)



Safe  
Patient  
Handling

# Introducing the *Aquatec® Ocean VIP*

This guide details popular configurations of the *Aquatec Ocean VIP* and shows how accessory selection can support a range of conditions and mild postural needs.

The *Aquatec Ocean VIP* is a modular, tilt-in-space shower chair commode that can be

easily adapted to suit a wide range of clients and conditions.

The seat tilt can be individually adjusted between 0 - 35° to help support posture and aid positioning during hygiene activities.

#### Included as standard are:

- Complete tilt-in-space shower chair
- Soft seat and insert
- Toilet pan, lid and holder
- Leg rests and footplate
- Headrest

Safe Working Load =  
150kg

Product Code: 1517060



# Aquatec® **Ocean VIP** with:

## Ergonomic Soft Seat

Application of the ergonomic soft seat helps to provide pressure relieving properties and promote immersion for pelvic stability.

### Physical needs:

- Increased risk of pressure ulcers
- Mild postural support
- Increased comfort whilst seated

### Ergonomic Soft Seat



Product Code: 1531259

Configuration product code: 1638498



# Aquatec® **Ocean VIP** with:

Ergonomic Soft Seat | Soft Backrest Cushion | Hip Belt

Introduction of the soft backrest cushion provides additional stability and positioning support to individuals with limited trunk control.

The multi-adjustable hip safety belt ensures safety, stability and reassurance during the positioning process.

## Physical needs:

- Trunk support
- Pelvic stability
- Re-assurance and sense of safety



### Soft Backrest Cushion



Product code: 1558769

### Hip Belt



Product code: 1470082

Configuration product code: 1638499

# Aquatec® **Ocean VIP** with:

Ergonomic Soft Seat | Hip Belt | Lateral Supports

Utilisation of the lateral supports helps to reduce the backrest width and enhance the position of an individual's trunk.

The lateral supports can be used either individually, or as a pair, to help provide additional assistance.

## Physical needs:

- Trunk support
- Scoliosis
- Bi-lateral weakness



## Lateral Supports



Product code:  
1535077 (Each)

Configuration product  
code: 1638500

## Aquatec® **Ocean VIP** with:

Ergonomic Soft Seat | Hip Belt | Lateral Supports |  
Multi-adjustable Head Support Stem | Ergonomic Headrest Pad

Introduction of the ergonomic headrest support offers exceptional management of an individual's head position.

The height and angle adjustment positions allow for safe and effective management of an individual's head position.

### Physical needs:

- Pelvic obliquity
- Scoliosis
- Trunk weakness



### Ergonomic Headrest pad and Stem



Product code:  
Stem = 1531613  
Headrest pad = 1531465

Configuration product  
code: 1638501

## Aquatec® **Ocean VIP** Accessories

There is a wide range of additional accessories available that enhance the modularity of the *Aquatec® Ocean VIP*:

**Calf and amputation support (1556894 & 1556893)**



Padded. Both height and angle adjustable.

**Short Foot Rests (1591003 & 1591006)**



Can be used to support abnormal posture or smaller/young adults (12cm shorter than standard foot rests).

**Ocean Skirt (1627829)**



Provides a protective barrier between the commode aperture and automatic toilet systems.

**XL Armrest Kit (1470086)**



Increase the distance between the armrests by 40mm on each side. These can be retrofitted.

**Variable Soft Seat (1544469)**



A 4-way position seat allows hygiene access from front, back and side.

**Universal Soft Seat (1526032)**



Closed seat surface. Smooth and comfortable.

**Visit [www.safepatienthandling.co.uk](http://www.safepatienthandling.co.uk)  
to view the entire range.**



# Meet the *Aquatec® Ocean* Family:

## *Aquatec® Ocean* Shower Chair Commode



A flexible, modular shower chair commode designed to meet a range of clients and conditions.



## *Aquatec® Ocean VIP* and *Dual VIP Shower Chair* Commodes



Modular tilt-in-space models that can compensate for a range of mild postural presentations and physical conditions.



## *Aquatec® Ocean E-VIP* Shower Chair Commode



Fully powered height and tilt adjustment to protect carers and support clients' hygiene needs.



## Scan the QR codes to watch the YouTube videos

Invacare Limited  
Pencoed Technology Park - Pencoed -  
Bridgend - CF35 5AQ -  
United Kingdom  
Tel: +44 1656 776 222  
Fax: +44 1656 776 220  
E-mail:  
safepatienthandling@invacare.com  
Sales Order E-mail:  
ordersuk@invacare.com  
www.invacare.co.uk

© 2018 Invacare International Sàrl. All rights reserved.

All information quoted is believed to be correct at time of print.  
Invacare reserves the right to alter product specifications without prior consultation.

Ocean VIP - UK - 03/2018

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS - County  
Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: ireland@invacare.com  
www.invacare.ie



Invacare UK



@InvacareSPH



Invacare Ltd - UK



Safe  
Patient  
Handling