



Yes, you can:

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
www.invacare.ie

AUGUST 2022
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Azalea Base

Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- Transit Propulsion, Attendant hub brakes and Solid black tyres
- 4 Seat Widths 340mm to 490mm (+50mm, +100mm with Transit)
- 2 Seat Heights 400mm or 450mm
- 2 seat depths 380-450mm (with 340mm width) or 430-500mm
- Seat plate or Trapeze Interface (without Interface option)
- Manually settable backrest recline
- Assistant operated gas piston -1° to +23° seat tilt
- Tilt lock - Pre-settable angle
- Fixed castor attachment
- Pearl Silver frame
- User manoeuvred tyre brakes
- Anti Tippers, Front tiedown brackets

Standard Base price: €1,520



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
BQA0010	Rea Azalea Base	€1,520
BQA8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
PROPULSION*		
BQA0020	Transit 16" wheels with hub brakes (NC 0705, 0715, 0724, 0810, 0850, 0863, 0867)	◆
BQA0030	Self Propel 24"/22" wheels with user brakes (NC 0100, 0430, 0480, 0702, 0726)	€36
BQA0040	Self Propel 24"/22" wheels with assistant hub brakes (NC 0100, 0430, 0480, 0724)	€167

SEAT WIDTH*				
BQA0050		340 mm/13" (Extra small chassis)	(NC 0120, 0164, 0168, 0169, 0180, 0186, 0210, 0260, 0600, 0665, 0675, 0863, 0867, 0900, 0902)	○
BQA0060		390 mm/15" (Small chassis)	(NC 0110)	○
BQA0070		440 mm/17" (Medium chassis)	(NC 0110)	○
BQA0080		490 mm/19" (Large chassis)	(NC 0110)	○
SEAT WIDTH EXTENSION*				
BQA0090		+ 50 mm seat width pre-adjustment on armrests and legrests		○
BQA0100		+ 100 mm seat width pre-adjustment on armrests and legrests	(NC 0863, 0867)	○
BQA0105		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required		○
SEAT DEPTH*				
BQA0110		380-450 mm (15" - 18") only possible with seat width 340mm		○
BQA0120		430-500 mm (17" - 20")		○
SEAT HEIGHT*				
BQA0130		400 mm / 16"	(NC 0180, 0330, 0340, 0390, 0665, 0675, 0715)	○
BQA0140		450 mm / 18"	(NC 0705)	◆
SEAT ANGLE*				
BQA0000		Standard Seat Angle +1° to +25°		◆
BQA0145		Optional Seat Angle -1 to +23°)		○
SEAT INTERFACE*				
BQA0150		Seat Plate standard		◆
BQA0155		Seat Plate with wedges (Seat Angle +/- 5°)		€254
BQA0157		H-Bracket		€91
BQA0170		Without interface (allows mounting of trapeze interface)		○
BACKREST REGULATION*				
BQA0160		Assistant operated gas piston Backrest Recline		€102
BQA0161		Backrest regulation with rigging screw		◆
TILT LOCK*				
BQA0932		Tilt Lock (locking of pre-set angle)		◆
BACKREST*				
BQA0164		Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box)		€777
BQA0168		Matrx Elite (includes hardware mounted to frame, backrest supplied in box)		€599
BQA0169		Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)		€599
BQA0167		Without Backrest	(NC 0172, 0173, 0567)	○
MATRX BACKREST HEIGHT				
BQA0172		Matrx Backrest Height 410 mm / 16"		○
BQA0173		Matrx Backrest Height 510 mm / 20"		○
PUSH HANDLE*				
BQA0174		Push Handles connected		◆
CENTRAL LEGREST				
BQA0180		Central Legrest, includes calfpad and footrest. Also select options BQA0230, BQA0280, BQA1111 and BQA2222	(NC 0186, 0188, 0189, 0190, 0200, 0210, 0220, 0240, 0250, 0260, 0270, 0290, 0295, 0300, 0310, 0330, 0340, 0360, 0370, 0390, 0400, 0580, 0590, 0960)	€253

LEGREST HANGER ATTACHMENT*			
BQA0185		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function	€116
BQA0186		High position 50mm higher for leg length compensation	○
BQA0188		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest (NC 0200, 0210, 0220, 0250, 0260, 0270)	○
LEGREST HANGER LH*			
BQA0190		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH (NC 0300, 0310, 0595, 3333)	€113
BQA0200		Legrest Hanger fixed 80 deg, swing in/swing out LH (NC 0576, 0577, 0580)	€68
BQA0210		Legrest Hanger fixed 90 deg, swing in/swing out LH (NC 0300, 0330, 0340, 0576, 0577, 0580, 0600, 0665, 0675)	€68
BQA0220		Legrest Amputee LH (NC 0290, 0295, 0300, 0310, 0330, 0340, 0580, 0576, 0577, 0595)	€163
BQA0230		Without Legrest, LH (NC 0290, 0295, 0300, 0310, 0330, 0340, 0576, 0577, 0580, 0595)	○
LEGREST HANGER RH*			
BQA0240		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH (NC 0360, 0370, 0595, 3333)	€113
BQA0250		Legrest Hanger fixed 80 deg, swing in/swing out RH (NC 0586, 0587, 0590)	€68
BQA0260		Legrest Hanger fixed 90 deg, swing in/swing out RH (NC 0360, 0390, 0400, 0586, 0587, 0590, 0600, 0665, 0675)	€68
BQA0270		Legrest Amputee RH (NC 0290, 0295, 0360, 0370, 0390, 0400, 0586, 0587, 0590, 0595)	€163
BQA0280		Without Legrest, RH (NC 0290, 0295, 0360, 0370, 0390, 0400, 0586, 0587, 0590, 0595)	○
FOOTREST LH*			
BQA0290		Footrest one piece, depth and angle adjustable (NC 0295, 0300, 0310, 0330, 0340, 0360, 0370, 0390, 0400, 0580, 0590)	€237
BQA0295		Footbox, angle adjustable, padded with short tube (NC 0300, 0310, 0330, 0340, 0360, 0370, 0390, 0400, 0580, 0590, 0595)	€653
BQA0300		Footrest lower part, fixed with short tube and heel strap, LH (NC 0580)	€40
BQA0310		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH	€40
BQA0330		Footrest lower part, fixed with long tube and heel strap, LH. (NC 0580)	€40
BQA0340		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH	€40
BQA1111		Without footrest. Only select with options BQA0180, Central legrest, BQA0220 Amputee legrest and BQA0230 Without legrest	○
FOOTREST RH*			
BQA0360		Footrest lower part, fixed with short tube and heel strap, RH (NC 0590)	€40
BQA0370		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH	€40
BQA0390		Footrest lower part, fixed with long tube and heel strap, RH (NC 0590)	€40
BQA0400		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH	€40
BQA2222		Without footrest. Only select with options BQA0180, Central legrest, BQA0270 Amputee legrest, BQA0280 Without legrest, BQA0290 One piece footrest and BQA0295 Footbox	○

ARMREST LH*			
BQA0420		Armrest Ducktail Pad LH, height adjustable with autolock	€98
BQA0430		Armrest Wide Pad LH, height adjustable with autolock Compatible with 16" wheels only	€98
BQA0435		Armrests Reinforced Wide Pad, LH, height adjustable with autolock	(NC 0900, 0902, 0930) €106
BQA0440		Armrest Padded LH, height adjustable with autolock	(NC 0520, 0530, 0900, 0902, 0920, 0930) €75
BQA0450		Armrest Hemiplegic LH, height adjustable with autolock	(NC 0900, 0902, 0910, 0915, 0920, 0930) €211
BQA0460		Without Armrest LH	(NC 0520, 0530, 0900, 0902, 0910, 0915, 0920, 0930) ○
ARMREST RH*			
BQA0470		Armrest Ducktail Pad RH, height adjustable with autolock	€98
BQA0480		Armrest Wide Pad RH, height adjustable with autolock compatible with 16" wheels only	€98
BQA0485		Armrest Reinforced Wide Pad, RH, height adjustable with autolock	(NC 0900, 0902, 0930) €106
BQA0490		Armrest Padded, RH, height adjustable with autolock.	(NC 0520, 0530, 0900, 0902, 0920, 0930) €75
BQA0500		Armrest Hemiplegic, RH, height adjustable with autolock	(NC 0900, 0902, 0910, 0915, 0920, 0930) €211
BQA0510		Without armrest RH	(NC 0520, 0530, 0900, 0902, 0910, 0915, 0920, 0930) ○
ARMREST OPTIONS. Side support cover 0520, must be selected if 0530 Foam for side support is required			
BQA0520		Side Support Cover, removable with tools (wipeable cover for the sideguard)	€60
BQA0530		Foam for side support cover (internal cushioned pad, which is only possible in combination with BNA0680)	○
BACKREST OPTIONS			
BQA0540		Headrest, Sideway Adjustable	(NC 0560, 0567) €130
BQA0550		Headrest with cheek support, sideway adjustable	(NC 0545, 0560, 0567) €162
BQA0545		Cover for headrest	€37
BQA0560		Neckrest	(NC 0567) €130
BQA0567		Headrest Matrix, 250 mm Soft Pad	(NC 0570) €362
BQA0570		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	€48
ELEVATING LEGREST OPTIONS*			
BQA0576		Calf pad cover, for LH elevating legrest	○
BQA0577		Knee pad for LH elevating legrest	◆
BQA0586		Calf pad cover, for RH elevating legrest	○
BQA0587		Knee pad for RH elevating legrest	◆
BQA3333		Without knee pad, select only with Fixed legrests, Amputee legrests and without legrests	○
LEGREST OPTIONS			
BQA0580		Footplate Extender LH	€67
BQA0590		Footplate Extender RH	€67
BQA0595		Calf Strap (for One piece footrest)	○
CHASSIS COLOUR*			
BQA0591		Pearl Grey	◆
BQA0592		Azurite	○
BQA0593		Sand	○
BQA0594		Electric Green	○
BQA0596		Happy Red	○
CHASSIS SAFETY			
BQA0600		Chassis extender, castors extend 50mm forward for added stability	○

CASTOR ATTACHMENT*			
BQA0610	<input type="checkbox"/>	Fixed castor attachment	◆
FRONT CASTORS*			
BQA0638	<input type="checkbox"/>	140x45 mm (5 1/2" x 1 3/4") Medium Hard Black tyre	€24
BQA0645	<input type="checkbox"/>	150x30 mm (6" x 1 1/4") Black Medium Hard tyre	€24
BQA0665	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre	○
BQA0675	<input type="checkbox"/>	200x45 mm (8" X 1 3/4") Black Solid tyre	◆
REAR WHEEL SIZE*			
BQA0702	<input type="checkbox"/>	16" / 405 mm Transit rear wheel with fixed attachment (NC 0724, 0745, 0790, 0810, 0850, 0863, 0867)	◆
BQA0705	<input type="checkbox"/>	22" / 560 mm Self Propel rear wheel with QR attachment	○
BQA0715	<input type="checkbox"/>	24" / 600 mm Self Propel rear wheel with QR attachment	○
REAR WHEEL TYPE*			
BQA0724	<input type="checkbox"/>	Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes	○
BQA0726	<input type="checkbox"/>	Hub Brake Self Propel or Transit 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes	◆
REAR WHEEL TYRE*			
BQA0745	<input type="checkbox"/>	Pneumatic Profile Black 1" 3/8 tyres	○
BQA0775	<input type="checkbox"/>	Solid Black tyres	◆
HANDRIM			
BQA0790	<input type="checkbox"/>	Aluminium handrim	○
BQA0810	<input type="checkbox"/>	Plastic coated handrim	€58
SPOKEGUARD			
BQA0850	<input type="checkbox"/>	Frost spokeguard	€195
BRAKE*			
BQA0860	<input type="checkbox"/>	User manoeuvred tyre brakes	◆
BQA0863	<input type="checkbox"/>	One Arm Brake LH	€146
BQA0867	<input type="checkbox"/>	One Arm Brake RH	€146
SAFETY OPTIONS			
BQA0870	<input type="checkbox"/>	Pelvic belt	○
TRANSPORTATION BRACKETS			
BQA0875	<input type="checkbox"/>	Tie down labels (Important Note selection of this option will require that a suitable risk assessment is carried out on The final Set up of the chair. Please refer to the Instructions for Use for further information regarding using the Rea Azalea Base as a seat in transit) (NC 0895)	◆
BQA0895	<input type="checkbox"/>	Transportation Brackets, mounted to the front of the wheelchair. (Important note selection of this option will require that a suitable risk assessment is carried out on the final seat up of the chair. Please refer to the Instructions for Use for further information regarding using the Rea Azalea Base as a seat in transit)	€15

OTHER OPTIONS			
BQA0900		Table Tray without elbow pads includes mounting kit (NC 0930)	€164
BQA0902		Table Tray with elbow Pads includes mounting kit (NC 0930)	€164
BQA0910		Swing away table without elbow pads (NC 0920, 0930)	€230
BQA0915		Swing away table with elbow pads (NC 0920, 0930)	€230
BQA0905		Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)	€60
BQA0920		Locking device for table tray (NC 0930)	€16
BQA0930		Mounting kit for table (if ordered separately)	€9
BQA0960		Clinic kit	€64
BQA0970		Harness attachment , must be selected in combination with Head/Neckrest	€34

TOTAL PRICE	
--------------------	--

Non Conforming Options:

BQA0100	+ 100 mm seat width pre-adjustment on armrests and legrests
BQA0110	380-450 mm (15" - 18") only possible with seat width 340mm
BQA0120	430-500 mm (17" - 20")
BQA0164	Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box)
BQA0168	Matrx Elite (includes hardware mounted to frame, backrest supplied in box)
BQA0169	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BQA0172	Matrx Backrest Height 410 mm / 16"
BQA0173	Matrx Backrest Height 510 mm / 20"
BQA0180	Central Legrest, includes calfpad and footrest. Also select options BQA0230, BQA0280, BQA1111 and BQA2222
SWO0189	Standard position with 0-30° Legrest Abduction function
BQA0186	High position 50mm higher for leg length compensation
BQA0188	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BQA0190	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH
BQA0200	Legrest Hanger fixed 80 deg, swing in/swing out LH
BQA0210	Legrest Hanger fixed 90 deg, swing in/swing out LH
BQA0220	Legrest Amputee LH
BQA0240	Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support with locking clip RH
BQA0250	Legrest Hanger fixed 80 deg, swing in/swing out RH
BQA0260	Legrest Hanger fixed 90 deg, swing in/swing out RH
BQA0270	Legrest Amputee RH
BQA0290	Footrest one piece, depth and angle adjustable
BQA0295	Footbox, angle adjustable, padded with short tube
BQA0300	Footrest lower part, fixed with short tube and heel strap, LH
BQA0310	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH
BQA0330	Footrest lower part, fixed with long tube and heel strap, LH.
BQA0340	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH
BQA0360	Footrest lower part, fixed with short tube and heel strap, RH
BQA0370	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH
BQA0390	Footrest lower part, fixed with long tube and heel strap, RH
BQA0400	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH
BQA0430	Armrest Wide Pad LH, height adjustable with autolock Compatible with 16" wheels only
BQA0480	Armrest Wide Pad RH, height adjustable with autolock compatible with 16" wheels only
BQA0520	Side Support Cover, removable with tools (wipeable cover for the sideguard)
BQA0530	Foam for side support cover (internal cushioned pad, which is only possible in combination with BNA0680)
BQA0545	Cover for headrest
BQA0560	Neckrest
BQA0567	Headrest Matrx, 250 mm Soft Pad
BQA0570	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BQA0576	Calf pad cover, for LH elevating legrest
BQA0577	Knee pad for LH elevating legrest
BQA0586	Calf pad cover, for RH elevating legrest
BQA0587	Knee pad for RH elevating legrest
BQA3333	Without knee pad, select only with Fixed legrests, Amputee legrests and without legrests

BQA0580	Footplate Extender LH
BQA0590	Footplate Extender RH
BQA0595	Calf Strap (for One piece footrest)
BQA0600	Chassis extender, castors extend 50mm forward for added stability
BQA0665	200x32 mm (8" x 1 1/4") Black Solid Tyre
BQA0675	200x45 mm (8" X 1 3/4") Black Solid tyre
BQA0702	16" / 405 mm Transit rear wheel with fixed attachment
BQA0705	22" / 560 mm Self Propel rear wheel with QR attachment
BQA0715	24" / 600 mm Self Propel rear wheel with QR attachment
BQA0724	Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes
BQA0726	Hub Brake Self Propel or Transit 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes
BQA0745	Pneumatic Profile Black 1" 3/8 tyres
BQA0790	Aluminium handrim
BQA0810	Plastic coated handrim
BQA0850	Frost spokeguard
BQA0863	One Arm Brake LH
BQA0867	One Arm Brake RH
BQA0895	Transportation Brackets, mounted to the front of the wheelchair. (Important note selection of this option will require that a suitable risk assessment is carried out on the final seat up of the chair. Please refer to the Instructions for Use for further information regarding using the Rea Azalea Base as a seat in transit)
BQA0900	Table Tray without elbow pads includes mounting kit
BQA0902	Table Tray with elbow Pads includes mounting kit
BQA0910	Swing away table without elbow pads
BQA0915	Swing away table with elbow pads
BQA0920	Locking device for table tray
BQA0930	Mounting kit for table (if ordered separately)
BQA0960	Clinic kit



For a more complex customer requirement
please refer to our Invacare Unique special order form.
E-mail: ordersireland@invacare.com or Fax: 01 810
7085

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):