



**Yes, you can:**

Tel: 01656 776222  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

# Manual Invacare Ireland Prescription Form

07 JANUARY 2026  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Rea Azalea Minor

### Standard Features

- Max User Weight 75KG/12 Stone
- Crash Tested ISO 7176-Part 19
- Transit Propulsion with Attendant Hub Brakes and Black Solid Tyres
- Seat Width 340mm (+50mm or +100mm Width Adjustment)
- Choice of 2 Seat Heights 400mm or 450mm
- Depth Adjustable Seat 380mm-450mm
- Backrest Plate with Braced Push Handles
- Assistant Operated Gas Piston Seat Tilt +1° - +24° and Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support
- Height Adjustable Wide Pad Armrests
- Pearl Silver Frame Colour
- User Manoeuvred Tyre Brakes
- Pelvic Belt
- Anti Tippers

**Standard Base price: €3,056**



Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA    special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
BNC0014	Azalea Minor	€3,056
BNC8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
<b>PROPULSION*</b>		
BNC0020	Transit 16" wheels with hub brakes	(NC 1055, 1062, 1074, 1100, 1120) <input checked="" type="radio"/>
BNC0025	Self Propel 24"/22" wheels with user brakes	(NC 0075, 0605, 0635, 1047, 1076, 1305) <input type="radio"/>
BNC0030	Self Propel 24"/22" wheels with assistant hub brakes	(NC 0075, 0605, 0635, 1047, 1074, 1305) €162

SEAT DIMENSIONS*			
BNC0037	<input type="checkbox"/>	Minor (- 50 mm In seat depth and - 30 mm In chassis length)	◆
SEAT WIDTH*			
BNC0038	<input type="checkbox"/>	340 mm/13" (Minor) (NC 0722)	◆
SEAT WIDTH EXTENSION*			
BNC0065	<input type="checkbox"/>	+ 50 mm seat width pre-adjustment on armrests and legrests	○
BNC0075	<input type="checkbox"/>	+ 100 mm seat width pre-adjustment on armrests and legrests (NC 0110)	○
BNC0085	<input type="checkbox"/>	Without seat width extension, set at 0mm adjustment but can be adjusted wider if required	○
SEAT DEPTH*			
BNC0095	<input type="checkbox"/>	380-450 mm / 15" - 18" (Azalea Minor)	◆
SEAT HEIGHT*			
BNC0110	<input type="checkbox"/>	400 mm / 16" (NC 1062)	○
BNC0120	<input type="checkbox"/>	450 mm / 18" (NC 1055)	○
SEAT ANGLE			
BNC0127	<input type="checkbox"/>	Negative Seat Angle -1 To +23° (For attendant operated gas piston in negative tilt)	○
SEAT TYPE*			
BNC0130	<input type="checkbox"/>	Depth Adjustable	◆
SEAT CUSHION*			
BNC0165	<input type="checkbox"/>	Cushion Flosape seat cushion (NC 1111)	€72
BNC0180	<input type="checkbox"/>	Without Seat Cushion, Seat Plate only (NC 0191, 0201, 0202)	◆
SEAT UPHOLSTERY*			
BNC0191	<input type="checkbox"/>	Black Plush TR35	○
BNC0201	<input type="checkbox"/>	Black Multi-stretch Black PU coated fabric	◆
BNC0202	<input type="checkbox"/>	Black Spacer TR36	○
BNC1111	<input type="checkbox"/>	No upholstery, select with option 0180 Without seat cushion	○
BACKREST*			
BNC0217	<input type="checkbox"/>	Flex 2 Backrest Tension adjustable (Azalea Minor)	€204
BNC0218	<input type="checkbox"/>	Flex 2 Backrest Fixed wooden plate (Azalea Minor) (NC 0240)	◆
BACKREST CUSHION*			
BNC0240	<input type="checkbox"/>	Cover for tension adjustable upholstery	○
BNC0250	<input type="checkbox"/>	Cushion, Laguna (Lateral Support)	○
BNC0260	<input type="checkbox"/>	Cushion, Mistral 2 (Waist Support) (NC 0322)	€179
BACKREST UPHOLSTERY*			
BNC0311	<input type="checkbox"/>	Black Plush TR35	○
BNC0321	<input type="checkbox"/>	Black Multi-stretch Black PU coated fabric	◆
BNC0322	<input type="checkbox"/>	Black Spacer TR36	○
PUSH HANDLES*			
BNC0340	<input type="checkbox"/>	Push Handles braced (Azalea Minor)	◆
BNC0345	<input type="checkbox"/>	Push Bar, angle adjustable	€123
SEAT TILT AND BACKREST REGULATION*			
BNC0350	<input type="checkbox"/>	Assistant operated gas piston seat tilt and backrest recline	◆
TILT LOCK*			
BNC1263	<input type="checkbox"/>	Tilt Lock (locking of pre-set angle)	◆

LEGREST HANGER ATTACHMENT*			
BNC0375		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function	€116
LEGREST HANGER LH*			
BNC0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH (NC 0475, 0480, 0490, 0942, 0945, 5555)	◆
BNC0390		Legrest Hanger fixed 80 deg, swing in/swing out LH (NC 0928, 0929, 0930)	○
BNC0410		Legrest Amputee LH (NC 0475, 0476, 0480, 0490, 0510, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956)	€171
LEGREST HANGER RH*			
BNC0430		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH (NC 0475, 0540, 0550, 0942, 0945, 5555)	◆
BNC0440		Legrest Hanger fixed 80 deg, swing in/swing out RH (NC 0931, 0932, 0940)	○
BNC0460		Legrest Amputee RH (NC 0475, 0476, 0540, 0550, 0570, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	€171
FOOTREST LH*			
BNC0475		Footrest one piece, depth and angle adjustable short tube (NC 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940, 0952, 0954, 0956)	€145
BNC0476		Footrest one piece, depth and angle adjustable long tube (NC 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940, 0952, 0954, 0956)	€145
BNC0480		Footrest lower part, fixed with short tube and heel strap, LH (NC 0930, 0941, 0952, 0954, 0956)	○
BNC0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.	○
BNC0510		Footrest lower part, fixed with long tube and heel strap, LH. (NC 0930, 0941, 0952, 0954, 0956)	○
BNC0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNC3333		Without footrest, only select with option 0410 Amputee legrest	○
FOOTREST RH*			
BNC0540		Footrest lower part, fixed with short tube and heel strap, RH (NC 0940, 0941, 0952, 0954, 0956)	○
BNC0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.	○
BNC0570		Footrest lower part, fixed with long tube and heel strap, RH (NC 0940, 0941, 0952, 0954, 0956)	○
BNC0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNC4444		Without footrest, only select with options 0460 Amputee legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BNC0600		Armrest Ducktail Pad LH, height adjustable with autolock (NC 1232)	○
BNC0605		Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only (NC 1232)	◆
BNC0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	€112
BNC0610		Armrest Padded LH, height adjustable with autolock (NC 0680, 0695, 1230, 1231, 1232, 1250, 1260, 9999)	€79
BNC0615		Armrest Hemiplegic LH, height adjustable with autolock (NC 0695, 1230, 1231, 1232, 1250, 1260)	€222
BNC0620		Without Armrest LH Select when choosing BNC1232 Half Tray RH, Narrow armrest is included (NC 0695, 1230, 1231, 1250, 1260)	○

ARMREST RH*				
BNC0630		Armrest Ducktail Pad RH, height adjustable with autolock	(NC 1233)	○
BNC0635		Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only	(NC 1233)	◆
BNC0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock	(NC 1230, 1231, 1233, 1250, 1260)	€112
BNC0640		Armrest Padded RH, height adjustable with autolock.	(NC 0680, 0695, 1230, 1231, 1233, 1250, 1260, 9999)	€79
BNC0645		Armrest Hemiplegic, RH, height adjustable with autolock	(NC 0695, 1230, 1231, 1233, 1240, 1250, 1260)	€222
BNC0650		Without Armrest RH Select when choosing BNC1233 Half Tray LH, Narrow armrest is included	(NC 0695, 1230, 1231, 1240, 1250, 1260)	○
ARMREST OPTIONS				
BNC0680		Side Support Cover, PU fabric cover for the sideguards (Removable with tools)	(NC 0695, 9999)	€63
BNC9999		Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools)	(NC 0695)	€63
BNC0695		Pad for sideguard, removable with hook and loop fastening		€63
SEAT OPTIONS				
BNC0710		Incontinence Cover		€48
BNC0720		Pommel Abduction Cushion	(NC 0722)	€145
BNC0722		2 cm Increased seat height plate		€89
BNC0725		Tilt Scale for seat tilt		€53
HEADREST OPTIONS				
BNC0730		Headrest, Sideway Adjustable	(NC 0750, 0770)	◆
BNC0750		Headrest with cheek support, sideway adjustable	(NC 0740, 0770)	€103
SWO0190		Junior Headrest, sideway adjustable. Pad size 230mm x 100mm	(NC 0730, 0740, 0750, 0770)	€103
BACKREST OPTIONS				
BNC0740		Cover for headrest		€39
BNC0770		Neckrest		○
BNC0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)		€61
BNC0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included		○
BNC0915		Tilt Scale for backrest recline		€53
BNC0920		Wedges (One Pair)		€14
TRUNK SUPPORT LH				
BNC0820		Trunk Support Multi Adjustable with extended height arm LH	(NC 0900)	○
BNC0880		NEW!! Trunk Support Swingaway LH, multi adjustable	(NC 0900)	€132
TRUNK SUPPORT RH				
BNC0830		Trunk Support Multi Adjustable with extended height arm RH	(NC 0900)	○
BNC0890		NEW!! Trunk Support Swingaway RH, multi adjustable	(NC 0900)	€132

ELEVATING LEGREST OPTIONS*			
BNC0928		Calf pad cover, for LH elevating legrest	○
BNC0929		Knee pad for LH elevating legrest	◆
BNC0931		Calf pad cover, for RH elevating legrest	○
BNC0932		Knee pad for RH elevating legrest	◆
BNC5555		Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)	○
LEGREST OPTIONS			
BNC0930		Footplate Extender LH (NC 0941)	€40
BNC0940		Footplate Extender RH (NC 0941)	€40
BNC0941		Cover Footplate, Black Plush Fabric	€66
BNC0942		Calf Strap Padded	€88
BNC0945		Calf Strap (for One Piece Footrest)	○
BODYPOINT ANKLE HUGGERS			
BNC0952		Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)	€189
BNC0954		Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)	€189
BNC0956		Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)	€189
CHASSIS COLOUR*			
BNC0944		Silver Pearl	◆
BNC0946		Azurite (Dark Blue)	○
BNC0947		Sand (Beige)	○
BNC0948		Electric Green	○
BNC0949		Happy Red	○
CASTOR ATTACHMENT*			
BNC0960		Fixed castor attachment	◆
FRONT CASTORS*			
BNC0982		140x45 mm (5 1/2" x 1 3/4") Medium Hard Black tyre	€21
BNC0992		150x30 mm (6" x 1 1/4") Medium Hard Black tyre	◆
REAR WHEEL SIZE*			
BNC1047		16" / 405 mm Transit rear wheel with fixed attachment (NC 1074, 1082, 1100, 1120, 1160)	◆
SWO1401		20" / 510 mm Self Propel rear wheel with QR attachment	€85
BNC1055		22" / 560 mm Self Propel rear wheel with QR attachment	○
BNC1062		24" / 600 mm Self Propel rear wheel with QR attachment	○
REAR WHEEL TYPE*			
BNC1074		Standard self propel with colourless hub, clear anodized rim and Stainless steel spokes	○
BNC1076		Hub Brake 16" Transit or 22" or 24" Self Propel Transit with Black Plastic Hub and Rim Self Propel with colourless hub, clear anodized rim and Stainless steel spokes	○
REAR WHEEL TYRE*			
BNC1082		Pneumatic Profile Black 1 3/8" tyres	○
BNC1092		Solid Black tyres	◆
HANDRIM			
BNC1100		Aluminium handrim	◆
BNC1120		Plastic coated handrim	€61

SPOKEGUARD			
BNC1160		Frost spokeguard	€258
BRAKE*			
BNC1170		User manoeuvred tyre brakes	◆
PELVIC BELT*			
BNC1200		Pelvic belt	◆
TRANSPORTATION BRACKETS*			
BNC1225		Transportation Brackets, mounted to the front of the wheelchair	€16
OTHER OPTIONS			
BNC1230		Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	€173
BNC1231		Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	€173
BNC1232		Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	€138
BNC1233		Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	€138
BNC1240		Swing away table without elbow pads (NC 1237, 1250, 1260)	€315
BNC1235		Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)	€63
BNC1237		Half table tray padded cushion	€63
BNC1250		Locking device for table tray	€15
BNC1260		Mounting kit for table (if ordered separately)	€10
BNC1290		Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)	€66
BNC1300		Harness attachment , must be selected in combination with Head/Neckrest	€98
BNC1305		Front castor locks, includes front tiedown brackets. Allows castors to remain parallel for easy access into confined spaces such as lifts or vehicles	€240
BNC1310		Attachment for positioning belt	€39

TOTAL PRICE			
-------------	--	--	--



Non Conforming Options:	
BNC0075	+ 100 mm seat width pre-adjustment on armrests and legrests
BNC0110	400 mm / 16"
BNC0191	Black Plush TR35
BNC0201	Black Multi-stretch Black PU coated fabric
BNC0202	Black Spacer TR36
BNC1111	No upholstery, select with option 0180 Without seat cushion
BNC0240	Cover for tension adjustable upholstery
BNC0322	Black Spacer TR36
BNC0475	Footrest one piece, depth and angle adjustable short tube
BNC0476	Footrest one piece, depth and angle adjustable long tube
BNC0480	Footrest lower part, fixed with short tube and heel strap, LH
BNC0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.
BNC0510	Footrest lower part, fixed with long tube and heel strap, LH.
BNC0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNC0540	Footrest lower part, fixed with short tube and heel strap, RH
BNC0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.
BNC0570	Footrest lower part, fixed with long tube and heel strap, RH
BNC0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNC0605	Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only
BNC0635	Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only
BNC0680	Side Support Cover, PU fabric cover for the sideguards (Removable with tools)
BNC9999	Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools)
BNC0695	Pad for sideguard, removable with hook and loop fastening
BNC0722	2 cm Increased seat height plate
BNC0730	Headrest, Sideway Adjustable
BNC0750	Headrest with cheek support, sideway adjustable
BNC0740	Cover for headrest
BNC0770	Neckrest
BNC0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNC0928	Calf pad cover, for LH elevating legrest
BNC0929	Knee pad for LH elevating legrest
BNC0931	Calf pad cover, for RH elevating legrest
BNC0932	Knee pad for RH elevating legrest
BNC5555	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)
BNC0930	Footplate Extender LH
BNC0940	Footplate Extender RH
BNC0941	Cover Footplate, Black Plush Fabric
BNC0942	Calf Strap Padded
BNC0945	Calf Strap (for One Piece Footrest)
BNC0952	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)
BNC0954	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)
BNC0956	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)

BNC1047	16" / 405 mm Transit rear wheel with fixed attachment
BNC1055	22" / 560 mm Self Propel rear wheel with QR attachment
BNC1062	24" / 600 mm Self Propel rear wheel with QR attachment
BNC1074	Standard self propel with colourless hub, clear anodized rim and Stainless steel spokes
BNC1076	Hub Brake 16" Transit or 22" or 24" Self Propel Transit with Black Plastic Hub and Rim Self Propel with colourless hub, clear anodized rim and Stainless steel spokes
BNC1082	Pneumatic Profile Black 1 3/8" tyres
BNC1100	Aluminium handrim
BNC1120	Plastic coated handrim
BNC1160	Frost spokeguard
BNC1230	Fixed Table Tray without elbow pads includes mounting kit
BNC1231	Fixed Table Tray with elbow Pads includes mounting kit
BNC1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNC1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNC1240	Swing away table without elbow pads
BNC1235	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)
BNC1237	Half table tray padded cushion
BNC1250	Locking device for table tray
BNC1260	Mounting kit for table (if ordered separately)
BNC1305	Front castor locks, includes front tiedown brackets. Allows castors to remain parallel for easy access into confined spaces such as lifts or vehicles





For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):