



Yes, you can:

Manual Ireland Prescription Form

01 2 9 (2023) 5

3 0 H D 0 H 0 0 U H 0 0 H 0 0 U 0 0 W R 0 0 0 0 0 0 L 0 0 D F 0 D 0 U 0 H 0 0 L 0 H 0 0 0 0 R 0 U 0 0 0 0 R 0 F 0 D 0 O 0 Customer Ref 0 L 0 E 0 X

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Azalea Transit Spec 1

Standard Features

- Max user weight 165kg (2 stone)
- Crash test ISO 7176-Part 19
- Transit Propulsion 16" Black solid tyre attendant hub brakes
- Small, medium or large base without seat width adjustment
- Depth adjustable seat 430mm-500mm/17"-20"
- Flex 3 Backrest
- Black Plush Laguna backrest cushion
- Black Plush Flo-shape seat cushion
- Connected push handles
- Assistant operated gas piston seat and backrest recline tilt lock
- Height adjustable wide pad armrests
- Removable, angle adjustable legrests with calf pads and angle adjustable footplates
- Sideways adjustable headrest
- Trunk supports, multi-adjustable extended height
- 200mm 45mm Black medium hard castors
- User manoeuvre brakes
- Pelvic belt
- Anti-tippers
- Pearl Silver grey frame

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
Rea Azalea Transit		
1536796	W 390mm x D 430-500mm	€3,267
1536797	W 440mm x D 430-500mm	€3,267
1536798	W 490mm x D 430-500mm	€3,267

TOTAL PRICE	
--------------------	--

Non Conforming Options:



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):