



Yes, you can:

Manual Ireland Prescription Form

NOVEMBER 2023
Customer Ref:

Please refer to www.invacare.ie for Local Distributor

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Azalea Transit Spec 3

Standard Features

- Max user weight: 135kg (21.5 stone)
- Crash tested ISO 7176-Part 19
- Transit Propulsion 16" Black solid tyre, attendant hub brakes
- Small, medium or large chassis, +50mm seat width adjustment
- Seat height 450mm/18"
- Depth adjustable seat 430mm-500mm/17"-20"
- Flex 3 Tension adjustable backrest
- Black Dartex Passad 2 backrest cushion
- Black Dartex Flo-shape seat cushion
- Connected push handles
- Assistant operated gas piston seat and backrest recline with tilt lock
- Height adjustable wide pad armrests
- Removable, angle adjustable legrests with calf pads and angle adjustable footplates
- Sideways adjustable headrest
- 200 X 45mm Black medium hard front castors
- User manoeuvred tyre brakes
- Pelvic belt
- Anti-tippers
- Pearl Silver grey frame



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
Rea Azalea Transit		
1542419	<input type="checkbox"/> W 390mm x D 430-500mm	€3,171
1542420	<input type="checkbox"/> W 440mm x D 430-500mm	€3,171
1542421	<input type="checkbox"/> W 490mm x D 430-500mm	€3,171

TOTAL PRICE	
--------------------	--

Non Conforming Options:



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):