



Yes, you can:

# Manual Ireland Prescription Form

NOVEMBER 2023

Customer Ref:

Please refer to [www.invacare.ie](http://www.invacare.ie) for Local Distributor

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Rea Azalea Transit Spec 5

### Standard Features

- Max user weight: 135kg (21.5 stone)
- Crash tested ISO 7176-Part 19
- Transit Propulsion 16" Black solid tyre, attendant hub brakes
- Small, medium or large chassis, without seat width adjustment
- Depth adjustable seat 430mm-500mm/17"-20"
- Seat height 450mm/18"
- Flex 3 Backrest Plate
- Black Dartex Passad 2 backrest cushion
- Black Dartex Flo-shape seat cushion
- Connected push handles
- Assistant operated gas piston seat and backrest recline with tilt lock
- Height adjustable wide pad armrests
- Removable, angle adjustable legrests with calf pads and angle adjustable footplates
- Sideways adjustable headrest
- Trunk supports, multi-adjustable extended height
- 200 X 45mm Black medium hard front castors
- User manoeuvred tyre brakes
- Pelvic belt
- Anti-tippers
- Pearl Silver grey frame



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>Rea Azalea Transit</b>		
1542580	W 390mm x D 430-500mm	€3,171
1542581	W 440mm x D 430-500mm	€3,171
1542582	W 490mm x D 430-500mm	€3,171

**TOTAL PRICE**

**Non Conforming Options:**



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):