



Yes, you can:

Tel: 01656 776222  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

**TDX SP2 MAX**

**Standard Features**



- Max User Weight 180kg (28 stone)
- Crash tested ISO7176 - 19
- Seat widths from 490mm - 610mm
- Seat depths from 560mm - 630mm
- Tension Adjustable backrest
- Armrest with reinforced depth adjustable bracket with short armpads
- Manual settable Seat Tilt
- Manual settable Backrest Recline
- VARI F hangers with fixed footplates
- LiNX DLX-REM110 with standard bracket
- Selection of 5 shroud and rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Stability Technology
- 6 km/h (4mph) 4 Pole Motors
- 8A (200 - 250V) Charger
- 73.5 Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: €5,513**



Q'STRAIT  
QLK-150  
Compatible

Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE MODEL*</b>		
ARH0000	TDX SP2 MAX	€5,513
ARZ0001	Special TDX SP2 MAX (Extended Lead Time). Please provide details in comments section.	POA
<b>USER WEIGHT*</b>		
ARH1924	User weight up to 180 kg	◆
<b>SEATING TYPE*</b>		
ARH1145	MAX Seating System	◆

<b>SEAT WIDTH*</b>			
ARH0378	<input type="checkbox"/>	Seat Width of Seat Frame 550 mm (+100 mm adjustment in arms)	<input type="radio"/>
ARH0380	<input type="checkbox"/>	Seat Width of Seat Frame 610 mm (+100 mm adjustment in arms) (NC 0533)	<input type="radio"/>
<b>SEAT DEPTH*</b>			
ARH0396	<input type="checkbox"/>	560 mm	<input type="radio"/>
ARH0397	<input type="checkbox"/>	595 mm	<input type="radio"/>
ARH0398	<input type="checkbox"/>	630 mm	<input type="radio"/>
<b>SEAT CUSHION*</b>			
ARH1484	<input type="checkbox"/>	MAX Seat Cushion	<input type="radio"/>
ARH1114	<input type="checkbox"/>	Without Seat Cushion (10 mm seat board only)	<input type="radio"/>
<b>ADAPTOR TO REDUCE THE SEAT DEPTH BY 130MM</b>			
ARH0399	<input type="checkbox"/>	Adaptor to shorten the seat depth by 130 mm. (Required if seat width 490 mm is selected in combination with seat depth 595 mm or 630 mm)	<input type="radio"/>
<b>BACKREST UPHOLSTERY*</b>			
ARH1106	<input type="checkbox"/>	Tension Adjustable	<input type="checkbox"/>
<b>UPHOLSTERY COLOUR*</b>			
ARH0278	<input type="checkbox"/>	Black Multi Stretch PU Coated Wipe Down Fabric TR26	<input type="checkbox"/>
<b>ARMRESTS*</b>			
ARH0635	<input type="checkbox"/>	Fixed Armrest, 250mm - 340mm tool adjustable in height, detachable. (Includes reinforced depth adjustable bracket).	<input type="checkbox"/>
<b>ARMPADS*</b>			
ARH0647	<input type="checkbox"/>	Short armpad (325mm)	<input type="checkbox"/>
ARH0649	<input type="checkbox"/>	Long armpad (375mm)	<input type="radio"/>
<b>PUSH BAR*</b>			
ARH1307	<input type="checkbox"/>	Push bar	<input type="checkbox"/>
<b>HEADREST / NECKREST (Adequate headrest option must be selected)*</b>			
ARH1432	<input type="checkbox"/>	Headrest sideways adjustable	<input type="checkbox"/>
ARH1427	<input type="checkbox"/>	Headrest with cheek support	€204
ARH1495	<input type="checkbox"/>	Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	<input type="radio"/>
<b>MANUAL / POWERED SEATING OPTIONS*</b>			
ARH1091	<input type="checkbox"/>	Manually settable Seat Tilt (max user weight 180kg)	<input type="checkbox"/>
ARH1095	<input type="checkbox"/>	Powered Heavy Duty Basic Tilt (0° - 20°). Select either 470mm or 495mm Seat to Ground Height. (max user weight 180kg)	€394
<b>SEAT TO GROUND HEIGHT*</b>			
ARH1024	<input type="checkbox"/>	Fixed seat height 470 mm (18.5")	<input type="radio"/>
ARH1046	<input type="checkbox"/>	Fixed seat height 495 mm (19.5")	<input type="radio"/>
<b>BACKREST RECLINER*</b>			
ARH0958	<input type="checkbox"/>	Manually settable recliner by spindle (0° - 25°)	<input type="checkbox"/>
ARH0947	<input type="checkbox"/>	Electrical recliner (0° - 25°)	€455

LEG REST TYPE*			
ARH0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0504, 0513, 0599) €335
ARH0589		Centre Mounted: LiNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 0504, 0513, 0599) €1,477
ARH0564		Side Mounted: Reinforced hangers, fixed angle (not detachable without tools). Includes calf pads.	○
ARH0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	◆
ARH0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	(NC 0559) €211
ARH0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 0513, 0599, 0559) €602
FOOTPLATE TYPE*			
ARH0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	(NC 0503, 0533, 0531) ○
ARH0512		Fixed Footplate pair	(NC 0533, 0556, 0589) ○
ARH0504		Angle Adjustable Footplate pair	(NC 0589) €66
ARH0513		Angle and depth adjustable Footplate pair	(NC 0589) €107
ARH0599		Angle adjustable flip-back, heavy-duty footboard (Recommended for SW 550mm)	(NC 0503, 0531, 0533, 0556, 0589) €301
LiNX CONTROLS*			
ARH2201		DLX-REM110 (Drive Only)	(NC 1437, 1095, 0947, 0589, 0531, 2219, 2245) ◆
ARH2203		DLX-REM211 (Drive & Actuators)	(NC 1437) ○
ARH2204		DLX-REM216 (Drive, Actuators & Lights)	€137
CONTROL BRACKETRY*			
ARH1207		Standard Bracket RH	○
ARH1208		Standard Bracket LH	○
ARH1707		Swing Away RH Bracket	€119
ARH1708		Swing Away LH Bracket	€119
ARH1709		Matrix Resolve Swing in/out Bracket RH	€281
ARH1710		Matrix Resolve Swing in/out Bracket LH	€281
ARH2227		Bracket RH mounted for Attendant	€58
ARH2228		Bracket LH mounted for Attendant	€58
DUAL CONTROLS			
ARH2220		LiNX DLX-REM050 (Dual Control used as secondary remote)	(NC 1221, 1222) €356
ARH2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1222, 2201) €356
ARH2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2219, 2227, 2228, 1207, 1208) €676
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)			
ARH1207		Dual Control Standard Bracket RH	○
ARH1208		Dual Control Standard Bracket LH	○
ARH1221		IDC Dual Control Detachable Bracket RH	(NC 2219, 2227, 2228) ○
ARH1222		IDC Dual Control Detachable Bracket LH	(NC 2219, 2227, 2228) ○
LiNX 10-WAY SWITCH			
ARH2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	€455

<b>LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)</b>			
ARH2217		Bracket RH, mounted for user	○
ARH2218		Bracket LH, mounted for user	○
ARH2227		Bracket RH, mounted for attendant	€58
ARH2228		Bracket LH, mounted for attendant	€58
<b>LiNX ACCESSORY CABLE</b>			
ARH2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	€90
<b>FURTHER LiNX COMPONENTS / ACCESSORIES</b>			
ARH2302		Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
ARH2300		To be enclosed only	○
ARZ2301		To be connected and programmed (Special)	POA
<b>BUDDY BUTTONS</b>			
ARH1685		Red Buddy Button	€80
ARH1686		Yellow Buddy Button	€80
ARH1687		Green Buddy Button	€80
ARH1690		Blue Buddy Button	€80
<b>EGG SWITCH</b>			
ARH2254		Egg Switch, black	€80
<b>PIKO BUTTONS</b>			
ARH1711		Red Piko Button	€142
ARH1712		Yellow Piko Button	€142
ARH1713		Green Piko Button	€142
ARH1714		Orange Piko Button	€142
ARH1715		Blue Piko Button	€142
ARH1716		Black Piko Button	€142
<b>TOGGLES / BUTTONS</b>			
ARZ2261		Stereo Toggle - direct control for Tilt only	€243
ARZ2262		Stereo Button - direct control for Tilt only	€232
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline	€302
ARZ2264		Four-Way Button - direct control for Tilt & Recline	€323
ARZ2207		Toggle / Button RH Mount	○
ARZ2208		Toggle / Button LH Mount	○
ARZ2225		Without Mounting	○
<b>BODYPOINT JOYSTICK KNOBS</b>			
ARH2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	€65
ARH2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	€65
ARH2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	€65
ARH2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	€65
<b>REMOTE STOP SWITCH</b>			
ARZ2258		Remote Stop Switch (ASL504)	€946

<b>SHROUD COLOUR*</b>			
ARH2640	<input type="checkbox"/>	Obsidian Black Frame Shrouds	○
ARH2641	<input type="checkbox"/>	Ruby Red Frame Shrouds	○
ARH2642	<input type="checkbox"/>	Purple Passion Frame Shrouds	○
ARH2643	<input type="checkbox"/>	Peacock Blue Frame Shrouds	○
ARH2644	<input type="checkbox"/>	Arctic Blue Frame Shrouds	○
<b>RIM COLOUR*</b>			
ARH2650	<input type="checkbox"/>	Obsidian Black Rim Inserts	○
ARH2651	<input type="checkbox"/>	Ruby Red Rim Inserts	○
ARH2652	<input type="checkbox"/>	Purple Passion Rim Inserts	○
ARH2653	<input type="checkbox"/>	Peacock Blue Rim Inserts	○
ARH2654	<input type="checkbox"/>	Arctic Blue Rim Inserts	○
<b>FORKS*</b>			
ARH2631	<input type="checkbox"/>	Single sided fork (New) includes Black castor	€317
ARH2632	<input type="checkbox"/>	Double sided fork	○
<b>CENTRE WHEELS (DRIVE WHEELS)*</b>			
ARH0704	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARH0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam	◆
ARH0706	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Puncture protected, gel-sealed	○
<b>MOTOR TYPE*</b>			
ARH1294	<input type="checkbox"/>	3 km/h (2mph)	○
ARH1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole)	◆
<b>BATTERY*</b>			
ARH1260	<input type="checkbox"/>	73.5 Ah	◆
<b>G TRAC</b>			
ARH1215	<input type="checkbox"/>	G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	€742
<b>CHARGER*</b>			
ARH1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
ARH1313	<input type="checkbox"/>	OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€137
<b>LIGHTS*</b>			
ARH1437	<input type="checkbox"/>	LED Light System with Indicators (30 LUX)	€433
ARH1254	<input type="checkbox"/>	Without light system	◆
<b>TRANSPORTATION KIT</b>			
ARH1450	<input type="checkbox"/>	TEMPORARILY UNAVAILABLE ON TDX MODELS UNTIL FURTHER NOTICE - Preparation for DAHL-Docking station. (NC 0001)	€527
ARH1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory.	◆
<b>FIXATION BELT*</b>			
ARH1329	<input type="checkbox"/>	MAX Lapbelt (adjustable both sides) fitted to base of backrest	◆
<b>REAR VIEW MIRROR</b>			
ARH1415	<input type="checkbox"/>	Rear-view Mirror RH	€38
ARH1339	<input type="checkbox"/>	Rear-view Mirror LH	€38

CANE HOLDER			
-------------	--	--	--

ARH1312		Cane holder RH	€38
ARH1363		Cane holder LH	€38

			<b>TOTAL PRICE</b>	
--	--	--	--------------------	--

**Non Conforming Options:**

ARZ0001	Special TDX SP2 MAX (Extended Lead Time). Please provide details in comments section.
ARH1095	Powered Heavy Duty Basic Tilt (0° - 20°). Select either 470mm or 495mm Seat to Ground Height. (max user weight 180kg)
ARH0947	Electrical recliner (0° - 25°)
ARH0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
ARH0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
ARH0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
ARH0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
ARH0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
ARH0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
ARH0504	Angle Adjustable Footplate pair
ARH0513	Angle and depth adjustable Footplate pair
ARH0599	Angle adjustable flip-back, heavy-duty footboard (Recommended for SW 550mm)
ARH2201	DLX-REM110 (Drive Only)
ARH1207	Standard Bracket RH
ARH1208	Standard Bracket LH
ARH2227	Bracket RH mounted for Attendant
ARH2228	Bracket LH mounted for Attendant
ARH2219	LiNX DLX-ACU200 (Compact Dual Control)
ARH2245	LiNX DLX-IDC (Intuitive Dual Control)
ARH1207	Dual Control Standard Bracket RH
ARH1208	Dual Control Standard Bracket LH
ARH1221	IDC Dual Control Detachable Bracket RH
ARH1222	IDC Dual Control Detachable Bracket LH
ARH2227	Bracket RH, mounted for attendant
ARH2228	Bracket LH, mounted for attendant
ARH1686	Yellow Buddy Button
ARH1687	Green Buddy Button
ARH1690	Blue Buddy Button
ARH2254	Egg Switch, black
ARH1711	Red Piko Button
ARH1712	Yellow Piko Button
ARH1713	Green Piko Button
ARH1714	Orange Piko Button
ARH1715	Blue Piko Button
ARH1716	Black Piko Button
ARH1437	LED Light System with Indicators (30 LUX)



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):