



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Power Invacare Ireland Prescription Form

07 JANUARY 2026
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 NB Modulite

Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash tested ISO7176 - 19
- Highly adjustable Modulite Seating System
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Tension Adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 20°)
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- Standard 80° hangers with fixed footplate
- LiNX DLX-REM110 with standard bracket
- Selection of 5 shroud and rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Stability Technology
- 6 km/h (4mph) 4 Pole Motors
- 8A (200 - 250V) Charger
- 50 Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



Standard Base price: €7,624

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ARM0000	TDX SP2 Narrow Base	€7,624
ARZ0001	Special TDX SP2 Narrow Base (Extended Lead Time). Please provide details in comments section.	POA
USER WEIGHT*		
ARM1920	User weight up to 125kg	<input type="radio"/>
ARM1921	User weight 126kg to 136kg	<input type="radio"/>
ARM1922	User weight 137kg to 150kg (NC 1092, 2052, 2053)	<input type="radio"/>

SEATING TYPE*			
ARM2100		Invacare® Modulite Seat Unit	◆
SEAT FRAME*			
ARM2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 530 mm. Additional 50 mm width by arm rest adjustment.	◆
SEAT WIDTH*			
ARZ1001		Seat width 330 mm. Seat width achieved on a 380 mm seat frame (2121) with inrigged armrests (2183)	(NC 2181, 2193, 2192, 2194) ○
ARM2121		380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.	(NC 2143, 2701, 2703, 2705) ○
ARM2122		430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.	(NC 2703, 2705) ○
ARM2123		480 mm Seat Plate	(NC 2702, 2704, 2706) ○
ARM2124		530 mm Seat Plate	(NC 2174, 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706) ○
ARM2125		550 mm Seat Plate	(NC 2132, 2134, 2142, 2143, 2147, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178, 2168, 2169, 2170, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2152, 2155, 2183, 2197, 2184, 2186, 2188, 2185, 2187, 2189, 1496, 0556, 0589, 2701, 2702, 2703, 2704, 2705, 1404, 1405) ○
BODY FIT KIT			
UNZ2106		Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side)	(NC 2124, 2192, 2194) ○
SEAT DEPTH RAIL LENGTH*			
ARM2129		Standard Depth Package (410mm - 510mm. Simple stepless adjustment from the side)	(NC 2138) ○
ARM2130		Extended Depth Package (510mm - 560mm. Simple stepless adjustment from the side)	(NC 2132, 2134) ○
SEAT DEPTH*			
ARM2132		410 mm	(NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706) ○
ARM2134		460 mm	○
ARM2136		510 mm	○
ARM2138		560 mm	(NC 2143) ○
SEAT OPTION			
ARZSWO		"H" Frame Interface Kit	€507
SEAT CUSHION*			
ARM2142		Matrx Flo-tech Cushion (Lite - 75mm up to 127kg) or (Contour - 100mm up to 191kg)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983) €126
ARM2143		Matrx Flo-tech Contour Visco; (2-way stretch material, 430 mm x 430 mm, up to 125kg)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983) €126
ARM2147		Matrx Libra (90mm - standard up to 136kg or HD up to 227kg). Extended lead time.	€469
ARM2140		Without seat cushion	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983) ○
POSITIONING ACCESORIES FOR MATRX LIBRA			
ARM1980		Tapered Wedges (Pair)	€28
ARM1981		Abductor	€58
ARM1982		Lateral Wedge	€46
ARM1983		Anterior Seat Wedges	€43
BACKREST SUPPORT*			
ARM2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.	(NC 2172, 2173, 2052, 2053) ○
ARM2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.	◆

PUSH HANDLES FOR BACKPOSTS*			
ARM2107		Backposts WITH push handles	◆
ARM2108		Backposts WITHOUT push handles. Push handles will be cut down on the vertical section.	(NC 1332, 2219, 2220, 2227, 2228, 2245) €48
BACKREST UPHOLSTERY / CUSHION*			
ARM2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.	(NC 1445, 1496, 2155, 2230, 2247, 2052, 2053) ◆
ARM2179		Flex3 Tension Adjustable (540 mm or 480 mm)	(NC 2151, 2192, 2193, 2194, 2199, 2155, 2230, 2247, 2128) €182
ARM2172		Flex3 adjustable Backrest Plates (without Backrest Cushion)	(NC 2163, 2128) €252
ARM2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion	(NC 2163, 2128) €411
ARM2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) €1,038
ARM2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) €1,038
ARM2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) €668
ARM2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) €572
ARM2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) €739
ARM2168		Matrx MAC Back - Foam Cushion [Mesh fabric]	(NC 1445) €650
ARM2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]	(NC 1445) €750
ARM2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 1496, 2230, 2247, 2151, 2152, 2155) ○
MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)			
ARM2380		405 mm (16 inch) MAC BH	(NC 1445) ○
ARM2381		455 mm (18 inch) MAC BH	(NC 1445) ○
ARM2382		505 mm (20 inch) MAC BH	(NC 1445) ○
ARM2383		560 mm (22 inch) MAC BH	(NC 1445) ○
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)			
ARM2385		90 mm (3.5 inch) MAC Contour Depth LH	(NC 1445) €130
ARM2386		150 mm (6 inch) MAC Contour Depth LH	(NC 1445) €130
ARM2387		230 mm (9 inch) MAC Contour Depth LH	(NC 1445) €130
MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)			
ARM2389		90 mm (3.5 inch) MAC Contour Depth RH	(NC 1445) €130
ARM2390		150 mm (6 inch) MAC Contour Depth RH	(NC 1445) €130
ARM2391		230 mm (9 inch) MAC Contour Depth RH	(NC 1445) €130
MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)			
ARM2393		Backrest -25 mm (-1 inch) to seat width	(NC 1445) ○
ARM2394		Backrest matches seat width	(NC 1445) ○
ARM2395		Backrest +25 mm (+1 inch) to seat width	(NC 1445) ○
UPHOLSTERY COLOUR*			
ARM2151		Black Nylon Fabric TR18	◆
ARM2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26	(NC 2171) ○
ARM2155		Black Spacer Fabric TR36	(NC 2171, 2179) ○
ARX1002		No Upholstery Colour required - No backrest upholstery required	(NC 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178) ○

ARMRESTS*			
ARM2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm	◆
ARM2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm	€48
ARM2183		Inrigged Armrests. Reduces SW by 25 mm each side	(NC 2185, 2186, 2187, 2188, 2189) €28
ARM2193		Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).	(NC 1445, 2179) €41
ARM2192		Following Armrest (height adjustable from 230mm - 300mm).	(NC 1445, 2179) €200
ARM2194		Following Armrest (height adjustable from 300mm - 360mm).	(NC 1445, 2179) €200
ARMPADS*			
ARM2195		Short Armpad (325mm)	◆
ARM2196		Long Armpad (375mm)	○
ARM2197		Long Ergonomic Armpad (397mm)	(NC 2230, 2247) ○
ARM2199		Soft Armpad (REA Azalea Wide Pad)	(NC 2179) €103
HIP SUPPORTS LH			
ARM2184		Hip Support LH, pad size 140 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) €272
ARM2186		Hip Support LH, pad size 230 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) €272
ARM2188		Hip Support LH, pad size 330 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) €272
HIP SUPPORTS RH			
ARM2185		Hip Support RH, pad size 140 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) €272
ARM2187		Hip Support RH, pad size 230 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) €272
ARM2189		Hip Support RH, pad size 330 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) €272
MOUNTING HARDWARE FOR HIP SUPPORTS			
ARM2164		Mounting hardware for hip supports (6.1" / 155mm high)	○
ARM2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ARM1350		Neckrest depth adjustable	(NC 2247) €38
ARM1432		Headrest sideways adjustable	(NC 2247) €228
ARM1427		Headrest with cheek support	(NC 2247) €239
ARM2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware	(NC 2247) €406
ARM2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware	(NC 2247) €406
ARM2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware	(NC 2247) €548
ARM1496		Non Invacare Headrest but profile mount fitting Flex3 hard back	(NC 2171, 2174, 2175, 2176, 2177, 2178) ○
ARM1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
LATERAL SUPPORT			
ARM1445		Lateral Support with extended arm, multi-adjustable pad	(NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194) €170

MANUAL / POWERED SEATING OPTIONS*			
ARM1091		Manually settable Seat Tilt 0° - 20°	○
ARM1092		Powered "fixed pivot" Tilt 20° (max user weight 136kg) (NC 2201)	€372
ARM1095		Powered "fixed pivot" Tilt 20°, heavy duty (max user weight >136kg) (NC 2201)	€372
ARM1093		Powered Tilt (0° - 30°) with CoG shift (Seat height 460mm) (NC 0354, 0364, 2201)	€374
ARM2052		Powered Tilt (0°-50°) with CoG shift (Seat height 460mm, 136kg max user weight) (NC 0364, 2201)	€1,735
ARM1094		Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (Seat height 460mm - 490mm) (NC 0354, 0364, 2201)	€1,419
ARM2053		Seat Lifter with Powered Tilt (300mm, 45°) with CoG shift (Seat height 460mm, 136kg max user weight) (NC 0364, 2201)	€3,096
SEAT TO GROUND HEIGHT*			
ARM0354		450mm (17.75") (from floor to front of seat plate at 0° tilt) (NC 1094, 2052, 2053)	○
ARM1003		460mm (18") (from floor to front of seat plate at 0° tilt) For GoG Tilt and/or Lift options (NC 1092, 1095)	○
ARM0364		475mm (18.75") (from floor to front of seat plate at 0° tilt) (NC 1094, 2052, 2053, 1091)	○
SEAT TILT OPTION			
ARM2054		5° Fixed Anterior Tilt, Assist mount. Will reduce max tilt angle by 5° (NC 1092, 1095)	€476
BACKREST RECLINER*			
ARM0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆
ARM0977		Powered Recliner (0° - 30°) Foldable (NC 2201)	€338
LEG REST TYPE*			
ARM0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm (NC 0517, 0512, 0504, 0513)	€321
ARM0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm (NC 2201, 0504, 0512, 0504, 0513)	€1,914
ARM0526		Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm (NC 2365, 0559, 0516, 0517)	◆
ARM0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools (NC 0559, 0516, 0517)	○
ARM0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad (NC 2365, 0559, 0516, 0517)	€133
ARM0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad (NC 2201, 2365, 0559, 0516, 0517)	€686
ARM0501		Without Leg Rests (NC 0504, 0513, 2365, 0559, 0516, 0517, 0559, 0516, 0517)	○
FOOTPLATE TYPE*			
ARM2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5") (NC 0503, 0533, 0531)	€129
ARM0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5") (NC 0503, 0533, 0531)	○
ARM0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5") (NC 0503, 0533, 0531)	€167
ARM0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75") (NC 0526, 0503, 0533, 0531)	€376
ARM0512		Fixed Footplate pair (NC 0533, 0531)	◆
ARM0504		Angle Adjustable Footplate pair (NC 0526)	€75
ARM0513		Angle and depth adjustable Footplate pair	€122
ARM0501		Without footplates for option 0501 without Leg Rests (NC 0503, 0533, 0531)	○
AMPUTEE SUPPORT			
ARM0554		Amputee Support RH (NC 0556, 0589)	€202
ARM0555		Amputee Support LH (NC 0556, 0589)	€202

LINX CONTROLS*				
ARM2201		DLX-REM110 (Drive Only)	(NC 1437, 1092, 1095, 1093, 1094, 0977, 0589, 0531, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2219, 2245, 2270, 1908, 2209, 2210)	◆
ARM2203		DLX-REM211 (Drive & Actuators)	(NC 1437, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2270, 1908, 2209, 2210)	○
ARM2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2270, 2209, 2210)	€159
ARM2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	(NC 2209, 2210)	€473
ARM2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 1208, 1707, 1708, 1709, 1710, 2213, 2214, 2227, 2228, 2220, 2270)	€909
REM400 TOGGLE OPTION				
ARM2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2203, 2204, 2206)	€329
CONTROL BRACKETRY*				
ARM1207		Standard Bracket RH		○
ARM1208		Standard Bracket LH		○
ARM1707		Swing Away RH Bracket		€151
ARM1708		Swing Away LH Bracket		€151
ARM1709		Matrx Resolve Swing in/out Bracket RH	(NC 2213, 2275, 2235, 2213, 2241, 1404)	€281
ARM1710		Matrx Resolve Swing in/out Bracket LH	(NC 2214, 2276, 2236, 2214, 2242, 1405)	€281
ARM2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 1405, 1415)	€739
ARM2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 1405, 1339)	€739
ARM2209		REM500 Swing Away RH Multi-Adjustable Bracket	(NC 1415)	○
ARM2210		REM500 Swing Away LH Multi-Adjustable Bracket	(NC 1339)	○
ARM2227		Bracket RH mounted for Attendant	(NC 2220, 2108, 2128)	€71
ARM2228		Bracket LH mounted for Attendant	(NC 2220, 2108, 2128)	€71
ARZ2239		Individual mounting option per request (Special). Requires additional information/specification.		POA
DUAL CONTROLS				
ARM2220		LiNX DLX-REM050 (Dual Control used as secondary remote)	(NC 2108, 2205, 2206, 1221, 1222, 2230, 2247, 2257, 2233, 2240)	€489
ARM2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1222, 2108, 2201)	€489
ARM2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2201)	€748
UNZ9259		mo-vis Scoot Control (Only available via Special)	(NC 2108, 0000, 1221, 1222)	€3,020
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
ARM1207		Dual Control Standard Bracket RH		○
ARM1208		Dual Control Standard Bracket LH		○
ARM1221		IDC Dual Control Detachable Bracket RH	(NC 2220, 2227)	○
ARM1222		IDC Dual Control Detachable Bracket LH	(NC 2220, 2228)	○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification.		POA

COMPACT REMOTES				
ARM2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2220)	€845
ARM2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2220)	€1,199
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€3,430
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€3,560
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€3,060
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€3,285
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€1,820
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€2,630
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€4,160
ARM2233		ASL MEC (Proportional MEC Joystick ASL130)	(NC 2203, 2204, 2220, 2230, 2247, 2257, 2240, 2128)	€4,230
ARM2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133)	(NC 2203, 2204, 2220, 2230, 2247, 2257, 2233, 2128)	€2,683
COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
ARM2207		RH Mount	(NC 2201, 2203, 2204)	○
ARM2208		LH Mount	(NC 2201, 2203, 2204)	○
ARM2213		Centre Mount Swing Away to RH (Nucleus)	(NC 1405, 1415, 2201, 2203, 2204)	€739
ARM2214		Centre Mount Swing Away to LH (Nucleus)	(NC 1405, 1339, 2210, 2203, 2204)	€739
ARM2235		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 1415, 2201, 2203, 2204)	○
ARM2236		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1339, 2201, 2203, 2204)	○
ARM2241		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1405, 1415, 2201, 2203, 2204)	€739
ARM2242		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1405, 1339, 2201, 2203, 2204)	€739
ARZ2215		Individual mounting option per request (Special). Requires additional information/specification.		POA
CHIN CONTROLS				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220)	€3,430
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	€3,560
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220)	€3,060
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	€3,285
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	€1,820
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	€2,630
CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)				
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)		€1,255
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)		€1,820
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)		€3,005
HEAD CONTROLS				
ARM2230		ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2171, 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 2327, 1350, 1432, 1427, 2257, 2240)	€4,497
ARM2247		ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2171, 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 2327, 1350, 1432, 1427, 2257, 2240)	€4,497
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0000, 2220)	€9,135
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0000, 2220)	€5,435

FOOT CONTROL			
UNZ9264	<input type="checkbox"/>	mo-vis Foot Control LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220) €3,800
UNZ9265	<input type="checkbox"/>	mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220) €3,615
FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)			
UNZ9266	<input type="checkbox"/>	Individual mounting option per request (Special). Requires additional information/specification.	POA
TRAY CONTROL			
ARM2257	<input type="checkbox"/>	Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2220, 2230, 2247, 2233, 2240, 2128) €8,205
TRAY CONTROL BRACKETRY (Only applicable if Tray Control is selected)			
ARM2275	<input type="checkbox"/>	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 2203, 2204, 1404, 1405, 1415) €586
ARM2276	<input type="checkbox"/>	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 2203, 2204, 1404, 1405, 1339) €586
ARZ2215	<input type="checkbox"/>	Individual mounting option per request (Special). Requires additional information/specification.	POA
WIRELESS ACCESSORIES FOR ASL CONTROLS			
ARM2249	<input type="checkbox"/>	Wireless Triple Switch Receiver (ASL557-3)	€422
ARM2250	<input type="checkbox"/>	Wireless Mouse Emulator (ASL558)	€1,013
LiNX 10-WAY SWITCH			
ARM2251	<input type="checkbox"/>	10-Way Switch DLX-SW10-A for direct Powered Seating Access	€513
LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)			
ARM2217	<input type="checkbox"/>	Bracket RH, mounted for user	(NC 2213, 2275, 2235, 2241, 1415) ○
ARM2218	<input type="checkbox"/>	Bracket LH, mounted for user	(NC 2214, 2276, 2236, 2242, 1339) ○
ARM2227	<input type="checkbox"/>	Bracket RH, mounted for attendant	(NC 2128, 1221) ○
ARM2228	<input type="checkbox"/>	Bracket LH, mounted for attendant	(NC 2128, 1222) ○
LiNX ACCESSORY CABLE			
ARM2252	<input type="checkbox"/>	Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	€105
LiNX USB CHARGER PORT			
ARM2253	<input type="checkbox"/>	USB Charger Port DLX-USB02-A	(NC 2701, 2702, 2703, 2704, 2705, 2706) €162
LiNX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)			
ARM2207	<input type="checkbox"/>	RH Mount	○
ARM2208	<input type="checkbox"/>	LH Mount	○
REMOTE STOP SWITCH			
ARM2258	<input type="checkbox"/>	Remote Stop Switch (ASL504)	(NC 2203, 2204) €1,019
FURTHER LiNX COMPONENTS / ACCESSORIES			
ARM2302	<input type="checkbox"/>	Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
ARM2300	<input type="checkbox"/>	To be enclosed only	○
ARZ2301	<input type="checkbox"/>	To be connected and programmed (Special)	POA
BUDDY BUTTONS			
ARM1685	<input type="checkbox"/>	Red Buddy Button	€90
ARM1686	<input type="checkbox"/>	Yellow Buddy Button	€90
ARM1687	<input type="checkbox"/>	Green Buddy Button	€90
ARM1690	<input type="checkbox"/>	Blue Buddy Button	€90
EGG SWITCH			
ARM2254	<input type="checkbox"/>	Egg Switch, black	€90

PIKO BUTTONS			
ARM1711		Red Piko Button	€161
ARM1712		Yellow Piko Button	€161
ARM1713		Green Piko Button	€161
ARM1714		Orange Piko Button	€161
ARM1715		Blue Piko Button	€161
ARM1716		Black Piko Button	€161
TOGGLES / BUTTONS			
ARZ2261		Stereo Toggle - direct control for Tilt only	€273
ARZ2262		Stereo Button - direct control for Tilt only	€263
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline	€342
ARZ2264		Four-Way Button - direct control for Tilt & Recline	€365
ARZ2207		Toggle / Button RH Mount	○
ARZ2208		Toggle / Button LH Mount	○
ARZ2225		Without Mounting	○
LINX INPUT MODULE			
UNZ2273		DLX-IN200 Module (Jack port for external power button) (NC 2220)	€303
UNZ9249		DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)	€814
BODYPOINT JOYSTICK KNOBS			
ARM2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	€75
ARM2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	€75
ARM2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	€75
ARM2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	€75
SHROUD COLOUR*			
ARM2640		Obsidian Black Frame Shrouds	○
ARM2641		Ruby Red Frame Shrouds	○
ARM2642		Purple Passion Frame Shrouds	○
ARM2643		Peacock Blue Frame Shrouds	○
ARM2644		Arctic Blue Frame Shrouds	○
RIM COLOUR*			
ARM2650		Obsidian Black Rim Inserts	○
ARM2651		Ruby Red Rim Inserts	○
ARM2652		Purple Passion Rim Inserts	○
ARM2653		Peacock Blue Rim Inserts	○
ARM2654		Arctic Blue Rim Inserts	○
FORKS*			
ARM2631		Single sided fork (New) includes Black castor	€358
ARM2632		Double sided fork	◆
CENTRE WHEELS (DRIVE WHEELS)*			
ARM0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARM0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	◆

MOTOR TYPE*			
ARM1294		3 km/h (2mph)	○
ARM1282		6 km/h (4mph)	◆
ARM1908		10 km/h (6mph) (includes light system and indicators (NC 2203, 1437, 1254)	€527
BATTERY*			
ARM1257		50 Ah Gel	◆
G TRAC			
ARM1215		G-Trac (G-trac ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	€865
CHARGER*			
ARM1343		OptiCharge Charger 8A (200 - 250V)	◆
ARM1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€174
LIGHTS*			
ARM1437		LED Light System with Indicators (30 LUX) (NC 2203, 1908)	€394
ARX1006		Lights Included with Motor Package 10 & 12kmh (NC 1282)	○
ARM1254		Without Light system	◆
TRANSPORTATION KIT*			
ARM1450		TEMPORARILY UNAVAILABLE ON TDX MODELS UNTIL FURTHER NOTICE - Preparation for DAHL-Docking station. (NC 0001)	€593
ARM1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory.	◆
FIXATION BELT*			
ARM1329		Lap Belt with Auto Lock (NC 2702, 2703, 2704, 2705, 2706)	◆
ARM2113		Lap belt with auto lock (flexible mount in side rails) (NC 2701, 2702, 2703, 2704, 2705, 2706)	○
ARM1332		H-Belt (NC 2108)	€75
ARM2701		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	€117
ARM2702		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm) (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	€117
ARM2703		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	€108
ARM2704		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm) (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	€108
ARM2705		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	€108
ARM2706		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm) (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	€108
REAR VIEW MIRROR			
ARM1415		Rear view Mirror RH - Must be selected with 10 km/h Motor (NC 2213, 2217, 2275, 2235, 2241)	○
ARM1339		Rear view Mirror LH - Must be selected with 10 km/h Motor (NC 2214, 2218, 2276, 2236, 2242)	○
BAG HOOKS			
ARM2128		Bag Hooks (Max load up to 6kg = impact on max user weight) (NC 2173, 2227, 2228, 2257, 2233, 2240)	€95
CANE HOLDER			
ARM1312		Cane holder RH	€25
ARM1363		Cane holder LH	€25
TRAY			
ARM1404		Tray, Swing Away RH (NC 1707, 2209, 2213, 2214, 2235, 2236, 2241, 2242, 2207, 2275, 2276)	€373
ARM1405		Tray, Swing Away LH (NC 1708, 2210, 2213, 2214, 2235, 2236, 2241, 2242, 2208, 2275, 2276)	€373

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARM0000	TDX SP2 Narrow Base
ARZ0001	Special TDX SP2 Narrow Base (Extended Lead Time). Please provide details in comments section.
ARM2121	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.
ARM2122	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.
ARM2123	480 mm Seat Plate
ARM2124	530 mm Seat Plate
ARM2132	410 mm
ARM2134	460 mm
ARM2138	560 mm
ARM2142	Matrx Flo-tech Cushion (Lite - 75mm up to 127kg) or (Contour - 100mm up to 191kg)
ARM2143	Matrx Flo-tech Contour Visco; (2-way stretch material, 430 mm x 430 mm, up to 125kg)
ARM2147	Matrx Libra (90mm - standard up to 136kg or HD up to 227kg). Extened lead time.
ARM1980	Tapered Wedges (Pair)
ARM1981	Abductor
ARM1982	Lateral Wedge
ARM1983	Anterior Seat Wedges
ARM2163	Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.
ARM2108	Backposts WITHOUT push handles. Push handles will be cut down on the vertical section.
ARM2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
ARM2179	Flex3 Tension Adjustable (540 mm or 480 mm)
ARM2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ARM2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
ARM2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen
ARM2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARM2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARM2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARM2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARM2168	Matrx MAC Back - Foam Cushion [Mesh fabric]
ARM2169	Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]
ARM2170	Without Backrest Upholstery (Back Posts/Back Frame only)
ARM2380	405 mm (16 inch) MAC BH
ARM2381	455 mm (18 inch) MAC BH
ARM2382	505 mm (20 inch) MAC BH
ARM2383	560 mm (22 inch) MAC BH
ARM2385	90 mm (3.5 inch) MAC Contour Depth LH
ARM2386	150 mm (6 inch) MAC Contour Depth LH
ARM2387	230 mm (9 inch) MAC Contour Depth LH
ARM2389	90 mm (3.5 inch) MAC Contour Depth RH

ARM2390	150 mm (6 inch) MAC Contour Depth RH
ARM2391	230 mm (9 inch) MAC Contour Depth RH
ARM2393	Backrest -25 mm (-1 inch) to seat width
ARM2394	Backrest matches seat width
ARM2395	Backrest +25 mm (+1 inch) to seat width
ARM2151	Black Nylon Fabric TR18
ARM2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
ARM2155	Black Spacer Fabric TR36
ARM2181	T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm
ARM2183	Inrigged Armrests. Reduces SW by 25 mm each side
ARM2193	Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).
ARM2192	Following Armrest (height adjustable from 230mm - 300mm).
ARM2194	Following Armrest (height adjustable from 300mm - 360mm).
ARM2196	Long Armpad (375mm)
ARM2197	Long Ergonomic Armpad (397mm)
ARM2199	Soft Armpad (REA Azalea Wide Pad)
ARM2184	Hip Support LH, pad size 140 x 100 mm
ARM2186	Hip Support LH, pad size 230 x 100 mm
ARM2188	Hip Support LH, pad size 330 x 100 mm
ARM2185	Hip Support RH, pad size 140 x 100 mm
ARM2187	Hip Support RH, pad size 230 x 100 mm
ARM2189	Hip Support RH, pad size 330 x 100 mm
ARM1350	Neckrest depth adjustable
ARM1432	Headrest sideways adjustable
ARM1427	Headrest with cheek support
ARM2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
ARM2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware
ARM2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware
ARM1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
ARM1445	Lateral Support with extended arm, multi-adjustable pad
ARM1091	Manually settable Seat Tilt 0° - 20°
ARM1092	Powered "fixed pivot" Tilt 20° (max user weight 136kg)
ARM1095	Powered "fixed pivot" Tilt 20°, heavy duty (max user weight >136kg)
ARM1093	Powered Tilt (0° - 30°) with CoG shift (Seat height 460mm)
ARM2052	Powered Tilt (0°-50°) with CoG shift (Seat height 460mm, 136kg max user weight)
ARM1094	Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (Seat height 460mm - 490mm)
ARM2053	Seat Lifter with Powered Tilt (300mm, 45°) with CoG shift (Seat height 460mm, 136kg max user weight)
ARM0354	450mm (17.75") (from floor to front of seat plate at 0° tilt)
ARM0364	475mm (18.75") (from floor to front of seat plate at 0° tilt)
ARM0977	Powered Recliner (0° - 30°) Foldable
ARM0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
ARM0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
ARM0526	Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm

ARM0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
ARM0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
ARM0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
ARM2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
ARM0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
ARM0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)
ARM0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
ARM0512	Fixed Footplate pair
ARM0504	Angle Adjustable Footplate pair
ARM0513	Angle and depth adjustable Footplate pair
ARM2201	DLX-REM110 (Drive Only)
ARM2203	DLX-REM211 (Drive & Actuators)
ARM2204	DLX-REM216 (Drive, Actuators & Lights)
ARM2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
ARM2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
ARM2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
ARM1208	Standard Bracket LH
ARM1707	Swing Away RH Bracket
ARM1708	Swing Away LH Bracket
ARM1709	Matrx Resolve Swing in/out Bracket RH
ARM1710	Matrx Resolve Swing in/out Bracket LH
ARM2213	Centre Mount Swing Away to RH Nucleus Bracket
ARM2214	Centre Mount Swing Away to LH Nucleus Bracket
ARM2209	REM500 Swing Away RH Multi-Adjustable Bracket
ARM2210	REM500 Swing Away LH Multi-Adjustable Bracket
ARM2227	Bracket RH mounted for Attendant
ARM2228	Bracket LH mounted for Attendant
ARM2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
ARM2219	LiNX DLX-ACU200 (Compact Dual Control)
ARM2245	LiNX DLX-IDC (Intuitive Dual Control)
ARM1208	Dual Control Standard Bracket LH
ARM1221	IDC Dual Control Detachable Bracket RH
ARM1222	IDC Dual Control Detachable Bracket LH
ARM2231	LiNX DLX-CR400 (Compact Remote)
ARM2232	LiNX DLX-CR400LF (Compact Remote Low Force)
ARM2233	ASL MEC (Proportional MEC Joystick ASL130)
ARM2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
ARM2207	RH Mount
ARM2208	LH Mount
ARM2213	Centre Mount Swing Away to RH (Nucleus)
ARM2214	Centre Mount Swing Away to LH (Nucleus)

ARM2235	Tray Mount RH (ASL603, 100mm x 150mm)
ARM2236	Tray Mount LH (ASL603, 100mm x 150mm)
ARM2241	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARM2242	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARM2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
ARM2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
ARM2257	Four-Switch Proximity Array (ASL106)
ARM2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
ARM2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
ARM2217	Bracket RH, mounted for user
ARM2218	Bracket LH, mounted for user
ARM2227	Bracket RH, mounted for attendant
ARM2228	Bracket LH, mounted for attendant
ARM2253	USB Charger Port DLX-USB02-A
ARM2207	RH Mount
ARM2208	LH Mount
ARM1686	Yellow Buddy Button
ARM1687	Green Buddy Button
ARM1690	Blue Buddy Button
ARM2254	Egg Switch, black
ARM1711	Red Piko Button
ARM1712	Yellow Piko Button
ARM1713	Green Piko Button
ARM1714	Orange Piko Button
ARM1715	Blue Piko Button
ARM1716	Black Piko Button
ARZ2207	Toggle / Button RH Mount
ARZ2208	Toggle / Button LH Mount
ARM1282	6 km/h (4mph)
ARM1908	10 km/h (6mph) (includes light system and indicators)
ARM1437	LED Light System with Indicators (30 LUX)
ARM1254	Without Light system
ARM1329	Lap Belt with Auto Lock
ARM2113	Lap belt with auto lock (flexible mount in side rails)
ARM1332	H-Belt
ARM2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ARM2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ARM2703	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ARM2704	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ARM2705	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)
ARM2706	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ARM1415	Rear view Mirror RH - Must be selected with 10 km/h Motor
ARM1339	Rear view Mirror LH - Must be selected with 10 km/h Motor
ARM2128	Bag Hooks (Max load up to 6kg = impact on max user weight)

ARM1404	Tray, Swing Away RH
---------	---------------------

ARM1405	Tray, Swing Away LH
---------	---------------------



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):