



Powerchairs Ireland

Prescription Form

Yes, you can:

Tel: 01 810 7084 Fax: 01 810 7085

Jan-22

Email: ordersireland@invacare.com

Customer Ref:

Online Spares available at: www.invacare.ie

Customer Ref :	Account No :
Contact no:	Order no:
Delivery Address:	Order Qty:
	Order Date :

TDX2 Ultra NB

Standard Features

Max user weight 136Kg (21.4 Stone)
 Crash test ISO7176-Part 19
 Ultra Low Maxx seating by Motion Concepts
 Compact, manoeuvrable, CWD TDX SP2 chassis
 Seat widths of 405mm - 610mm (16" - 24"), adjustable between 405mm - 505mm (16" - 20") and 505mm - 610mm (20" - 24")
 Seat depth of 380mm - 580mm (15" - 23") adjustable between 380mm - 480mm (15" - 19") and 480mm - 580mm (19 - 23")
 Backrest heights of 560mm - 660mm in 25mm increments
 Cantilever flip-back arm, width, height and angle adjustable with desk armpad with startex covered visco foam
 Elan multi positioning headrest
 ULTRA Low Maxx CG power tilt (0° to 50°) seat heights of 420mm - 480mm
 Powered Recline 90° to 170° with extended shear reduction (ESR)
 Vari F hangers with fixed footplates
 LiNX control REM211 with standard bracket
 Selection of 10 stunning shroud and Rim colours
 Centre (driving wheels) 355mm x 76mm (14" x 3") Black tyres
 10km/h (6mph) 4 pole motors and 50 Ah batteries
 G-Trac
 4 Point tie down system
 Lapbelt with autolock (adjustable both sides)



Please make your selection with a cross in the box			☐ Free Of Charge	☒ Standard Feature	☐ POA Special Request
Customer weight=	KG				
Code	Description				
BASE PRODUCT *					
O0001	☐	TDX SP2 ULTRA NARROW BASE	☐		
CHAIR TYPE *					

ARP0080		TDX SP2 ULTRA Narrow Base	◆
ARX0080		Special TDX SP2 ULTRA Narrow Base (Please provide details in special section)	POA
SEATING TYPE *			
ARP2020		Ultra Low Maxx Powered Positioning System	€13,186
SEAT ADJUSTMENTS *			
ARP2052		Centre of gravity Tilt (0 to 50°)	◆
ARP2053		300mm seat lifter (increased STF height by 4mm) with centre of gravity tilt (0 to 45°) (NC 0339)	€2,203
SEAT HEIGHT SETTING *			
ARP0337		Seat-to -Floor Height 440mm	○
ARP0331		Seat-to Floor Height 460mm	○
ARP0339		Seat-to Floor Height 480mm (NC 2053)	○
ANTERIOR AND POSTERIOR TILT *			
ARP2057		Tilt Standard Mount (No anterior or posterior tilt)	◆
ARP2054		5° Fixed Anterior Assist Mount (will reduce max tilt angle by 5°)	€298
ARP2055		10° Fixed Anterior Assist Mount (will reduce max tilt angle by 10°)	€372
ARP2056		5° Pre-Tilt (will increase max tilt angle by 5°)	€298
BACKREST RECLINER *			
ARP2060		Powered Recliner (90° - 170°) with Extended Shear Reduction (ESR), Standard Mount	◆
ARP2063		Powered Recliner 8° Recline mount 82° to 162° with Extended Shear Reduction (ESR)	€32
ARP2066		Powered Recline, 30° Recline mount, 60° to 140° with Extended Shear Reduction (ESR)	€32
SEAT CUSHION *			
ARP2144		Matrx Vi Cushion (Mesh fabric)	€315
ARP2146		Matrx Flo-tech Contour Visco	€158
ARP2147		Matrx Libra 90mm thick	€483
ARP2140		Solid seat Base, For alternative cushion options	◆
BACKREST SUPPORT *			
ARP2050		High Back, height 410mm (Mesh fabric) (NC 2073, 2074)	◆
ARP2051		High Back height 460mm (Mesh fabric) (NC 2070, 2071)	○
ARP2320		High Back, height 505mm / 20" (Mesh fabric)	○
ARP2550		Matrx Elite Back, height 355mm / 14"	€519
ARP2551		Matrx Elite Back, height 405mm / 15"	€519

ARP2552		Matrx Elite Back, height 455mm / 18"	€519
ARP2553		Matrx Elite Back, height 505mm / 20"	€519
ARP2554		Matrx Elite HD Back, height 405mm / 16"	€672
ARP2555		Matrx Elite HD Back, height 505mm / 20"	€672
ARP2556		Matrx Elite TR Back, height 405mm / 16" (Mesh fabric)	€519
ARP2557		Matrx Elite TR Back, height 455mm / 18" (Mesh fabric)	€519
ARP2177		Matrx Elite TR Back height 505mm (Mesh fabric)	€519
ARP2563		Matrx Elite TR HD Back, height 505mm / 20" (Mesh fabric)	€672
ARP2558		Matrx Elite Deep Back, height 405mm / 16" (Mesh fabric)	€672
ARP2559		Matrx Elite Deep Back, height 455mm / 18" (Mesh fabric)	€672
ARP2178		Matrx Elite Deep Back, height 500mm (Mesh fabric)	€672
ARP2174		Matrx Posture Back, height 405mm (Mesh fabric max dimension 455mm wide)	€519
ARP2175		Matrx Posture Back height 505mm (Mesh fabric max dimension 455mm wide)	€519
ARP2560		Matrx PB Deep Back, height 405mm / 16"	€672
ARP2561		Matrx PB Deep Back, height 505mm / 20"	€672
ARP2562		Matrx Elite Deep HD Back, height 505mm / 20"	€672
BACKREST HEIGHT *			
ARP2506		505mm (20") measured from seat pan to top of back shell <i>(NC 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 21)</i>	○
ARP2501		535mm (21") measured from seat pan to top of back shell <i>(NC 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 21)</i>	○
ARP2070		560mm (22") measured from seat pan to top of back shell <i>(NC 2050, 2551, 2553, 2554, 2555, 2177, 2558, 2178, 2174, 21)</i>	○
ARP2071		585mm (23") measured from seat pan to top of back shell <i>(NC 2050, 2551, 2553, 2554, 2555, 2177, 2558, 2178, 2174, 21)</i>	○
ARP2072		610mm (24") measured from seat pan to top of back shell <i>(NC 2050, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 25)</i>	○
ARP2073		635mm (25") measured from seat pan to top of back shell <i>(NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2558, 2559, 21)</i>	○
ARP2074		660mm (26") measured from seat pan to top of back shell <i>(NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 25)</i>	○

ARP2502		685mm (27") measured from seat pan to top of back shell (NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 25	○
ARMREST *			
ARP2095		Cantilever Flip-Back Arm, width, height and angle adjustable (230mm - 330mm) (9" - 13) (NC 2308, 2309)	◆
ARP2096		Cantilever Flip-Back Arm, width, height and angle adjustable (320mm - 405mm) (12.5" - 16") (NC 2308, 2309)	○
ARP2097		Two-post Recline Arm, Flip-Back, adjustable height (240mm - 330mm)	○
ARP2098		Two-post Recline Arm, Flip-Back, adjustable height (330mm - 405mm)	○
ARP0690		Without Armrest	○
ARMREST ACCESSORIES			
ARP2308		1" Offset Armrest Spacer (pair), for two-post recline only (NC 2095, 2096)	€123
ARP2309		1/2" Offset Armrest Spacer (pair), for two recline only (NC 2095, 2096)	€123
ARMPAD LH *			
ARP2012		Desk arm pad with startex cover visco foam (255mm) 40m thick	◆
ARP2013		Full length armpad with startex cover foam (355mm) 40mm thick	○
ARP2014		Ergonomic arm though 30mm thick (NC 1404, 1405, 2347, 2348, 2349, 2350)	€100
ARP2015		Without LH arm pad	○
ARMPAD RH *			
ARP2016		Desk arm pad with startex cover visco foam (255mm) 40mm thick	◆
ARP2017		Full length armpad with startex cover foam (355mm) 40mm thick	○
ARP2018		Ergonomic arm though 30mm thick (NC 1404, 1405, 2351, 2352, 2353, 2354)	€100
ARP2019		Without RH arm pad	○
ARMPAD ACCESSORIES & MOUNTING OPTIONS LH			
ARP2080		Flat Hand Pad Left (in combination with arm trough only): to be mounted on the opposite side of remote	€100
ARP2081		Ergonomic Mounting Hardware to mount Otto Bock (EMO) Hand pad options Left (in combination with arm trough only): to be mounted on the opposite side of remote	€44
ARMPAD ACCESSORIES & MOUNTING OPTIONS RH			
ARP2083		Flat Hand Pad Right (in combination with arm trough only): to be mounted on the opposite side of remote	€100
ARP2084		Ergonomic Mounting Hardware to mount Otto Bock (EMO) Hand pad options Right (in combination with arm trough only): to be mounted on the opposite side of remote	€44
MULTI-AXIX UPPER EXTREMITY SUPPORTS LH			
ARP2082		Multi Axis Upper Extremity Support (MACES) Left, raises arm pad by 43mm	€252
MULTI-AXIX UPPER EXTREMITY SUPPORTS RH			
ARP2085		Multi Axis Upper Extremity Support (MACES) Right, raises arm pad by 43mm	€252

ELBOW BLOCK LH PADS			
ARP2347		Elbow Block XS (110mm wide x 83mm high) (NC 2014)	€472
ARP2348		Elbow Block S (135mm wide x 90mm high) (NC 2014)	€472
ARP2349		Elbow Block M (145mm wide x 95mm high) (NC 2014)	€472
ARP2350		Elbow Block L (165mm wide x 110mm high) (NC 2014)	€472
ELBOW BLOCK LH MOUNTING HARDWARE			
ARP2361		Vertical 3" / Horizontal 5"	○
ARP2362		Vertical 3" / Horizontal 7"	○
ARP2363		Vertical 5" / Horizontal 5"	○
ARP2364		Vertical 5" / Horizontal 7"	○
ELBOW BLOCKS RH PADS			
ARP2351		Elbow Block XS (110mm wide x 83mm high) (NC 2018)	€472
ARP2352		Elbow Block S (135mm wide x 90mm high) (NC 2018)	€472
ARP2353		Elbow Block M (145mm wide x 95mm high) (NC 2018)	€472
ARP2354		Elbow Block L (165mm wide x 110mm high) (NC 2018)	€472
ELBOW BLOCK RH MOUNTING HARDWARE			
ARP2371		Vertical 3" / Horizontal 5"	○
ARP2372		Vertical 3" / Horizontal 7"	○
ARP2373		Vertical 5" / Horizontal 5"	○
ARP2374		Vertical 5" / Horizontal 7"	○
HIP SUPPORTS LH (Not available with 410mm seat depth and Manual Centre legrests)			
ARP2184		Hip Support LH, pad size 140 x 100 mm	€233
ARP2114		Hip Support LH, pad size 180 x 150 mm	€233
ARP2186		Hip Support LH, pad size 230 x 100 mm	€233
ARP2188		Hip Support LH, pad size 330 x 100 mm	€233
HIP SUPPORTS RH (Not available with 410mm seat depth and Manual Centre legrests)			
ARP2185		Hip Support RH, pad size 140 x 100 mm	€233
ARP2115		Hip Support RH, pad size 180 x 150 mm	€233
ARP2187		Hip Support RH, pad size 230 x 100 mm	€233
ARP2189		Hip Support RH, pad size 330 x 100 mm	€233
MOUNTING HARDWARE FOR HIP SUPPORTS			
ARP2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
HEADREST *			
ARP2092		Onyx Auto Style headrest (NC 2255, 2230, 2247)	€114
ARP2323		Elan Standard Pad 150mm / 6" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	€303
ARP2324		Elan Standard Pad 255mm / 10" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	◆

ARP2322		Elan Standard Pad 355mm / 14" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	€303
ARP2325		Elan Occipital Pad 230mm x 180mm / 9" x 7" w/standard cover (EOP23) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	€303
ARP2326		Elan Occipital Pad 305mm x 200mm / 12" x 8" w/standard cover complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	€303
ARP2327		Elan 4 Point Headrest Pad 280mm x 255mm / 11" x 10" w/standard cover complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	€408
ARP1495		Head Control or equivalent Headrest, provided independently from this order (ACCORDING TO THE REQUIREMENTS OF ISO 7176-19, A HEADREST IS COMPULSORY FOR USING THE WHEELCHAIR AS VEHICLE SEAT) (NC 2229)	○
LATERAL SUPPORTS LH			
ARP2343		110mm wide x 83mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€365
ARP2345		135mm wide x 90mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€365
ARP2043		145mm wide x 95mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€365
ARP2045		165mm wide x 110mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€365
LATERAL SUPPORTS RH			
ARP2344		110mm wide x 83mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€365
ARP2346		135mm wide x 90mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€365
ARP2044		145mm wide x 95mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€365
ARP2046		165mm wide x 110mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€365
LATERAL SUPPORT ACCESSORIES			
ARP2049		Additional telescopic link for Lateral support, 25mm, pair (NC 2558, 2559, 2178, 2560, 2561, 2562)	€145
CENTRE MOUNTED LEGREST			
ARP0586		Hanger, centre mounted, angle settable (0° - 90°): knee to heel length 310mm - 410mm	€294
ARP0589		Hangers, centre mounted, length compensating, powered, LNX (0° - 90°): (NC 2261, 2262, 2263, 2264) KNEE TO HEL LENGTH 340mm - 410mm	€1,637
FOOTPLATES			
ARP2365		Angle & width adjustable footplates Small (127 width x 190 depth / 5" x 7.5")	€96
ARP0559		Angle & width adjustable footplate (240mm depth)	○
ARP0516		Angle & width adjustable footplates (290mm depth)	€125
FOOT PLATFORM			
ARP2366		Angle adjustable foot platform Small (305mm width x 215mm depth / 12" x 8.5") (NC 0586)	€249

ARP0517		Angle adjustable foot platform (240mm depth, 305mm width) (NC 0586)	€249
SIDE MOUNTED HANGERS LH			
ARP0536		Vari F: Angle adjustable with tools 0° - 70°, length 290mm - 460mm (NC 0505)	◆
ARP0537		Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm (NC 0505)	€122
ARP0538		ADM MANUAL: Angle adjustable with quick release lever 0° - 80°, length 290mm - 460mm (NC 0505)	€116
ARP0539		ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm (NC 0505, 2261, 2262, 2263, 2264)	€302
ARP0572		Pivot Plus Manual Elevating Legrest (NC 0505, 0509)	€655
ARP0577		Pivot Plus Power Elevating Legrest (NC 0505, 0509, 2261, 2262, 2263, 2264)	€1,070
ARP0535		Without Hanger and footplates (NC 0505, 0509, 0511)	○
SIDE MOUNTED HANGERS RH			
ARP0546		Vari F: Angle adjustable with tools 0° - 70°, length 290mm - 460mm (NC 0502)	◆
ARP0547		Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm (NC 0502)	€122
ARP0548		ADM MANUAL: Angle adjustable with quick release lever 0° - 80°, length 290mm - 460mm (NC 0502)	€116
ARP0549		ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm (NC 0502, 2261, 2262, 2263, 2264)	€302
ARP0579		Pivot Plus Manual Elevating Legrest (NC 0502, 0506)	€655
ARP0584		Pivot Plus Power Elevating Legrest (NC 0502, 0506, 2261, 2262, 2263, 2264)	€1,070
ARP0545		Without Hanger and footplates (NC 0502, 0506, 0508)	○
INDIVIDUAL FOOTPLATE TYPE LH			
ARP0505		Fixed Footplate (NC 0535, 0537, 0538, 0539, 0572, 0577)	◆
ARP0509		Angle Adjustable Footplate (NC 0535, 0572, 0577)	€65
ARP0511		Angle and depth adjustable Footplate (NC 0535)	○
INDIVIDUAL FOOTPLATE TYPE RH			
ARP0502		Fixed Footplate (NC 0545, 0547, 0548, 0549, 0579, 0584)	◆
ARP0506		Angle Adjustable Footplate (NC 0545, 0579, 0584)	€65
ARP0508		Angle and depth adjustable Footplate (NC 0545)	○

CONTROLS *				
ARP2203		REM211 Drive & Actuators	(NC 1437, 2231, 2232, 2229, 2233, 2230, 2247, 2257, 2240, 22	○
ARP2204		REM216 Drive, Actuators & Lights	(NC 2231, 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2243, 22	◆
ARP2205		REM400, touch-screen Interface, built in Bluetooth HID	(NC 2209, 2210, 2213, 2215, 2075, 2076)	€449
ARP2206		REM500, without joystick for alternative access, built in Bluetooth HID	(NC 1207, 1208, 1707, 1708, 2227, 2228, 2075, 2076)	€825
CONTROL BRACKETRY *				
ARP1207		Standard Bracket RH	(NC 1404, 2311)	○
ARP1208		Standard Bracket LH	(NC 1405, 2310)	○
ARP1707		Swing Away + height adjustable Bracket RH High end	(NC 1404, 2311)	€122
ARP1708		Swing Away + height adjustable Bracket LH High end	(NC 1405, 2310)	€122
ARP2075		Swing Away + height / width adjust mount on RH for user, Quad link design		€207
ARP2076		Swing Away + height / width adjust mount on LH for user, Quad link design		€207
ARP2227		Centre Mount swing-away Bracket to RH (Nucleus)	(NC 2219, 2245)	€65
ARP2228		Centre Mount swing-away Bracket to LH (Nucleus)	(NC 2219, 2245)	€65
ARP2209		REM500 Swing away RH multi-adjustable	(NC 1404, 2311)	○
ARP2210		REM500 Swing away LH multi-adjustable	(NC 1405, 2310)	○
ARP2213		REM500 Centre mount swing away to RH Nucleus	(NC 1404, 1405, 2310, 2311)	€670
ARP2214		REM500 Centre mount swing away to LH Nucleus	(NC 1404, 1405, 2310, 2311)	€670
ARP2075		Swing Away + height / width adjustable Bracket Mount on RH, Quad link design		€207
ARP2076		Swing Away + height / width adjustable Bracket Mount on RH Quad link design		€207
ARX2239		Individual mounting option per (special) - Requires additional information / specification		POA
REM400 TOGGLE OPTION				
ARP2270		Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial)	(NC 2203, 2204)	€298
CHIN CONTROL				

ARP2229		Proportional Chin Joystick (ASL138) (Requires Primary Remote REM400 or REM500) (Please select chin control below) <i>(NC 1495, 2203, 2204, 2230, 2247, 2233, 2257, 2240)</i>	€2,817
CHIN CONTROL BRACKETRY			
ARP1289		Chin Control Swing Away Bracket RH Manual <i>(NC 1495, 2203, 2204, 2230, 2247, 2233, 2257, 2240)</i>	○
ARP1288		Chin Control Swing Away Bracket LH Manual <i>(NC 1495, 2203, 2204, 2230, 2247, 2233, 2257, 2240)</i>	○
DUAL CONTROL			
ARP2219		Compact Dual Control (DLX-ACU200) (Please select Dual control bracket below) <i>(NC 1221, 1222, 2227, 2228, 2245)</i>	€364
DUAL CONTROL BRACKETRY			
ARP1207		Dual Control Standard Bracket RH <i>(NC 1221, 1222, 2227, 2228, 2245)</i>	○
ARP1208		Dual Control Standard Bracket LH <i>(NC 1221, 1222, 2227, 2228, 2245)</i>	○
INVACARE LINX IDC DUAL CONTROL			
ARP2245		IDC Dual control (Please select AQE1221 OR AQE1222 for bracket option below) <i>(NC 2219, 2227, 2228)</i>	€678
INVACARE LINX IDC DUAL CONTROL BRACKETRY			
ARP1221		IDC Dual control Detachable Bracket RH <i>(NC 2219, 2227, 2228)</i>	○
ARP1222		IDC Dual control Detachable Bracket LH <i>(NC 2219, 2227, 2228)</i>	○
LINX COMPACT REMOTE			
ARP2231		Proportional Compact Remote DLX-CR400 (Please select Proportional Bracket) <i>(NC 2203, 2204)</i>	€767
ARP2232		Proportional Compact Remote Low Force DLX-CR400LF (Please select Proportional Bracketry) <i>(NC 2203, 2204)</i>	€1,088
LINX COMPACT REMOTE BRACKET			
ARP2207		Proportional Compact RH Mount Bracket <i>(NC 2203, 2204)</i>	○
ARP2208		Proportional Compact LH Mount Bracket <i>(NC 2203, 2204)</i>	○
ARP2213		Proportional Compact Centre Mount swing away to RH (Nucleus) <i>(NC 2203, 2204)</i>	€670
ARP2214		Proportional Compact Centre Mount swing away to LH (Nucleus) <i>(NC 2203, 2204)</i>	€670
ARX2215		Individual mounting option per request (Special) Requires additional information / specification <i>(NC 2203, 2204)</i>	POA
LINX HEAD CONTROL			
ARP2230		Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M3 <i>(NC 2203, 2204)</i>	€4,080

ARP2247		Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings (NC 2203, 2204)	€4,759
LiNX MICRO EXTREMITY CONTROL			
ARP2233		MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry) (NC 2203, 2204, 2229, 2230, 2247, 2257, 2240)	€3,837
LiNX MICRO EXTREMITY CONTROL BRACKET			
ARP2235		MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please also select MEC Bracketry)	○
ARP2236		MEC Tray Mount LH (ASL603, 100mm x 150mm)	○
ARP2213		Centre Mount, swing away to RH, Nucleus (NC 1404, 1405, 2310, 2311)	€670
ARP2214		Centre Mount, swing away to LH, Nucleus (NC 1404, 1405, 2310, 2311)	€670
ARP2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm) (NC 1404, 1405, 2310, 2311)	€670
ARP2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm) (NC 1404, 1405, 2310, 2311)	€670
LiNX HEAVY DUTY REMOTE			
ARP2240		Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LiNX Heavy Duty Bracketry) (NC 2203, 2204, 2229, 2230, 2247, 2257, 2233)	€2,433
LiNX HEAVY DUTY REMOTE BRACKET			
ARP2235		Tray Mount RH, (ASL603, 100mm x 150mm)	○
ARP2236		Tray Mount LH, (ASL603, 100mm x 150mm)	○
ARP2213		Centre Mount, swing away to RH, Nucleus (NC 1404, 1405, 2310, 2311)	€670
ARP2214		Centre Mount, swing away to LH, Nucleus (NC 1404, 1405, 2310, 2311)	€670
ARP2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm) (NC 1404, 1405, 2310, 2311)	€670
ARP2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm) (NC 1404, 1405, 2310, 2311)	€670
ARX2215		Individual mounting option per request (Special) Requires additional information / specification	POA
LiNX 10-WAY SWITCH			
ARP2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch) (Please select 10 way switch bracket below)	€466
LiNX 10-WAY SWITCH BRACKET			
ARP2217		10-Way Switch Bracket RH, Mounted for User	○
ARP2218		10-Way Switch Bracket LH, Mounted for User	○
ARP2227		10-Way Switch Bracket RH, Mounted for Attendant	€65
ARP2228		10-Way Switch Bracket LH, Mounted for Attendant	€65

LINX FOUR -SWITCH PROXIMITY ARRAY			
ARP2257		Four-Switch Proximity Array (ASL106) NEW! (NC 2203, 2204, 2229, 2230, 2092, 2323, 2324, 2322, 2325, 23)	€6,119
ARP2275		Centre Mount, swing-away to RH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm) (NC 1404, 1405, 2310, 2311)	€670
ARP2276		Centre Mount, swing-away to LH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm) (NC 1404, 1405, 2310, 2311)	€670
LINX WIRELESS ACCESSORIES FOR ATOM			
ARP2249		Wireless Triple Switch Receiver (ASL557-3) (Sub D plug) NEW!	€315
ARP2250		Wireless Mouse Emulator (ASL558) NEW!	€755
LINX PEDIATRIC COMPACT			
ARP2243		Proportional Paediatric Compact (ASK132), Includes input Module NEW! (NC 2271, 2203, 2204, 2255, 2256)	€2,211
ARP2235		Tray Mount RH (ASL603, 100mm x 150mm)	○
ARP2236		Tray Mount LH (ASL603, 100mm x 150mm)	○
ARP2213		Centre Mount, swing-away to RH, Nucleus (NC 1404, 1405, 2310, 2311)	€670
ARP2214		Centre Mount, swing-away to LH, Nucleus (NC 1404, 1405, 2310, 2311)	€670
ARP2241		Centre Mount, swing-away to RH, Nucleus (incl Tray (ASL603, 100mm x 150mm) (NC 1404, 1405, 2310, 2311)	€670
ARP2242		Centre Mount, swing-away to RH, Nucleus (incl Tray (ASL603, 100mm x 150mm) (NC 1404, 1405, 2310, 2311)	€670
LINX SIP N PUFF HEAD CONTROL			
ARP2255		Sip'n Puff Head Control (ASL109) NEW! includes Input Module and Lip Switch (ASL 301 LIP) (Combines simple sip and puff controls with head movements) (NC 2271, 2203, 2204, 2092, 2323, 2324, 2322, 2325, 2326, 23)	€5,731
LINX SIP N PUFF CONTROL & INPUT MODULE			
ARP2256		Sip'n Puff Control NEW! (including Sip'n Puff kit, Lip Switch ASL301LIP and Input module) (NC 2271, 2203, 2204, 2243)	€796
LINX INTERFACE OPTIONS			
ARP2271		Input Module DLX-IN500-A (sub D, Sip'n Puff) NEW! (NC 2243, 2255, 2256, 2203, 2204)	€737
LINX USB CHARGER PORT			
ARP2253		USB Charger Port DLX-USB02-A	€147
LINX USB CHARGER PORT BRACKET			
ARP2207		USB Charger Port RH Mount	○
ARP2208		USB Charger Port LH Mount	○
OUTPUT MODULE			
ARP2272		LiNX Output Module (OUT500)	€755
FURTHER LINX COMPONENTS / ACCESSORIES			
ARP2300		To be enclosed only	○

BUDDY BUTTONS			
ARP1685		Red Buddy Button	€82
ARP1686		Yellow Buddy Button	€82
ARP1687		Green Buddy Button	€82
ARP1690		Blue Buddy Button	€82
EGG SWITCH			
ARP2254		Egg Switch, Black	€82
PIKO BUTTONS			
ARP1711		Red Piko Button	€146
ARP1712		Yellow Piko Button	€146
ARP1713		Green Piko Button	€146
ARP1714		Orange Piko Button	€146
ARP1715		Blue Piko Button	€146
ARP1716		Black Piko Button	€146
TOGGLES / BUTTONS			
ARP2261		Stereo Toggle - direct control for Tilt only	€248
ARP2262		Stereo Button - direct control for Tilt only	€238
ARP2263		Four-Way Toggle - direct control for Tilt & Recline	€310
ARP2264		Four-Way Button - direct control for Tilt & Recline	€331
ARP2207		Toggle / Button RH Mount	○
ARP2208		Toggle / Button LH Mount	○
ARP2225		Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
ARP2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	€68
ARP2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	€68
ARP2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	€68
ARP2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	€68
REMOTE STOP SWITCH			
ARP2258		Remote Stop Switch (ASL504)	€970
SHROUD COLOUR *			
ARP2601		SpaceStation Silver	○
ARP2602		Island Blue	○
ARP2603		Grape JellyBean (Purple)	○
ARP2604		Monster Green	○
ARP2605		Tangy Orange	○
ARP2606		WhiteOut	○
ARP2607		Rockstar Red	○
ARP2608		Lights Out Black - Matt	○
ARP2609		Invacare Blue	○
RIM COLOUR *			
ARP2611		SpaceStation Silver	○
ARP2612		Island Blue	○
ARP2613		Grape JellyBean (Purple)	○
ARP2614		Monster Green	○
ARP2615		Tangy Orange	○
ARP2616		WhiteOut	○
ARP2617		Rockstar Red	○
ARP2618		Lights Out Black - Matt	○
ARP2619		Invacare Blue	○
ARP2620		Wet Black	○

ARP2621		Without Rim insert	○
FORKS *			
ARP2631		Single sided fork (New) includes Black castor	€325
ARP2632		Double sided fork	○
CENTRE WHEELS (DRIVING WHEELS) *			
ARP0704		355mm x 76mm (14" x 3") Black Tyres, Pneumatic	○
ARP0733		355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	◆
MOTOR TYPE *			
ARP1294		3 km/h (2mph)	○
ARP1282		6 km/h (4mph)	○
ARP1285		10 km/h (6mph) (includes light system and indicators - please select option ARM1367 or ARM1482)	◆
BATTERY *			
ARP1257		50 Ah	◆
PERFORMANCE ACCESSORIES			
ARP1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	◆
CHARGER *			
ARP1343		OptiCharge Charger 8A (200 - 250V)	◆
ARP1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€158
LIGHTS *			
ARP1437		Light system with indicators & LED Headlights (Only available with REM24SD or (NC 2203) G90A joystick controls)	€420
ARP1254		Without Light system	◆
TRANSPORTATION KIT *			
ARP1450		Preparation for DAHL Docking Station: Docking station for car implementation, please order docking station separately. (Adequate headrest is compulsory.	€541
ARP1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below) (NC 2053)	◆
ARP1401		6-Point-tie--down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19-2001) (Adequate headrest is compulsory) (NC 2052)	○
FIXATION BELT *			
ARP2113		Lap belt with auto lock (flexible mount in side rails)	◆
BODYPOINT PELVIC BELTS (includes mounting kit)			
ARP2701		Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	€106
ARP2702		Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	€106
ARP2703		2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	€98

ARP2704		2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	€98
ARP2705		2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)	€98
ARP2706		2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)	€98
TRAY			
ARP1404		Tray Swing Away RH <i>(NC 1207, 1707, 2014, 2018, 2209, 2207, 2213, 2214, 2241, 22</i>	€180
ARP1405		Tray Swing Away LH <i>(NC 1208, 1708, 2014, 2018, 2210, 2208, 2213, 2214, 2241, 22</i>	€180
ARP2310		Flip and fold away 1/2 lap tray LH <i>(NC 1208, 1708, 2014, 2018, 2210, 2208, 2213, 2214, 2241, 22</i>	€762
ARP2311		Flip and fold away 1/2 lap tray RH <i>(NC 1207, 1707, 2014, 2018, 2209, 2207, 2213, 2214, 2241, 22</i>	€762
TOTAL PRICE			

Non Conforming Options:

ARP2052	Centre of gravity Tilt (0 to 50°)
ARP2053	300mm seat lifter (increased STF height by 4mm) with centre of gravity tilt (0 to 45°)
ARP0339	Seat-to Floor Height 480mm
ARP2050	High Back, height 410mm (Mesh fabric)
ARP2051	High Back height 460mm (Mesh fabric)
ARP2320	High Back, height 505mm / 20" (Mesh fabric)
ARP2550	Matrx Elite Back, height 355mm / 14"
ARP2551	Matrx Elite Back, height 405mm / 15"
ARP2552	Matrx Elite Back, height 455mm / 18"
ARP2553	Matrx Elite Back, height 505mm / 20"
ARP2554	Matrx Elite HD Back, height 405mm / 16"
ARP2555	Matrx Elite HD Back, height 505mm / 20"
ARP2556	Matrx Elite TR Back, height 405mm / 16" (Mesh fabric)
ARP2557	Matrx Elite TR Back, height 455mm / 18" (Mesh fabric)
ARP2177	Matrx Elite TR Back height 505mm (Mesh fabric)
ARP2563	Matrx Elite TR HD Back, height 505mm / 20" (Mesh fabric)
ARP2558	Matrx Elite Deep Back, height 405mm / 16" (Mesh fabric)
ARP2559	Matrx Elite Deep Back, height 455mm / 18" (Mesh fabric)
ARP2178	Matrx Elite Deep Back, height 500mm (Mesh fabric)
ARP2174	Matrx Posture Back, height 405mm (Mesh fabric max dimension 455mm wide)
ARP2175	Matrx Posture Back height 505mm (Mesh fabric max dimension 455mm wide)
ARP2560	Matrx PB Deep Back, height 405mm / 16"
ARP2561	Matrx PB Deep Back, height 505mm / 20"
ARP2562	Matrx Elite Deep HD Back, height 505mm / 20"
ARP2070	560mm (22") measured from seat pan to top of back shell
ARP2071	585mm (23") measured from seat pan to top of back shell
ARP2073	635mm (25") measured from seat pan to top of back shell
ARP2074	660mm (26") measured from seat pan to top of back shell

ARP2095	Cantilever Flip-Back Arm, width, height and angle adjustable (230mm - 330mm) (9" - 13)
ARP2096	Cantilever Flip-Back Arm, width, height and angle adjustable (320mm - 405mm) (12.5" - 16")
ARP2308	1" Offset Armrest Spacer (pair), for two-post recline only
ARP2309	1/2" Offset Armrest Spacer (pair), for two recline only
ARP2014	Ergonomic arm though 30mm thick
ARP2018	Ergonomic arm though 30mm thick
ARP2347	Elbow Block XS (110mm wide x 83mm high)
ARP2348	Elbow Block S (135mm wide x 90mm high)
ARP2349	Elbow Block M (145mm wide x 95mm high)
ARP2350	Elbow Block L (165mm wide x 110mm high)
ARP2351	Elbow Block XS (110mm wide x 83mm high)
ARP2352	Elbow Block S (135mm wide x 90mm high)
ARP2353	Elbow Block M (145mm wide x 95mm high)
ARP2354	Elbow Block L (165mm wide x 110mm high)
ARP2092	Onyx Auto Style headrest
ARP2323	Elan Standard Pad 150mm / 6" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW)
ARP2324	Elan Standard Pad 255mm / 10" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW)
ARP2322	Elan Standard Pad 355mm / 14" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW)
ARP2325	Elan Occipital Pad 230mm x 180mm / 9" x 7" w/standard cover (EOP23) complete with Elan Headrest Hardware (MEHW)
ARP2326	Elan Occipital Pad 305mm x 200mm / 12" x 8" w/standard cover complete with Elan Headrest Hardware (MEHW)
ARP2327	Elan 4 Point Headrest Pad 280mm x 255mm / 11" x 10" w/standard cover complete with Elan Headrest Hardware (MEHW)
ARP1495	Head Control or equivalent Headrest, provided independently from this order (ACCORDING TO THE REQUIREMENTS OF ISO 7176-19, A HEADREST IS COMPULSORY FOR USING THE WHEELCHAIR AS VEHICLE SEAT)
ARP0586	Hanger, centre mounted, angle settable (0° - 90°): knee to heel length 310mm - 410mm
ARP0537	Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARP0538	ADM MANUAL: Angle adjustable with quick release lever 0° - 80°, length 290mm - 460mm
ARP0539	ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm
ARP0572	Pivot Plus Manual Elevating Legrest
ARP0577	Pivot Plus Power Elevating Legrest
ARP0535	Without Hanger and footplates
ARP0547	Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARP0548	ADM MANUAL: Angle adjustable with quick release lever 0° - 80°, length 290mm - 460mm
ARP0549	ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm

ARP0579	Pivot Plus Manual Elevating Legrest
ARP0584	Pivot Plus Power Elevating Legrest
ARP0545	Without Hanger and footplates
ARP0505	Fixed Footplate
ARP0509	Angle Adjustable Footplate
ARP0511	Angle and depth adjustable Footplate
ARP0502	Fixed Footplate
ARP0506	Angle Adjustable Footplate
ARP0508	Angle and depth adjustable Footplate
ARP2203	REM211 Drive & Actuators
ARP2204	REM216 Drive, Actuators & Lights
ARP1207	Standard Bracket RH
ARP1208	Standard Bracket LH
ARP1707	Swing Away + height adjustable Bracket RH High end
ARP1708	Swing Away + height adjustable Bracket LH High end
ARP2075	Swing Away + height / width adjust mount on RH for user, Quad link design
ARP2076	Swing Away + height / width adjust mount on LH for user, Quad link design
ARP2227	Centre Mount swing-away Bracket to RH (Nucleus)
ARP2228	Centre Mount swing-away Bracket to LH (Nucleus)
ARP2209	REM500 Swing away RH multi-adjustable
ARP2210	REM500 Swing away LH multi-adjustable
ARP2213	REM500 Centre mount swing away to RH Nucleus
ARP2214	REM500 Centre mount swing away to LH Nucleus
ARP2075	Swing Away + height / width adjustable Bracket Mount on RH, Quad link design
ARP2076	Swing Away + height / width adjustable Bracket Mount on RH Quad link design
ARP2229	Proportional Chin Joystick (ASL138) (Requires Primary Remote REM400 or REM500) (Please select chin control below)
ARP2219	Compact Dual Control (DLX-ACU200) (Please select Dual control bracket below)
ARP1207	Dual Control Standard Bracket RH
ARP1208	Dual Control Standard Bracket LH
ARP2245	IDC Dual control (Please select AQE1221 OR AQE1222 for bracket option below)
ARP1221	IDC Dual control Detachable Bracket RH
ARP1222	IDC Dual control Detachable Bracket LH
ARP2231	Proportional Compact Remote DLX-CR400 (Please select Proportional Bracket)
ARP2232	Proportional Compact Remote Low Force DLX-CR400LF (Please select Proportional Bracketry)
ARP2207	Proportional Compact RH Mount Bracket
ARP2208	Proportional Compact LH Mount Bracket
ARP2213	Proportional Compact Centre Mount swing away to RH (Nucleus)
ARP2214	Proportional Compact Centre Mount swing away to LH (Nucleus)
ARX2215	Individual mounting option per request (Special) Requires additional information / specification

ARP2230	Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M3
ARP2247	Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings
ARP2233	MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry)
ARP2213	Centre Mount, swing away to RH, Nucleus
ARP2214	Centre Mount, swing away to LH, Nucleus
ARP2241	Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARP2242	Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARP2240	Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LiNX Heavy Duty Bracketry)
ARP2213	Centre Mount, swing away to RH, Nucleus
ARP2214	Centre Mount, swing away to LH, Nucleus
ARP2241	Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARP2242	Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARX2215	Individual mounting option per request (Special) Requires additional information / specification
ARP2227	10-Way Switch Bracket RH, Mounted for Attendant
ARP2228	10-Way Switch Bracket LH, Mounted for Attendant
ARP2257	Four-Switch Proximity Array (ASL106) NEW!
ARP2275	Centre Mount, swing-away to RH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm)
ARP2276	Centre Mount, swing-away to LH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm)
ARP2243	Proportional Paediatric Compact (ASK132), Includes input Module NEW!
ARP2213	Centre Mount, swing-away to RH, Nucleus
ARP2214	Centre Mount, swing-away to LH, Nucleus
ARP2241	Centre Mount, swing-away to RH, Nucleus (incl Tray (ASL603, 100mm x 150mm)
ARP2242	Centre Mount, swing-away to RH, Nucleus (incl Tray (ASL603, 100mm x 150mm)
ARP2255	Sip'n Puff Head Control (ASL109) NEW! includes Input Module and Lip Switch (ASL 301 LIP) (Combines simple sip and puff controls with head movements)
ARP2256	Sip'n Puff Control NEW! (including Sip'n Puff kit, Lip Switch ASL301LIP and Input module)
ARP2271	Input Module DLX-IN500-A (sub D, Sip'n Puff) NEW!
ARP2207	USB Charger Port RH Mount
ARP2208	USB Charger Port LH Mount
ARP2261	Stereo Toggle - direct control for Tilt only
ARP2262	Stereo Button - direct control for Tilt only
ARP2263	Four-Way Toggle - direct control for Tilt & Recline
ARP2264	Four-Way Button - direct control for Tilt & Recline

ARP2207	Toggle / Button RH Mount
ARP2208	Toggle / Button LH Mount
ARP1437	Light system with indicators & LED Headlights (Only available with REM24SD or G90A joystick controls)
ARP1404	Tray Swing Away RH
ARP1405	Tray Swing Away LH
ARP2310	Flip and fold away 1/2 lap tray LH
ARP2311	Flip and fold away 1/2 lap tray RH



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions - Please supply the following information:

Name :

Email :

Phone :

Fax :

Description of special request (please include relevant option code or part number):