



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 Ultra

Standard Features

- Max User Weight 136kg (21.4 Stone)
- Crash tested ISO7176- 19
- Ultra Low Maxx Seating System
- Seat widths of 405mm - 610mm (16" - 24"), adjustable between 405mm - 505mm (16" - 20") and 505mm - 610mm (20" - 24")
- Seat depths of 380mm - 580mm (15" - 23") adjustable between 380mm - 480mm (15" - 19") and 480mm - 580mm (19 - 23")
- Backrest heights of 560mm - 660mm in 25mm increments
- Cantilever flip-back arm, width, height and angle adjustable with desk armpad with startex covered visco foam
- Elan multi positioning headrest
- Ultra Low Maxx CG power tilt (0° to 50°) seat heights of 420mm - 480mm
- Powered Recline 90° to 170° with extended shear reduction (ESR)
- Vari F hangers with fixed footplates
- LiNX DLX-REM211 with standard bracket
- Selection of 5 shroud and rim colours
- Centre (driving wheels) 355mm x 76mm (14" x 3") Black tyres
- 10km/h (6mph) 4 pole motors and 50 Ah batteries
- G-Trac
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: €18,703



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code **Description**

BASE MODEL*

ARO0080	☐	TDX SP2 Ultra	◆
ARZ1300	☐	Special TDX SP2 Ultra (Extended Lead Time). Please provide details in comments section.	POA

SEATING TYPE*

ARO2020	☐	Ultra Low Maxx Powered Positioning System	€18,703
---------	--------------------------	---	---------

SEAT ADJUSTMENTS*				
ARO2052		Centre of gravity Tilt (0 to 50°)		◆
ARO2053		300mm seat lifter (increased STF height by 4mm) with centre of gravity tilt (0 to 45°)	(NC 0339)	€2,787
SEAT TO GROUND HEIGHT*				
ARO0337		440mm		○
ARO0331		460mm		○
ARO0339		480mm	(NC 2053)	○
ANTERIOR AND POSTERIOR TILT*				
ARO2057		Tilt Standard Mount (No anterior or posterior tilt)		◆
ARO2054		5° Fixed Anterior Assist Mount (will reduce max tilt angle by 5°)		€481
ARO2055		10° Fixed Anterior Assist Mount (will reduce max tilt angle by 10°)	(NC 2035)	€481
ARO2056		5° Pre-Tilt (will increase max tilt angle by 5°)		€163
BACKREST RECLINER*				
ARO2060		Powered Recliner (90° - 170°) with Extended Shear Reduction (ESR), Standard Mount		◆
ARO2063		Powered Recliner 8° Recline mount 82° to 162° with Extended Shear Reduction (ESR)		€38
ARO2066		Powered Recline, 30° Recline mount, 60° to 140° with Extended Shear Reduction (ESR)		€38
SEAT FRAME WIDTH RANGE*				
ARO2021		405mm - 505mm (increments of 25mm) (16" - 20")	(NC 2030, 2031, 2032, 2033)	○
ARZ2003		480mm - 555mm (Increments of 25mm) 19" - 22") - Additional build time of 15+days	(NC 2025, 2026, 2027, 2032, 2033)	○
ARZ2022		505mm - 610mm (increments of 25mm) (20" - 24") - Additional build time of 15+days	(NC 2025, 2026, 2027, 2174, 2175)	○
SEAT FRAME DEPTH RANGE*				
ARO2024		480mm - 580mm (increments of 25mm) (19" - 23")	(NC 2035, 2036, 2037)	○
ARO2008		405mm - 505mm (increments of 25mm) (16" - 20")	(NC 2040, 2041, 2042)	○
SEAT WIDTH*				
ARO2025		405mm (16")	(NC 2003, 2022, 2041, 2042, 2701, 2703, 2705)	○
ARO2026		430mm (17")	(NC 2003, 2022, 2041, 2042, 2701, 2703, 2705)	○
ARO2027		455mm (18")	(NC 2003, 2022, 2701, 2703, 2705)	○
ARO2028		480mm (19")	(NC 2022, 2174, 2175, 2702, 2704, 2706)	○
ARO2029		505mm (20")	(NC 2174, 2175, 2702, 2704, 2706)	○
ARO2030		530mm (21")	(NC 2021, 2035, 2036, 2174, 2175, 2702, 2704, 2706)	○
ARO2031		555mm (22")	(NC 2021, 2035, 2036, 2174, 2175, 2702, 2704, 2706)	○
ARO2032		580mm (23")	(NC 2021, 2035, 2036, 2174, 2175, 2147, 2702, 2704, 2706)	○
ARO2033		640mm (24")	(NC 2021, 2035, 2036, 2174, 2175, 1404, 1405, 2147, 2702, 2704, 2706)	○
SEAT DEPTH*				
ARO2035		405mm (16")	(NC 2024, 2030, 2031, 2032, 2033, 2055, 2701, 2702, 2703, 2704, 2705, 2706)	○
ARO2036		430mm (17")	(NC 2024, 2030, 2031, 2032, 2033)	○
ARO2037		455mm (18")	(NC 2024)	○
ARO2038		480mm (19")		○
ARO2039		505mm (20")		○
ARO2040		530mm (21")		○
ARO2041		555mm (22")	(NC 2025, 2026)	○
ARO2042		580mm (23")	(NC 2025, 2026, 2146)	○

SEAT CUSHION*			
ARO2144		Matrx Vi Cushion (Mesh fabric) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983)	€348
ARO2146		Matrx Flo-tech Contour Visco (NC 2042, 2030, 2031, 2033, 1985)	€174
ARO2147		Matrx Libra 90mm thick	€532
ARO2140		Without seat cushion (NC 2033, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	◆
FIRST POSITIONING ACCESORIES FOR MATRX LIBRA AND FLO-TECH LITE			
ARO1980		Tapered Wedges (Pair)	€28
ARO1981		Abductor	€58
ARO1982		Lateral Wedge	€46
ARO1983		Anterior Seat Wedges	€43
POSITIONING ACCESSORIES FOR MATRX VI			
ARO1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)	€71
BACKREST SUPPORT*			
ARO2050		High Back, height 410mm (Mesh fabric) (NC 2073, 2074)	◆
ARO2051		High Back height 460mm (Mesh fabric) (NC 2070, 2071)	○
ARO2320		High Back, height 505mm / 20" (Mesh fabric)	○
ARO2550		Matrx Elite E2 Back, height 355mm / 14" (NC 2028, 2029, 2030, 2032, 2033)	€572
ARO2551		Matrx Elite E2 Back, height 405mm / 15"	€572
ARO2552		Matrx Elite E2 Back, height 455mm / 18"	€572
ARO2553		Matrx Elite E2 Back, height 505mm / 20"	€572
ARO2554		Matrx Elite HD Back, height 405mm / 16"	€741
ARO2555		Matrx Elite HD Back, height 505mm / 20"	€741
ARO2556		Matrx Elite TR Back, height 405mm / 16" (Mesh fabric)	€572
ARO2557		Matrx Elite TR Back, height 455mm / 18"(Mesh fabric)	€572
ARO2177		Matrx Elite TR Back height 505mm (Mesh fabric)	€572
ARO2563		Matrix Elite TR HD Back, height 505mm / 20" (Mesh fabric)	€741
ARO2558		Matrx Elite Deep Back, height 405mm / 16" (Mesh fabric)	€741
ARO2559		Matrx Elite Deep Back, height 455mm / 18" (Mesh fabric)	€741
ARO2178		Matrx Elite Deep Back, height 500mm (Mesh fabric)	€741
ARO2174		Matrx Posture Back, height 405mm (Mesh fabric max dimension 455mm wide)	€572
ARO2175		Matrx Posture Back height 505mm (Mesh fabric max dimension 455mm wide)	€572
ARO2560		Matrx PB Deep Back, height 405mm / 16"	€741
ARO2561		Matrx PB Deep Back, height 505mm / 20"	€741
ARO2562		Matrx Elite Deep HD Back, height 505mm / 20"	€741
ARO2168		Matrx MAC Back - Foam Cushion [Mesh fabric]	○
ARO2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]	○

BACKREST HEIGHT*				
ARO2506		505mm (20") measured from seat pan to top of back shell	(NC 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562)	○
ARO2501		535mm (21") measured from seat pan to top of back shell	(NC 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562)	○
ARO2070		560mm (22") measured from seat pan to top of back shell	(NC 2050, 2551, 2553, 2554, 2555, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563)	○
ARO2071		585mm (23") measured from seat pan to top of back shell	(NC 2050, 2551, 2553, 2554, 2555, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563)	○
ARO2072		610mm (24") measured from seat pan to top of back shell	(NC 2050, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	○
ARO2073		635mm (25") measured from seat pan to top of back shell	(NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2558, 2559, 2174, 2560)	○
ARO2074		660mm (26") measured from seat pan to top of back shell	(NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	○
ARO2502		685mm (27") measured from seat pan to top of back shell	(NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	◆
MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)				
ARO2380		405 mm (16 inch) MAC BH		○
ARO2381		455 mm (18 inch) MAC BH		○
ARO2382		505 mm (20 inch) MAC BH		○
ARO2383		560 mm (22 inch) MAC BH		○
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)				
ARO2385		90 mm (3.5 inch) MAC Contour Depth LH		○
ARO2386		150 mm (6 inch) MAC Contour Depth LH		○
ARO2387		230 mm (9 inch) MAC Contour Depth LH		○
MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)				
ARO2389		90 mm (3.5 inch) MAC Contour Depth RH		○
ARO2390		150 mm (6 inch) MAC Contour Depth RH		○
ARO2391		230 mm (9 inch) MAC Contour Depth RH		○
MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)				
ARO2393		Backrest -25 mm (-1 inch) to seat width		○
ARO2394		Backrest matches seat width		○
ARO2395		Backrest +25 mm (+1 inch) to seat width		○
ARMREST*				
ARO2095		Cantilever Flip-Back Arm, width, height and angle adjustable (240mm - 330mm)		○
ARO2096		Cantilever Flip-Back Arm, width, height and angle adjustable (330mm - 405mm)		○
ARO2097		Reclining Arm, Flip-Back, adjustable height (240mm - 330mm)		○
ARO0690		Without Armrest		○
ARMPAD LH*				
ARO2012		Desk arm pad with startex cover visco foam (255mm) 40m thick		○
ARO2013		Full length armpad with startex cover foam (355mm) 40mm thick		○
ARO2014		Ergonomic arm though 30mm thick	(NC 1404, 1405, 2347, 2348, 2349, 2350)	€223
ARO2015		Without LH arm pad		◆

ARMPAD RH*			
ARO2016		Desk arm pad with startex cover visco foam (255mm) 40mm thick	○
ARO2017		Full length armpad with startex cover foam (355mm) 40mm thick	○
ARO2018		Ergonomic arm trough 30mm thick (NC 1404, 1405, 2351, 2352, 2353, 2354)	€223
ARO2019		Without RH arm pad	€148
ARMPAD ACCESSORIES & MOUNTING OPTIONS			
ARO2080		Flat Hand Pad Left (in combination with arm trough only): to be mounted on the opposite side of remote	€60
ARO2083		Flat Hand Pad Right (in combination with arm trough only): to be mounted on the opposite side of remote	€60
MULTI-AXIS UPPER EXTREMITY SUPPORTS			
ARO2082		Multi Axis Upper Extremity Support (MACES) Left, raises arm pad by 43mm	€362
ARO2085		Multi Axis Upper Extremity Support (MACES) Right, raises arm pad by 43mm	€362
ARMREST HEIGHT SETTING LH*			
O0001		Factory Setting	○
O0002		Customised 230mm	○
O0003		Customised 240mm	○
O0004		Customised 250mm	○
O0005		Customised 260mm	○
O0006		Customised 270mm	○
O0007		Customised 280mm	○
O0008		Customised 290mm	○
O0009		Customised 300mm	○
O0010		Customised 310mm	○
O0011		Customised 320mm	○
O0012		Customised 330mm	○
O0013		Customised 340mm	○
O0014		Customised 350mm	○
O0015		Customised 360mm	○
O0016		Customised 370mm	○
O0017		Customised 380mm	○
O0018		Customised 390mm	○
O0019		Customised 405mm	○

ARMREST HEIGHT SETTING RH*				
O0020		Factory Setting		○
O0021		Customised 230mm		○
O0022		Customised 240mm		○
O0023		Customised 250mm		○
O0024		Customised 260mm		○
O0025		Customised 270mm		○
O0026		Customised 280mm		○
O0027		Customised 290mm		○
O0028		Customised 300mm		○
O0029		Customised 310mm		○
O0030		Customised 320mm		○
O0031		Customised 330mm		○
O0032		Customised 340mm		○
O0033		Customised 350mm		○
O0034		Customised 360mm		○
O0035		Customised 370mm		○
O0036		Customised 380mm		○
O0037		Customised 390mm		○
O0038		Customised 405mm		○
ELBOW BLOCK LH PADS				
ARO2347		Elbow Block XS (110mm wide x 83mm high)	(NC 2014)	€521
ARO2348		Elbow Block S (135mm wide x 90mm high)	(NC 2014)	€521
ARO2349		Elbow Block M (145mm wide x 95mm high)	(NC 2014)	€521
ARO2350		Elbow Block L (165mm wide x 110mm high)	(NC 2014)	€521
ELBOW BLOCK LH MOUNTING HARDWARE				
ARO2361		Vertical 3"/ Horizontal 5"		○
ARO2362		Vertical 3"/ Horizontal 7"		○
ARO2363		Vertical 5"/ Horizontal 5"		○
ARO2364		Vertical 5"/ Horizontal 7"		○
ELBOW BLOCKS RH PADS				
ARO2351		Elbow Block XS (110mm wide x 83mm high)	(NC 2018)	€521
ARO2352		Elbow Block S (135mm wide x 90mm high)	(NC 2018)	€521
ARO2353		Elbow Block M (145mm wide x 95mm high)	(NC 2018)	€521
ARO2354		Elbow Block L (165mm wide x 110mm high)	(NC 2018)	€521
ELBOW BLOCK RH MOUNTING HARDWARE				
ARO2371		Vertical 3"/ Horizontal 5"		○
ARO2372		Vertical 3"/ Horizontal 7"		○
ARO2373		Vertical 5"/ Horizontal 5"		○
ARO2374		Vertical 5"/ Horizontal 7"		○

HIP SUPPORTS LH (Not available with 410mm seat depth and Manual Centre legrests)			
ARO2184		Hip Support LH, Quick Release, pad size 140 x 100 mm (NC 2701, 2702, 2703, 2704, 2705, 2706)	€162
ARO2114		Hip Support LH, Quick Release, pad size 180 x 150 mm (NC 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARO2186		Hip Support LH, Quick Release, pad size 230 x 100 mm (NC 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARO2188		Hip Support LH, Quick Release, pad size 330 x 100 mm (NC 2701, 2702, 2703, 2704, 2705, 2706)	€272
HIP SUPPORTS RH (Not available with 410mm seat depth and Manual Centre legrests)			
ARO2185		Hip Support RH, Quick release, pad size 140 x 100 mm (NC 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARO2115		Hip Support RH, Quick Release, pad size 180 x 150 mm (NC 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARO2187		Hip Support RH, Quick release, pad size 230 x 100 mm	€272
ARO2189		Hip Support RH, Quick Release, pad size 330 x 100 mm	€272
HEADREST*			
ARO2092		Onyx Auto Style headrest (NC 2255, 2230, 2247)	€207
ARO2324		Elan Standard Pad 255mm / 10" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	€334
ARO2325		Elan Occipital Pad 230mm x 180mm / 9" x 7" w/standard cover (EOP23) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	€334
ARO1495		Head Control or equivalent Headrest, provided independently from this order (ACCORDING TO THE REQUIREMENTS OF ISO 7176-19, A HEADREST IS COMPULSORY FOR USING THE WHEELCHAIR AS VEHICLE SEAT)	○
LATERAL SUPPORTS LH			
ARO2343		110mm wide x 83mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€402
ARO2345		135mm wide x 90mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€402
ARO2043		145mm wide x 95mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€362
ARO2045		165mm wide x 110mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€428
LATERAL SUPPORTS RH			
ARO2344		110mm wide x 83mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€402
ARO2346		135mm wide x 90mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€402
ARO2044		145mm wide x 95mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€362
ARO2046		165mm wide x 110mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€427
LATERAL SUPPORT ACCESSORIES			
ARO2049		Additional telescopic link for Lateral support, 25mm, pair (NC 2558, 2559, 2178, 2560, 2561, 2562)	€170
CENTRE MOUNTED LEGREST			
ARO0586		Hanger, centre mounted, angle settable (0° - 90°): knee to heel length 310 - 410mm	€357
ARO0589		Hangers, centre mounted, length compensating, powered, LNX (0° - 90°): KNEE TO HEL LENGTH 340 - 410mm (NC 2261, 2262, 2263, 2264)	€1,985
FOOTPLATES / FOOT PLATFORM			
ARO2365		Angle & width adjustable footplates Small (127mm width x 190mm depth / 5" x 7.5")	€106
ARO0559		Angle & width adjustable footplate (240mm depth)	○
ARO0516		Angle & width adjustable footplates (290mm depth)	€274
ARO0517		Angle adjustable foot platform (240mm depth, 305mm width)	€274

CENTRE MOUNT LEGREST HEIGHT SETTINGS 0586 & 0589			
O0039		Factory Setting	○
O0040		Customised 310mm	○
O0041		Customised 325mm	○
O0042		Customised 335mm	○
O0043		Customised 355mm	○
O0044		Customised 355mm	○
O0045		Customised 364mm	○
O0046		Customised 375mm	○
O0047		Customised 385mm	○
O0048		Customised 410mm	○
CENTRE MOUNT W TELESCOPING FOOT PLATFORM HEIGHT SETTINGS			
O0049		Factory Setting	○
O0050		Customised 380mm	○
O0051		Customised 393mm	○
O0052		Customised 406mm	○
O0053		Customised 419mm	◆
SIDE MOUNTED HANGERS LH			
ARO0536		Vari F: Angle adjustable with tools 0° - 70°, length 290mm - 460mm (NC 0505)	€134
ARO0537		Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm (NC 0505)	€128
ARO0539		ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm (NC 0505, 2261, 2262, 2263, 2264)	€1,225
ARO0535		Without Hanger and footplates (NC 0505, 0509, 0511)	○
VARI HANGER LH HEIGHT SETTING			
O0054		Factory Setting	○
O0055		Customised 270mm	○
O0056		Customised 280mm	○
O0057		Customised 290mm	○
O0058		Customised 300mm	○
O0059		Customised 310mm	○
O0060		Customised 320mm	○
O0061		Customised 330mm	○
O0062		Customised 340mm	○
O0063		Customised 350mm	○
O0064		Customised 360mm	○
O0065		Customised 370mm	○
O0066		Customised 380mm	○
O0067		Customised 390mm	○
O0068		Customised 400mm	○
O0069		Customised 410mm	○
O0070		Customised 420mm	○
O0071		Customised 430mm	○

SIDE MOUNTED HANGERS RH				
ARO0546		Vari F: Angle adjustable with tools 0° - 70°, length 290mm - 460mm	(NC 0502)	€134
ARO0547		Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	(NC 0502)	€128
ARO0549		ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm	(NC 0502, 2261, 2262, 2263, 2264)	€1,225
ARO0545		Without Hanger and footplates	(NC 0502, 0506, 0508)	○
VARI HANGER RH HEIGHT SETTING				
O0080		Factory Setting		○
O0081		Customised 270mm		○
O0082		Customised 280mm		○
O0083		Customised 290mm		○
O0084		Customised 300mm		○
O0085		Customised 310mm		○
O0086		Customised 320mm		○
O0087		Customised 330mm		○
O0088		Customised 340mm		○
O0089		Customised 350mm		○
O0090		Customised 360mm		○
O0091		Customised 370mm		○
O0092		Customised 380mm		○
O0093		Customised 390mm		○
O0094		Customised 400mm		○
O0095		Customised 410mm		○
O0096		Customised 420mm		○
O0097		Customised 430mm		○
INDIVIDUAL FOOTPLATE TYPE LH				
ARO0505		Fixed Footplate	(NC 0535, 0537, 0539)	€75
ARO0509		Angle Adjustable Footplate	(NC 0535)	○
ARO0511		Angle and depth adjustable Footplate	(NC 0535)	◆
INDIVIDUAL FOOTPLATE TYPE RH				
ARO0502		Fixed Footplate	(NC 0545, 0547, 0549)	€75
ARO0506		Angle Adjustable Footplate	(NC 0545)	○
ARO0508		Angle and depth adjustable Footplate	(NC 0545)	○
CONTROLS*				
ARO2203		DLX-REM211 (Drive & Actuators)	(NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2209, 2210, 2213, 2214, 2075, 2076)	◆
ARO2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2209, 2210, 2213, 2214, 2075, 2076)	€495
ARO2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	(NC 2209, 2210, 2213, 2215, 2075, 2076)	€909
ARO2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 1207, 1208, 1707, 1708, 1709, 1710, 2227, 2228, 2075, 2076)	○
REM400 TOGGLE OPTION				
ARO2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2203, 2204)	€3,106

CONTROL BRACKETRY*				
ARO1207		Standard Bracket RH	(NC 1404)	○
ARO1208		Standard Bracket LH	(NC 1405)	○
ARO1707		Swing Away RH Bracket	(NC 1404)	€134
ARO1708		Swing Away LH Bracket	(NC 1405)	€134
ARO1709		Matrx Resolve Swing in/out Bracket RH		€281
ARO1710		Matrx Resolve Swing in/out Bracket LH		€281
ARO2209		REM500 Swing Away RH Multi-Adjustable Bracket	(NC 1404)	○
ARO2210		REM500 Swing Away LH Multi-Adjustable Bracket	(NC 1405)	○
ARO2075		Swing Away + height / width adjust. Mount on RH for user, Quad link design	(NC 1415)	€228
ARO2076		Swing Away + height / width adjust. Mount on LH for user, Quad link design	(NC 1339)	€228
ARO2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 1404, 1405)	€739
ARO2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 1404, 1405)	€739
ARO2227		Bracket RH mounted for Attendant	(NC 1221, 2206, 2219)	€84
ARO2228		Bracket LH mounted for Attendant	(NC 1222, 2206, 2219)	€84
ARZ2239		Individual mounting option per request (Special). Requires additional information/specification.		POA
DUAL CONTROLS				
ARO2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1221, 1222, 2227, 2228, 2245)	€401
ARO2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2219, 2227, 2228)	€748
UNZ9259		mo-vis Scoot Control (Only available via Special)	(NC 0080, 1221, 1222)	€3,020
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
ARO1207		Dual Control Standard Bracket RH	(NC 1221, 1222, 2227, 2228, 2245)	○
ARO1208		Dual Control Standard Bracket LH	(NC 1221, 1222, 2227, 2228, 2245)	○
ARO1221		IDC Dual Control Detachable Bracket RH	(NC 2219, 2227, 2228)	○
ARO1222		IDC Dual Control Detachable Bracket LH	(NC 2219, 2227, 2228)	○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification.		POA
COMPACT REMOTES				
ARO2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204)	€845
ARO2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204)	€1,199
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	€3,430
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	€3,560
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	€3,060
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	€3,285
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	€1,820
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	€2,630
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	€4,160
ARO2233		ASL MEC (Proportional MEC Joystick ASL130)	(NC 2203, 2204, 2230, 2247, 2257, 2240)	€4,230
ARO2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133)	(NC 2203, 2204, 2230, 2247, 2257, 2233)	€2,683
ARO2243		ASL Paediatric (Proportional Joystick ASL132, incl. DLX-IN500)	(NC 2271, 2203, 2204, 2255, 2256)	€2,438

COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
ARO2207		RH Mount	(NC 2203, 2204)	○
ARO2208		LH Mount	(NC 2203, 2204)	○
ARO2213		Centre Mount Swing Away to RH (Nucleus)	(NC 2203, 2204)	€739
ARO2214		Centre Mount Swing Away to LH (Nucleus)	(NC 2203, 2204)	€739
ARO2235		Tray Mount RH (ASL603, 100mm x 150mm)		○
ARO2236		Tray Mount LH (ASL603, 100mm x 150mm)		○
ARO2213		Centre Mount Swing Away to RH (Nucleus)	(NC 1404, 1405)	€739
ARO2214		Centre Mount Swing Away to LH (Nucleus)	(NC 1404, 1405)	€739
ARO2241		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405)	€739
ARO2242		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405)	€739
ARZ2215		Individual mounting option per request (Special). Requires additional information/specification.	(NC 2203, 2204)	POA
CHIN CONTROLS				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204)	€3,430
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	€3,560
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204)	€3,060
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	€3,285
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	€1,820
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	€2,630
CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)				
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)		€1,255
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)		€1,820
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)		€3,005
HEAD CONTROLS				
ARO2230		ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2203, 2204)	€4,498
ARO2247		ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2203, 2204)	€5,247
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2092, 2324, 2325, 2257, 2240, 2175, 0080)	€9,135
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2092, 2324, 2325, 2257, 2240, 2175, 0080)	€5,435
FOOT CONTROL				
UNZ9264		mo-vis Foot Control LiNX (Only available via Special)	(NC 2203, 2204)	€3,800
UNZ9265		mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	€3,615
FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)				
UNZ9266		Individual mounting option per request (Special). Requires additional information/specification.		POA
TRAY CONTROL				
ARO2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2230, 2092, 2324, 2325, 2257, 2240)	€6,746
TRAY CONTROL BRACKETRY (Only applicable if Tray Control is selected)				
ARO2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1404, 1405)	€739
ARO2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1404, 1405)	€739

SIP N PUFF HEAD CONTROL			
ARO2255	☐	Sip'n Puff Head Control (ASL109) includes Input Module and Lip Switch (ASL 301 LIP) (Combines simple sip and puff controls with head movements)	(NC 2271, 2203, 2204, 2092, 2324, 2325, 2243) €6,319
SIP N PUFF CONTROL & INPUT MODULE			
ARO2256	☐	Sip'n Puff Control (including Sip'n Puff kit, Lip Switch ASL301LIP and Input module)	(NC 2271, 2203, 2204, 2243) €878
WIRELESS ACCESSORIES FOR ASL CONTROLS			
ARO2249	☐	Wireless Triple Switch Receiver (ASL557-3)	€348
ARO2250	☐	Wireless Mouse Emulator (ASL558)	€833
LiNX INTERFACE OPTIONS			
ARO2271	☐	Input Module DLX-IN500-A (sub D, Sip'n Puff)	(NC 2243, 2255, 2256, 2203, 2204) €813
ARO2272	☐	LiNX Output Module (OUT500)	€833
LiNX 10-WAY SWITCH			
ARO2251	☐	10-Way Switch DLX-SW10-A for direct Powered Seating Access	€513
LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)			
ARO2217	☐	Bracket RH, mounted for user	☐
ARO2218	☐	Bracket LH, mounted for user	☐
ARO2227	☐	For attendant, mounted on back shroud	€71
ARO2228	☐	Without mounting	€71
LiNX USB CHARGER PORT			
ARO2253	☐	USB Charger Port DLX-USB02-A	€162
LiNX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)			
ARO2207	☐	RH Mount	☐
ARO2208	☐	LH Mount	☐
REMOTE STOP SWITCH			
ARO2258	☐	Remote Stop Switch (ASL504)	€1,070
ARO2258	☐	Remote Stop Switch (ASL504)	€1,070
FURTHER LiNX COMPONENTS / ACCESSORIES			
ARO2302	☐	Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
ARO2300	☐	To be enclosed only	☐
BUDDY BUTTONS			
ARO1685	☐	Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€90
ARO1686	☐	Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€90
ARO1687	☐	Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€90
ARO1690	☐	Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€90
EGG SWITCH			
ARO2254	☐	Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	€90

PIKO BUTTONS			
ARO1711		Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
ARO1712		Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
ARO1713		Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
ARO1714		Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
ARO1715		Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
ARO1716		Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
TOGGLES / BUTTONS			
ARO2261		Stereo Toggle - direct control for Tilt only	€263
ARO2262		Stereo Button - direct control for Tilt only	€263
ARO2263		Four-Way Toggle - direct control for Tilt & Recline	€342
ARO2264		Four-Way Button - direct control for Tilt & Recline	€365
ARO2207		Toggle / Button RH Mount	○
ARO2208		Toggle / Button LH Mount	○
ARO2225		Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
ARO2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	€75
ARO2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	€75
ARO2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	€75
ARO2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	€75
SHROUD COLOUR*			
ARO2640		Obsidian Black Frame Shrouds	○
ARO2641		Ruby Red Frame Shrouds	○
ARO2642		Purple Passion Frame Shrouds	○
ARO2643		Peacock Blue Frame Shrouds	○
ARO2644		Arctic Blue Frame Shrouds	○
RIM COLOUR*			
ARO2650		Obsidian Black Rim Inserts	○
ARO2651		Ruby Red Rim Inserts	○
ARO2652		Purple Passion Rim Inserts	○
ARO2653		Peacock Blue Rim Inserts	○
ARO2654		Arctic Blue Rim Inserts	○
FORKS*			
ARO2631		Single sided fork (New) includes Black castor	€358
ARO2632		Double sided fork	○
CENTRE WHEELS (DRIVING WHEELS)*			
ARO0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARO0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	◆

MOTOR TYPE*				
ARO1294		3km/h (2mph)		○
ARO1282		6km/h (4mph)		○
ARO1283		10 km/h (6mph) (includes light system and indicators)	(NC 2203)	◆
ARO1224		12 km/h (8mph) (4 pole) (includes LED light system and indicators)	(NC 2203)	€527
BATTERY*				
ARO1260		73.5 Ah Battery		◆
G TRAC				
ARO1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)		◆
CHARGER*				
ARO1343		OptiCharge Charger 8A (200 - 250V)		◆
ARO1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster		€174
LIGHTS*				
ARO1437		LED Light System with Indicators (30 LUX)	(NC 2203)	€463
ARX1006		Lights Included with Motor Package 10 & 12kmh	(NC 1282)	◆
ARO1254		Without Light system		○
TRANSPORTATION KIT*				
ARO1450		TEMPORARILY UNAVAILABLE ON TDX MODELS UNTIL FURTHER NOTICE - Preparation for DAHL-Docking station.	(NC 1300)	€596
ARO1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.	(NC 2053)	◆
ARO1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.		○
FIXATION BELT*				
ARO2113		Lap belt with auto lock (flexible mount in side rails)	(NC 2701, 2702, 2703, 2704, 2705, 2706)	◆
ARO2701		Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2025, 2027)	€117
ARO2702		Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2028, 2029, 2030, 2031, 2032, 2033)	€117
ARO2703		2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2025, 2027)	€108
ARO2704		2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2028, 2029, 2030, 2031, 2032, 2033)	€108
ARO2705		2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2025, 2027)	€108
ARO2706		2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2028, 2029, 2030, 2031, 2032, 2033)	€108
REAR VIEW MIRROR				
ARO1415		Rear view Mirror RH - Must be selected with 10 km/h Motor	(NC 1207, 1707, 2075)	○
ARO1339		Rear view Mirror LH - Must be selected with 10 km/h Motor	(NC 1208, 1708, 2076)	○
TRAY				
ARO1404		Tray Swing Away RH	(NC 1207, 1707, 2033, 2014, 2018, 2209, 2207, 2213, 2214, 2241, 2242, 2275, 2276)	€355
ARO1405		Tray Swing Away LH	(NC 1208, 1708, 2033, 2014, 2018, 2210, 2208, 2213, 2214, 2241, 2242, 2275, 2276)	€355

TOTAL PRICE	
--------------------	--

Non Conforming Options:	
ARO0080	TDX SP2 Ultra
ARZ1300	Special TDX SP2 Ultra (Extended Lead Time). Please provide details in comments section.
ARO2053	300mm seat lifter (increased STF height by 4mm) with centre of gravity tilt (0 to 45°)
ARO0339	480mm
ARO2055	10° Fixed Anterior Assist Mount (will reduce max tilt angle by 10°)
ARO2021	405mm - 505mm (increments of 25mm) (16" - 20")
ARZ2003	480mm - 555mm (Increments of 25mm) 19" - 22") - Additional build time of 15+days
ARZ2022	505mm - 610mm (increments of 25mm) (20" - 24") - Additional build time of 15+days
ARO2024	480mm - 580mm (increments of 25mm) (19" - 23")
ARO2025	405mm (16")
ARO2026	430mm (17")
ARO2027	455mm (18")
ARO2028	480mm (19")
ARO2029	505mm (20")
ARO2030	530mm (21")
ARO2031	555mm (22")
ARO2032	580mm (23")
ARO2033	640mm (24")
ARO2035	405mm (16")
ARO2036	430mm (17")
ARO2037	455mm (18")
ARO2040	530mm (21")
ARO2041	555mm (22")
ARO2042	580mm (23")
ARO2146	Matrx Flo-tech Contour Visco
ARO2147	Matrx Libra 90mm thick
ARO1980	Tapered Wedges (Pair)
ARO1981	Abductor
ARO1982	Lateral Wedge
ARO1983	Anterior Seat Wedges
ARO1985	Pelvic Obliquity Kit (Contains 2x positioning accessories)
ARO2050	High Back, height 410mm (Mesh fabric)
ARO2051	High Back height 460mm (Mesh fabric)
ARO2320	High Back, height 505mm / 20" (Mesh fabric)
ARO2550	Matrx Elite E2 Back, height 355mm / 14"
ARO2551	Matrx Elite E2 Back, height 405mm / 15"
ARO2552	Matrx Elite E2 Back, height 455mm / 18"
ARO2553	Matrx Elite E2 Back, height 505mm / 20"
ARO2554	Matrx Elite HD Back, height 405mm / 16"
ARO2555	Matrx Elite HD Back, height 505mm / 20"
ARO2556	Matrx Elite TR Back, height 405mm / 16" (Mesh fabric)
ARO2557	Matrx Elite TR Back, height 455mm / 18"(Mesh fabric)
ARO2177	Matrx Elite TR Back height 505mm (Mesh fabric)

ARO2563	Matrix Elite TR HD Back, height 505mm / 20" (Mesh fabric)
ARO2558	Matrx Elite Deep Back, height 405mm / 16" (Mesh fabric)
ARO2559	Matrx Elite Deep Back, height 455mm / 18" (Mesh fabric)
ARO2178	Matrx Elite Deep Back, height 500mm (Mesh fabric)
ARO2174	Matrx Posture Back, height 405mm (Mesh fabric max dimension 455mm wide)
ARO2175	Matrx Posture Back height 505mm (Mesh fabric max dimension 455mm wide)
ARO2560	Matrx PB Deep Back, height 405mm / 16"
ARO2561	Matrx PB Deep Back, height 505mm / 20"
ARO2562	Matrx Elite Deep HD Back, height 505mm / 20"
ARO2070	560mm (22") measured from seat pan to top of back shell
ARO2071	585mm (23") measured from seat pan to top of back shell
ARO2073	635mm (25") measured from seat pan to top of back shell
ARO2074	660mm (26") measured from seat pan to top of back shell
ARO2014	Ergonomic arm though 30mm thick
ARO2018	Ergonomic arm though 30mm thick
ARO2347	Elbow Block XS (110mm wide x 83mm high)
ARO2348	Elbow Block S (135mm wide x 90mm high)
ARO2349	Elbow Block M (145mm wide x 95mm high)
ARO2350	Elbow Block L (165mm wide x 110mm high)
ARO2351	Elbow Block XS (110mm wide x 83mm high)
ARO2352	Elbow Block S (135mm wide x 90mm high)
ARO2353	Elbow Block M (145mm wide x 95mm high)
ARO2354	Elbow Block L (165mm wide x 110mm high)
ARO2184	Hip Support LH, Quick Release, pad size 140 x 100 mm
ARO2186	Hip Support LH, Quick Release, pad size 230 x 100 mm
ARO2188	Hip Support LH, Quick Release, pad size 330 x 100 mm
ARO2185	Hip Support RH, Quick release, pad size 140 x 100 mm
ARO2187	Hip Support RH, Quick release, pad size 230 x 100 mm
ARO2189	Hip Support RH, Quick Release, pad size 330 x 100 mm
ARO2092	Onyx Auto Style headrest
ARO2324	Elan Standard Pad 255mm / 10" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW)
ARO2325	Elan Occipital Pad 230mm x 180mm / 9" x 7" w/standard cover (EOP23) complete with Elan Headrest Hardware (MEHW)
ARO0537	Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARO0539	ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm
ARO0535	Without Hanger and footplates
ARO0547	Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARO0549	ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm
ARO0545	Without Hanger and footplates
O0080	Factory Setting
ARO0505	Fixed Footplate
ARO0509	Angle Adjustable Footplate
ARO0511	Angle and depth adjustable Footplate
ARO0502	Fixed Footplate
ARO0506	Angle Adjustable Footplate

ARO0508	Angle and depth adjustable Footplate
ARO2203	DLX-REM211 (Drive & Actuators)
ARO2204	DLX-REM216 (Drive, Actuators & Lights)
ARO2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
ARO1207	Standard Bracket RH
ARO1208	Standard Bracket LH
ARO1707	Swing Away RH Bracket
ARO1708	Swing Away LH Bracket
ARO1709	Matrx Resolve Swing in/out Bracket RH
ARO1710	Matrx Resolve Swing in/out Bracket LH
ARO2209	REM500 Swing Away RH Multi-Adjustable Bracket
ARO2210	REM500 Swing Away LH Multi-Adjustable Bracket
ARO2075	Swing Away + height / width adjust. Mount on RH for user, Quad link design
ARO2076	Swing Away + height / width adjust. Mount on LH for user, Quad link design
ARO2213	Centre Mount Swing Away to RH Nucleus Bracket
ARO2214	Centre Mount Swing Away to LH Nucleus Bracket
ARO2227	Bracket RH mounted for Attendant
ARO2228	Bracket LH mounted for Attendant
ARO2219	LiNX DLX-ACU200 (Compact Dual Control)
ARO2245	LiNX DLX-IDC (Intuitive Dual Control)
ARO1207	Dual Control Standard Bracket RH
ARO1208	Dual Control Standard Bracket LH
ARO1221	IDC Dual Control Detachable Bracket RH
ARO1222	IDC Dual Control Detachable Bracket LH
ARO2231	LiNX DLX-CR400 (Compact Remote)
ARO2232	LiNX DLX-CR400LF (Compact Remote Low Force)
ARO2233	ASL MEC (Proportional MEC Joystick ASL130)
ARO2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
ARO2243	ASL Paediatric (Proportional Joystick ASL132, incl. DLX-IN500)
ARO2207	RH Mount
ARO2208	LH Mount
ARO2213	Centre Mount Swing Away to RH (Nucleus)
ARO2214	Centre Mount Swing Away to LH (Nucleus)
ARO2235	Tray Mount RH (ASL603, 100mm x 150mm)
ARO2236	Tray Mount LH (ASL603, 100mm x 150mm)
ARO2213	Centre Mount Swing Away to RH (Nucleus)
ARO2214	Centre Mount Swing Away to LH (Nucleus)
ARO2241	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARO2242	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARZ2215	Individual mounting option per request (Special). Requires additional information/specification.
ARO2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
ARO2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
ARO2257	Four-Switch Proximity Array (ASL106)
ARO2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)

ARO2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
ARO2255	Sip'n Puff Head Control (ASL109) includes Input Module and Lip Switch (ASL 301 LIP) (Combines simple sip and puff controls with head movements)
ARO2256	Sip'n Puff Control (including Sip'n Puff kit, Lip Switch ASL301LIP and Input module)
ARO2271	Input Module DLX-IN500-A (sub D, Sip'n Puff)
ARO2227	For attendant, mounted on back shroud
ARO2228	Without mounting
ARO2253	USB Charger Port DLX-USB02-A
ARO2207	RH Mount
ARO2208	LH Mount
ARO1686	Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)
ARO1687	Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)
ARO1690	Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)
ARO2254	Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)
ARO1711	Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)
ARO1712	Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)
ARO1713	Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)
ARO1714	Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)
ARO1715	Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)
ARO1716	Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)
ARO2261	Stereo Toggle - direct control for Tilt only
ARO2262	Stereo Button - direct control for Tilt only
ARO2263	Four-Way Toggle - direct control for Tilt & Recline
ARO2264	Four-Way Button - direct control for Tilt & Recline
ARO2207	Toggle / Button RH Mount
ARO2208	Toggle / Button LH Mount
ARO1282	6km/h (4mph)
ARO1437	LED Light System with Indicators (30 LUX)
ARO2113	Lap belt with auto lock (flexible mount in side rails)
ARO2701	Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ARO2702	Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ARO2703	2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ARO2704	2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ARO2705	2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)
ARO2706	2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ARO1415	Rear view Mirror RH - Must be selected with 10 km/h Motor
ARO1339	Rear view Mirror LH - Must be selected with 10 km/h Motor
ARO1404	Tray Swing Away RH
ARO1405	Tray Swing Away LH



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):