



Yes, you can:

# Powerchairs Ireland Prescription Form

NOVEMBER 2023  
Customer Ref:

Please refer to [www.invacare.ie](http://www.invacare.ie) for Local Distributor

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## TDX2 Ultra

### Standard Features

- Max user weight 136Kg (21.4 Stone)
- Crash tested ISO7176-Part 19
- Seat widths of 405mm - 610mm (16" - 24"), adjustable between 405mm - 505mm (16" - 20") and 505mm - 610mm (20" - 24")
- Seat depth of 380mm - 580mm (15" - 23") adjustable between 380mm - 480mm (15" - 19") and 480mm - 580mm (19 - 23")
- Backrest heights of 560mm - 660mm in 25mm increments
- Cantilever flip-back arm, width, height and angle adjustable with desk armpad with Startex covered visco foam
- Elan multi positioning headrest
- ULTRA Low Maxx CG power Tilt (0° to 50°) seat heights of 420mm - 480mm
- Powered Recline 90° to 170° with extended shear reduction (ESR)
- Angle adjustable legrest hanger (swing away and detachable)
- Vari F hangers with fixed footplates
- LiNX control REM211 with standard bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 355mm x 76mm (14" x 3") Black tyres
- 10km/h (6mph) 4 pole motors and 73Ah batteries
- G-Trac
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: €18,703**



Q'STRAIT  
QLK-150  
Compatible

Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ARO0080	<input type="checkbox"/> TDX SP2 Ultra	◆
ARX0080	<input type="checkbox"/> Special TDX SP2 Ultra (Please provide details in special section)	POA
<b>SEATING TYPE*</b>		
ARO2020	<input type="checkbox"/> Ultra Low Maxx Powered Positioning System	€18,703

SEAT ADJUSTMENTS*			
ARO2052		Centre of gravity Tilt (0 to 50° )	◆
ARO2053		300mm seat lifter with centre of gravity tilt (0 to 45° ) (NC 0339)	€2,787
SEAT HEIGHT SETTING*			
ARO0337		Seat-to -Floor Height 440mm	○
ARO0331		Seat-to Floor Height 460mm	○
ARO0339		Seat-to Floor Height 480mm (NC 2053)	○
ANTERIOR AND POSTERIOR TILT*			
ARO2057		Tilt Standard Mount	◆
ARO2054		5° Fixed Anterior Assist Mount	€481
ARO2055		10° Fixed Anterior Assist Mount (NC 2035)	€481
ARO2056		5° Pre-Tilt	€163
BACKREST RECLINER*			
ARO2060		Powered Recliner (90° - 170°) with Extended Shear Reduction (ESR), Standard Mount	◆
ARO2063		Powered Recliner 8° Recline mount 82° to 162° with Extended Shear Reduction (ESR)	€38
ARO2066		Powered Recline, 30° Recline mount, 60° to 140° with Extended Shear Reduction (ESR)	€38
SEAT FRAME WIDTH RANGE*			
ARO2021		405mm - 505mm (increments of 25mm) (16" - 20") (NC 2030, 2031, 2032, 2033)	○
ARX2003		480mm - 555mm (Increments of 25mm) 19" - 22") (NC 2025, 2026, 2027, 2032, 2033)	○
ARX2022		505mm - 610mm (increments of 25mm) (20" - 24") (NC 2025, 2026, 2027, 2174, 2175)	○
SEAT FRAME DEPTH RANGE*			
ARO2024		480mm - 580mm (increments of 25mm) (19" - 23") (NC 2035, 2036, 2037)	○
ARO2008		405mm - 505mm (increments of 25mm) (16" - 20") (NC 2040, 2041, 2042)	○
SEAT WIDTH*			
ARO2025		405mm (16") SEAT WIDTH (NC 2003, 2022, 2041, 2042, 2701, 2703, 2705)	○
ARO2026		430mm (17") SEAT WIDTH (NC 2003, 2022, 2041, 2042, 2701, 2703, 2705)	○
ARO2027		455mm (18") SEAT WIDTH (NC 2003, 2022, 2701, 2703, 2705)	○
ARO2028		480mm (19") SEAT WIDTH (NC 2022, 2174, 2175, 2702, 2704, 2706)	○
ARO2029		505mm (20") SEAT WIDTH (NC 2174, 2175, 2702, 2704, 2706)	○
ARO2030		530mm (21") SEAT WIDTH (NC 2021, 2035, 2036, 2174, 2175, 2702, 2704, 2706)	○
ARO2031		555mm (22") SEAT WIDTH (NC 2021, 2035, 2036, 2174, 2175, 2702, 2704, 2706)	○
ARO2032		580mm (23") SEAT WIDTH (NC 2021, 2035, 2036, 2174, 2175, 2147, 2702, 2704, 2706)	○
ARO2033		640mm (24") SEAT WIDTH (NC 2021, 2035, 2036, 2174, 2175, 1404, 1405, 2147, 2702, 2704, 2706)	○
SEAT DEPTH SETTING*			
ARO2035		405mm (16") SEAT WIDTH (NC 2024, 2030, 2031, 2032, 2033, 2055, 2701, 2702, 2703, 2704, 2705, 2706)	○
ARO2036		430mm (17") SEAT WIDTH (NC 2024, 2030, 2031, 2032, 2033)	○
ARO2037		455mm (18") SEAT WIDTH (NC 2024)	○
ARO2038		480mm (19") SEAT WIDTH	○
ARO2039		505mm (20") SEAT WIDTH	○
ARO2040		530mm (21") SEAT WIDTH	○
ARO2041		555mm (22") SEAT WIDTH (NC 2025, 2026)	○
ARO2042		580mm (23") SEAT WIDTH (NC 2025, 2026, 2146)	○
SEAT CUSHION*			
ARO2144		Matrx Vi Cushion (Mesh fabric)	€348
ARO2146		Matrx Flo-tech Contour Visco 43cm x 43cm (NC 2041, 2042, 2030, 2031, 2033)	€174
ARO2147		Matrx Libra Cushion 90mm thick (NC 2042)	€532
ARO2140		Solid seat Base, For alternative cushion options (NC 2032, 2033)	◆

BACKREST SUPPORT*				
ARO2050		High Back, height 410mm (Mesh fabric)	(NC 2073, 2074)	○
ARO2051		High Back height 460mm (Mesh fabric)	(NC 2070, 2071)	◆
ARO2320		High Back, height 505mm / 20" (Mesh fabric)		○
ARO2550		Matrx Elite Back, height 355mm / 14"	(NC 2028, 2029, 2030, 2032, 2033)	€572
ARO2551		Matrx Elite Back, height 405mm / 15"		€572
ARO2552		Matrx Elite Back, height 455mm / 18"		€572
ARO2553		Matrx Elite Back, height 505mm / 20"		€572
ARO2554		Matrx Elite HD Back, height 405mm / 16"		€741
ARO2555		Matrx Elite HD Back, height 505mm / 20"		€741
ARO2556		Matrx Elite TR Back, height 405mm / 16" (Mesh fabric)		€572
ARO2557		Matrx Elite TR Back, height 455mm / 18"(Mesh fabric)		€572
ARO2177		Matrx Elite TR Back height 505mm (Mesh fabric)		€572
ARO2563		Matrix Elite TR HD Back, height 505mm / 20" (Mesh fabric)		€741
ARO2558		Matrx Elite Deep Back, height 405mm / 16" (Mesh fabric)		€741
ARO2559		Matrx Elite Deep Back, height 455mm / 18" (Mesh fabric)		€741
ARO2178		Matrx Elite Deep Back, height 500mm (Mesh fabric)		€741
ARO2174		Matrx Posture Back, height 405mm (Mesh fabric max dimension 455mm wide)		€572
ARO2175		Matrx Posture Back height 505mm (Mesh fabric max dimension 455mm wide)		€572
ARO2560		Matrx PB Deep Back, height 405mm / 16"		€741
ARO2561		Matrx PB Deep Back, height 505mm / 20"		€741
ARO2562		Matrx Elite Deep HD Back, height 505mm / 20"		€741

BACKREST HEIGHT*				
ARO2506		505mm (20") measured from seat pan to top of back shell	(NC 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562)	○
ARO2501		535mm (21") measured from seat pan to top of back shell	(NC 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562)	○
ARO2070		560mm (22") measured from seat pan to top of back shell	(NC 2050, 2551, 2553, 2554, 2555, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563)	○
ARO2071		585mm (23") measured from seat pan to top of back shell	(NC 2050, 2551, 2553, 2554, 2555, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563)	○
ARO2072		610mm (24") measured from seat pan to top of back shell	(NC 2050, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	○
ARO2073		635mm (25") measured from seat pan to top of back shell	(NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2558, 2559, 2174, 2560)	○
ARO2074		660mm (26") measured from seat pan to top of back shell	(NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	○
ARO2502		685mm (27") measured from seat pan to top of back shell	(NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	◆

ARMREST*				
ARO2095		Cantilever Flip-Back Arm, width, height and angle adjustable (240mm - 330mm)		○
ARO2096		Cantilever Flip-Back Arm, width, height and angle adjustable (320mm - 405mm)		○
ARO2097		Reclining Arm, Flip-Back, adjustable height (240mm - 330mm)		○
ARO0690		Without Armrest		○

ARMPAD LH*			
ARO2012		Desk arm pad with startex cover visco foam (255mm) 40m thick	○
ARO2013		Full length armpad with startex cover foam (355mm) 40mm thick	○
ARO2014		Ergonomic arm though 30mm thick	(NC 1404, 1405, 2347, 2348, 2349, 2350) €223
ARO2015		Without LH arm pad	◆
ARMPAD RH*			
ARO2016		Desk arm pad with startex cover visco foam (255mm) 40mm thick	○
ARO2017		Full length armpad with startex cover foam (355mm) 40mm thick	○
ARO2018		Ergonomic arm though 30mm thick	(NC 1404, 1405, 2351, 2352, 2353, 2354) €223
ARO2019		Without RH arm pad	€148
ARMPAD ACCESSORIES & MOUNTING OPTIONS			
ARO2080		Flat Hand Pad Left (in combination with arm trough only): to be mounted on the opposite side of remote	€60
ARO2083		Flat Hand Pad Right (in combination with arm trough only): to be mounted on the opposite side of remote	€60
MULTI-AXIX UPPER EXTREMITY SUPPORTS			
ARO2082		Multi Axis Upper Extremity Support (MACES) Left, raises arm pad by 43mm	€362
ARO2085		Multi Axis Upper Extremity Support (MACES) Right, raises arm pad by 43mm	€362
ARMREST HEIGHT SETTING LH*			
O0001		Armrest height setting LH - Factory Setting	○
O0002		Customised 230mm - Armrest height setting LH	○
O0003		Customised 240mm - Armrest height setting LH	○
O0004		Customised 250mm - Armrest height setting LH	○
O0005		Customised 260mm - Armrest height setting LH	○
O0006		Customised 270mm - Armrest height setting LH	○
O0007		Customised 280mm - Armrest height setting LH	○
O0008		Customised 290mm - Armrest height setting LH	○
O0009		Customised 300mm - Armrest height setting LH	○
O0010		Customised 310mm - Armrest height setting LH	○
O0011		Customised 320mm - Armrest height setting LH	○
O0012		Customised 330mm - Armrest height setting LH	○
O0013		Customised 340mm - Armrest height setting LH	○
O0014		Customised 350mm - Armrest height setting LH	○
O0015		Customised 360mm - Armrest height setting LH	○
O0016		Customised 370mm - Armrest height setting LH	○
O0017		Customised 380mm - Armrest height setting LH	○
O0018		Customised 390mm - Armrest height setting LH	○
O0019		Customised 405mm - Armrest height setting LH	○

<b>ARMREST HEIGHT SETTING RH*</b>				
O0020		Armrest height setting RH - Factory Setting		○
O0021		Customised 230mm - Armrest height setting RH		○
O0022		Customised 240mm - Armrest height setting RH		○
O0023		Customised 250mm - Armrest height setting RH		○
O0024		Customised 260mm - Armrest height setting RH		○
O0025		Customised 270mm - Armrest height setting RH		○
O0026		Customised 280mm - Armrest height setting RH		○
O0027		Customised 290mm - Armrest height setting RH		○
O0028		Customised 300mm - Armrest height setting RH		○
O0029		Customised 310mm - Armrest height setting RH		○
O0030		Customised 320mm - Armrest height setting RH		○
O0031		Customised 330mm - Armrest height setting RH		○
O0032		Customised 340mm - Armrest height setting RH		○
O0033		Customised 350mm - Armrest height setting RH		○
O0034		Customised 360mm - Armrest height setting RH		○
O0035		Customised 370mm - Armrest height setting RH		○
O0036		Customised 380mm - Armrest height setting RH		○
O0037		Customised 390mm - Armrest height setting RH		○
O0038		Customised 405mm - Armrest height setting RH		○
<b>ELBOW BLOCK LH PADS</b>				
ARO2347		Elbow Block XS (110mm wide x 83mm high)	(NC 2014)	€521
ARO2348		Elbow Block S (135mm wide x 90mm high)	(NC 2014)	€521
ARO2349		Elbow Block M (145mm wide x 95mm high)	(NC 2014)	€521
ARO2350		Elbow Block L (165mm wide x 110mm high)	(NC 2014)	€521
<b>ELBOW BLOCK LH MOUNTING HARDWARE</b>				
ARO2361		Vertical 3"/ Horizontal 5"		○
ARO2362		Vertical 3"/ Horizontal 7"		○
ARO2363		Vertical 5"/ Horizontal 5"		○
ARO2364		Vertical 5"/ Horizontal 7"		○
<b>ELBOW BLOCKS RH PADS</b>				
ARO2351		Elbow Block XS (110mm wide x 83mm high)	(NC 2018)	€521
ARO2352		Elbow Block S (135mm wide x 90mm high)	(NC 2018)	€521
ARO2353		Elbow Block M (145mm wide x 95mm high)	(NC 2018)	€521
ARO2354		Elbow Block L (165mm wide x 110mm high)	(NC 2018)	€521
<b>ELBOW BLOCK RH MOUNTING HARDWARE</b>				
ARO2371		Vertical 3"/ Horizontal 5"		○
ARO2372		Vertical 3"/ Horizontal 7"		○
ARO2373		Vertical 5"/ Horizontal 5"		○
ARO2374		Vertical 5"/ Horizontal 7"		○
<b>HIP SUPPORTS LH (Not available with 410mm seat depth and Manual Centre legrests)</b>				
ARO2184		Hip Support LH pad size 140 x 100 mm	(NC 2701, 2702, 2703, 2704, 2705, 2706)	€162
ARO2114		Hip Support LH, pad size 180 x 150 mm	(NC 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARO2186		Hip Support LH, pad size 230 x 100 mm	(NC 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARO2188		Hip Support LH, pad size 330 x 100 mm	(NC 2701, 2702, 2703, 2704, 2705, 2706)	€272

<b>HIP SUPPORTS RH (Not available with 410mm seat depth and Manual Centre legrests)</b>			
ARO2185		Hip Support RH, pad size 140 x 100 mm (NC 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARO2115		Hip Support RH, , pad size 180 x 150 mm (NC 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARO2187		Hip Support RH, , pad size 230 x 100 mm	€272
ARO2189		Hip Support RH, pad size 330 x 100 mm	€272
<b>HEADREST*</b>			
ARO2092		Onyx Auto Style headrest (NC 2255, 2230, 2247)	€207
ARO2324		Elan Standard Pad 255mm / 10" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	€334
ARO2325		Elan Occipital Pad 230mm x 180mm / 9" x 7" w/standard cover (EOP23) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	€334
ARO1495		Head Control or equivalent Headrest, provided independently from this order (ACCORDING TO THE REQUIREMENTS OF ISO 7176-19, A HEADREST IS COMPULSORY FOR USING THE WHEELCHAIR AS VEHICLE SEAT) (NC 2229)	○
<b>LATERAL SUPPORTS LH</b>			
ARO2343		110mm wide x 83mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€402
ARO2345		135mm wide x 90mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€402
ARO2043		145mm wide x 95mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€362
ARO2045		165mm wide x 110mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€428
<b>LATERAL SUPPORTS RH</b>			
ARO2344		110mm wide x 83mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€402
ARO2346		135mm wide x 90mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€402
ARO2044		145mm wide x 95mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€362
ARO2046		165mm wide x 110mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	€427
<b>LATERAL SUPPORT ACCESSORIES</b>			
ARO2049		Additional telescopic link for Lateral support, 25mm, pair (NC 2558, 2559, 2178, 2560, 2561, 2562)	€170
<b>CENTRE MOUNTED LEGREST</b>			
ARO0586		Hanger, centre mounted, angle settable (0° - 90°): knee to heel length 310 - 410mm	€357
ARO0589		Hangers, centre mounted, length compensating, powered, LNX (0° - 90° ): KNEE TO HEL LENGTH 340 - 410mm (NC 2261, 2262, 2263, 2264)	€1,985
ARO0592		LNX Power centre mount legrest with telescoping foot platform (lowers to the ground) inc calf pads: Length 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting (NC 0559, 0516, 0517, 2056, 2263, 2264, 2365)	€3,022
<b>FOOTPLATES / FOOT PLATFORM</b>			
ARO2365		Angle & width adjustable footplates Small (127mm width x 190mm depth / 5" x 7.5") (NC 0592)	€106
ARO0559		Angle & width adjustable footplate (240mm depth) (NC 0592)	○
ARO0516		Angle & width adjustable footplates (290mm depth) (NC 0592)	€274
ARO0517		Angle adjustable foot platform (240mm depth, 305mm width) (NC 0592)	€274

**CENTRE MOUNT LEGREST HEIGHT SETTINGS 0586 & 0589**

O0039		Factory Setting	○
O0040		Customised 310mm - centre mount legrest height settings	○
O0041		Customised 325mm - centre mount legrest height settings	○
O0042		Customised 335mm - centre mount legrest height settings	○
O0043		Customised 355mm - centre mount legrest height settings	○
O0044		Customised 355mm - centre mount legrest height settings	○
O0045		Customised 364mm - centre mount legrest height settings	○
O0046		Customised 375mm - centre mount legrest height settings	○
O0047		Customised 385mm - centre mount legrest height settings	○
O0048		Customised 410mm - centre mount legrest height settings	○

**CENTRE MOUNT W TELESCOPING FOOT PLATFORM HEIGHT SETTINGS**

O0049		Factory Setting - centre mount with telescopic foot platform height settings	○
O0050		Customised 380mm - centre mount legrest height settings	○
O0051		Customised 393mm - centre mount legrest height settings	○
O0052		Customised 406mm - centre mount legrest height settings	○
O0053		Customised 419mm - centre mount legrest height settings	◆

**SIDE MOUNTED HANGERS LH**

ARO536		Vari F: Angle adjustable with tools 0° - 70°, length 290mm - 460mm (NC 0505)	€134
ARO537		Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm (NC 0505)	€128
ARO539		ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm (NC 0505, 2261, 2262, 2263, 2264)	€1,225
ARO535		Without Hanger and footplates (NC 0505, 0509, 0511)	○

**VARI HANGER LH HEIGHT SETTING**

O0054		Factory Setting - Vari hanger LH height setting	○
O0055		Customised 270mm - Vari hanger LH height setting	○
O0056		Customised 280mm - Vari hanger LH height setting	○
O0057		Customised 290mm - Vari hanger LH height setting	○
O0058		Customised 300mm - Vari hanger LH height setting	○
O0059		Customised 310mm - Vari hanger LH height setting	○
O0060		Customised 320mm - Vari hanger LH height setting	○
O0061		Customised 330mm - Vari hanger LH height setting	○
O0062		Customised 340mm - Vari hanger LH height setting	○
O0063		Customised 350mm - Vari hanger LH height setting	○
O0064		Customised 360mm - Vari hanger LH height setting	○
O0065		Customised 370mm - Vari hanger LH height setting	○
O0066		Customised 380mm - Vari hanger LH height setting	○
O0067		Customised 390mm - Vari hanger LH height setting	○
O0068		Customised 400mm - Vari hanger LH height setting	○
O0069		Customised 410mm - Vari hanger LH height setting	○
O0070		Customised 420mm - Vari hanger LH height setting	○
O0071		Customised 430mm - Vari hanger LH height setting	○

**SIDE MOUNTED HANGERS RH**

ARO0546		Vari F: Angle adjustable with tools 0° - 70°, length 290mm - 460mm (NC 0502)	€134
ARO0547		Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm (NC 0502)	€128
ARO0549		ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm (NC 0502, 2261, 2262, 2263, 2264)	€1,225
ARO0545		Without Hanger and footplates (NC 0502, 0506, 0508)	○

**VARI HANGER RH HEIGHT SETTING**

O0080		Factory Setting - Vari hanger RH height setting	○
O0081		Customised 270mm - Vari hanger RH height setting	○
O0082		Customised 280mm - Vari hanger RH height setting	○
O0083		Customised 290mm - Vari hanger RH height setting	○
O0084		Customised 300mm - Vari hanger RH height setting	○
O0085		Customised 310mm - Vari hanger RH height setting	○
O0086		Customised 320mm - Vari hanger RH height setting	○
O0087		Customised 330mm - Vari hanger RH height setting	○
O0088		Customised 340mm - Vari hanger RH height setting	○
O0089		Customised 350mm - Vari hanger RH height setting	○
O0090		Customised 360mm - Vari hanger RH height setting	○
O0091		Customised 370mm - Vari hanger RH height setting	○
O0092		Customised 380mm - Vari hanger RH height setting	○
O0093		Customised 390mm - Vari hanger RH height setting	○
O0094		Customised 400mm - Vari hanger RH height setting	○
O0095		Customised 410mm - Vari hanger RH height setting	○
O0096		Customised 420mm - Vari hanger RH height setting	○
O0097		Customised 430mm - Vari hanger RH height setting	○



INDIVIDUAL FOOTPLATE TYPE LH				
ARO0505		Fixed Footplate	(NC 0535, 0537, 0539)	€75
ARO0509		Angle Adjustable Footplate	(NC 0535)	○
ARO0511		Angle and depth adjustable Footplate	(NC 0535)	◆
INDIVIDUAL FOOTPLATE TYPE RH				
ARO0502		Fixed Footplate	(NC 0545, 0547, 0549)	€75
ARO0506		Angle Adjustable Footplate	(NC 0545)	○
ARO0508		Angle and depth adjustable Footplate	(NC 0545)	○
CONTROLS*				
ARO2203		REM211 Drive & Actuators	(NC 1437, 2231, 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2209, 2210, 2213, 2214, 2075, 2076)	◆
ARO2204		REM216 Drive, Actuators & Lights	(NC 2231, 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2209, 2210, 2213, 2214, 2075, 2076)	€495
ARO2205		REM400, touch-screen Interface, built in Bluetooth HID	(NC 2209, 2210, 2213, 2215, 2075, 2076)	€909
ARO2206		REM500, without joystick for alternative access, built in Bluetooth HID	(NC 1207, 1208, 1707, 1708, 2227, 2228, 2075, 2076)	○
CONTROL BRACKETRY*				
ARO1207		Standard Bracket RH	(NC 1404)	○
ARO1208		Standard Bracket LH	(NC 1405)	€151
ARO1707		Swing Away + height adjustable Bracket RH High end	(NC 1404)	€151
ARO1708		Swing Away + height adjustable Bracket LH High end	(NC 1405)	€244
ARO2075		Swing Away + height / width adjust mount on RH for user, Quad link design	(NC 1415)	€244
ARO2076		Swing Away + height / width adjust mount on LH for user, Quad link design	(NC 1339)	€71
ARO2227		Centre Mount swing-away Bracket to RH (Nucleus)	(NC 2219, 2245)	€71
ARO2228		Centre Mount swing-away Bracket to LH (Nucleus)	(NC 2219, 2245)	€71
ARO2209		REM500 Swing away RH multi-adjustable	(NC 1404)	○
ARO2210		REM500 Swing away LH multi-adjustable	(NC 1405)	€739
ARO2213		REM500 Centre mount swing away to RH Nucleus	(NC 1404, 1405)	€739
ARO2214		REM500 Centre mount swing away to LH Nucleus	(NC 1404, 1405)	€244
ARO2075		Swing Away + height / width adjustable Bracket Mount on RH, Quad link design	(NC 1415)	€244
ARO2076		Swing Away + height / width adjustable Bracket Mount on RH Quad link design	(NC 1339)	POA
ARX2239		Individual mounting option per (special) - Requires additional information / specification		€329
REM400 TOGGLE OPTION				
ARO2270		Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial)	(NC 2203, 2204)	€3,106
CHIN CONTROL				
ARO2229		Proportional Chin Joystick (ASL138) (Requires Primary Remote REM400 or REM500) (Please select chin control below)	(NC 1495, 2203, 2204, 2230, 2247, 2233, 2257, 2240, 1287, 1286)	€3,106
CHIN CONTROL BRACKETRY				
ARO1289		Chin Control Swing Away Bracket RH Manual	(NC 1495, 2203, 2204, 2230, 2247, 2233, 2257, 2240, 1287, 1286)	○
ARO1288		Chin Control Swing Away Bracket LH Manual	(NC 1495, 2203, 2204, 2230, 2247, 2233, 2257, 2240, 1287, 1286)	○
ARO1287		Chin Control Swing Away Bracket RH Powered	(NC 1495, 2203, 2204, 2230, 2247, 2233, 2257, 2240, 1287, 1286)	€2,500
ARO1286		Chin Control Swing Away Bracket LH Powered	(NC 1495, 2203, 2204, 2230, 2247, 2233, 2257, 2240, 1287, 1286)	€1,306
DUAL CONTROL				
ARO2219		Compact Dual Control (DLX-ACU200) (Please select Dual control bracket below)	(NC 1221, 1222, 2227, 2228, 2245)	€356
DUAL CONTROL BRACKETRY				
ARO1207		Dual Control Standard Bracket RH	(NC 1221, 1222, 2227, 2228, 2245)	○
ARO1208		Dual Control Standard Bracket LH	(NC 1221, 1222, 2227, 2228, 2245)	○

INVACARE LINX IDC DUAL CONTROL			
ARO2245		IDC Dual control (Please select AQE1221 OR AQE1222 for bracket option below)	(NC 2219, 2227, 2228) €748
INVACARE LINX IDC DUAL CONTROL BRACKETRY			
ARO1221		IDC Dual control Detachable Bracket RH	(NC 2219, 2227, 2228) ○
ARO1222		IDC Dual control Detachable Bracket LH	(NC 2219, 2227, 2228) ○
LINX COMPACT REMOTE			
ARO2231		Proportional Compact Remote DLX-CR400 (Please select Proportional Bracket)	(NC 2203, 2204) €1,199
ARO2232		Proportional Compact Remote Low Force DLX-CR400LF (Please select Proportional Bracketry)	(NC 2203, 2204) €1,199
LINX COMPACT REMOTE BRACKET			
ARO2207		Proportional Compact RH Mount Bracket	(NC 2203, 2204) ○
ARO2208		Proportional Compact LH Mount Bracket	(NC 2203, 2204) ○
ARO2213		Proportional Compact Centre Mount swing away to RH (Nucleus)	(NC 2203, 2204) €739
ARO2214		Proportional Compact Centre Mount swing away to LH (Nucleus)	(NC 2203, 2204) €739
ARX2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2203, 2204) POA
LINX HEAD CONTROL			
ARO2230		Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M4	(NC 2203, 2204) €4,498
ARO2247		Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings	(NC 2203, 2204) €5,247
LINX MICRO EXTREMITY CONTROL			
ARO2233		MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry)	(NC 2203, 2204, 2229, 2230, 2247, 2257, 2240) €4,230
LINX MICRO EXTREMITY CONTROL BRACKET			
ARO2235		MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please also select MEC Bracketry)	○
ARO2236		MEC Tray Mount LH (ASL603, 100mm x 150mm)	○
ARO2213		Centre Mount, swing away to RH, Nucleus	(NC 1404, 1405) €739
ARO2214		Centre Mount, swing away to LH, Nucleus	(NC 1404, 1405) €739
ARO2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405) €739
ARO2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405) €739
ARO2215		Individual mounting option per request (Special) Requires additional information / specification	€2,683
LINX HEAVY DUTY REMOTE			
ARO2240		Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LiNX Heavy Duty Bracketry)	(NC 2203, 2204, 2229, 2230, 2247, 2257, 2233) €2,683
LINX HEAVY DUTY REMOTE BRACKET			
ARO2235		Tray Mount RH, (ASL603, 100mm x 150mm)	○
ARO2236		Tray Mount LH, (ASL603, 100mm x 150mm)	○
ARO2213		Centre Mount, swing away to RH, Nucleus	(NC 1404, 1405) €739
ARO2214		Centre Mount, swing away to LH, Nucleus	(NC 1404, 1405) €739
ARO2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405) €739
ARO2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405) €739
ARX2215		Individual mounting option per request (Special) Requires additional information / specification	POA

LINX 10-WAY SWITCH			
ARO2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch) (Please select 10 way switch bracket below)	€513
LINX 10-WAY SWITCH BRACKET			
ARO2217		10-Way Switch Bracket RH, Mounted for User	○
ARO2218		10-Way Switch Bracket LH, Mounted for User	○
ARO2227		10-Way Switch Bracket RH, Mounted for Attendant	€71
ARO2228		10-Way Switch Bracket LH, Mounted for Attendant	€71
LINX FOUR -SWITCH PROXIMITY ARRAY			
ARO2257		Four-Switch Proximity Array (ASL106) NEW! (NC 2203, 2204, 2229, 2230, 2092, 2324, 2325, 2257, 2240)	€6,746
ARO2275		Centre Mount, swing-away to RH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm) (NC 1404, 1405)	€739
ARO2276		Centre Mount, swing-away to LH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm) (NC 1404, 1405)	€739
ARO2215		Individual mounting option per request (Special) - Requires additional information / specification	POA
LINX WIRELESS ACCESSORIES FOR ATOM			
ARO2248		Wireless Dual Switch Receiver (ASL557-2) (Jack plugs) NEW!	€348
ARO2249		Wireless Triple Switch Receiver (ASL557-3) (Sub D plug) NEW!	€833
ARO2250		Wireless Mouse Emulator (ASL558) NEW!	€736
LINX PEDIATRIC COMPACT			
ARO2243		Proportional Paediatric Compact (ASK132), Includes input Module NEW! (NC 2271, 2203, 2204, 2255, 2256)	€2,438
ARO2235		Tray Mount RH (ASL603, 100mm x 150mm)	○
ARO2236		Tray Mount LH (ASL603, 100mm x 150mm)	○
ARO2213		Centre Mount, swing-away to RH, Nucleus (NC 1404, 1405)	€739
ARO2214		Centre Mount, swing-away to LH, Nucleus (NC 1404, 1405)	€739
ARO2241		Centre Mount, swing-away to RH, Nucleus (incl Tray (ASL603, 100mm x 150mm) (NC 1404, 1405)	€739
ARO2242		Centre Mount, swing-away to RH, Nucleus (incl Tray (ASL603, 100mm x 150mm) (NC 1404, 1405)	€739
ARO2215		Individual mounting option per request (Special) - Requires additional information / specification	POA
LINX SIP N PUFF HEAD CONTROL			
ARO2255		Sip'n Puff Head Control (ASL109) NEW! includes Input Module and Lip Switch (ASL 301 LIP) (Combines simple sip and puff controls with head movements) (NC 2271, 2203, 2204, 2092, 2324, 2325, 2243)	€6,319
LINX SIP N PUFF CONTROL & INPUT MODULE			
ARO2256		Sip'n Puff Control NEW! (including Sip'n Puff kit, Lip Switch ASL301LIP and Input module) (NC 2271, 2203, 2204, 2243)	€878
LINX INTERFACE OPTIONS			
ARO2271		Input Module DLX-IN500-A (sub D, Sip'n Puff) NEW! (NC 2243, 2255, 2256, 2203, 2204)	€814
SPECIAL ALTERNATIVE JOYSTICKS & INPUTS			
ARX2216		Special alternative Joystick/Input and bracket as per request	POA
LINX USB CHARGER PORT			
ARO2253		USB Charger Port DLX-USB02-A (NC 2701, 2702, 2703, 2704, 2705, 2706)	€162
LINX USB CHARGER PORT BRACKET			
ARO2207		USB Charger Port RH Mount	○
ARO2208		USB Charger Port LH Mount	○
OUTPUT MODULE			
ARO2272		LiNX Output Module (OUT500)	€833

<b>FURTHER LiNX COMPONENTS / ACCESSORIES</b>		
ARO2300		To be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected) ○
ARZ2301		To be connected and programmed (special) POA
<b>BUDDY BUTTONS</b>		
ARO1685		Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) €90
ARO1686		Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) €90
ARO1687		Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) €90
ARO1690		Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) €90
<b>EGG SWITCH</b>		
ARO2254		Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected) €90
<b>PIKO BUTTONS</b>		
ARO1711		Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) €161
ARO1712		Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) €161
ARO1713		Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) €161
ARO1714		Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) €161
ARO1715		Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) €161
ARO1716		Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) €161
<b>TOGGLES / BUTTONS</b>		
ARO2261		Stereo Toggle - direct control for Tilt only €263
ARO2262		Stereo Button - direct control for Tilt only €263
ARO2263		Four-Way Toggle - direct control for Tilt & Recline €342
ARO2264		Four-Way Button - direct control for Tilt & Recline €365
ARO2207		Toggle / Button RH Mount ○
ARO2208		Toggle / Button LH Mount ○
ARO2225		Without Mounting ○
<b>BODYPOINT JOYSTICK KNOBS</b>		
ARO2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A €75
ARO2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A €75
ARO2723		Domed Rubber, Joystick Handle, Bodypoint PC107A €75
ARO2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A €75
<b>REMOTE STOP SWITCH</b>		
ARO2258		Remote Stop Switch (ASL504) €1,070
<b>SHROUD COLOUR*</b>		
ARO0028		Wet Black ○
ARO2601		SpaceStation Silver ○
ARO2602		Island Blue ○
ARO2603		Grape JellyBean (Purple) ○
ARO2607		Rockstar Red ○
ARO2608		Lights Out Black - Matt ○
ARO2609		Invacare Blue ○

<b>RIM COLOUR*</b>			
ARO2611		SpaceStation Silver	○
ARO2612		Island Blue	○
ARO2613		Grape JellyBean (Purple)	○
ARO2617		Rockstar Red	○
ARO2618		Lights Out Black - Matt	○
ARO2619		Invacare Blue	○
ARO2620		Wet Black	○
ARO2621		Without Rim insert	○
<b>FORKS*</b>			
ARO2631		Single sided fork (New) includes Black castor	€358
ARO2632		Double sided fork	○
<b>CENTRE WHEELS (DRIVING WHEELS)*</b>			
ARO0704		355mm x 76mm (14" x 3") Black Tyres, Pneumatic	◆
ARO0733		355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	○
<b>MOTOR TYPE*</b>			
ARO1294		3km/h (2mph)	○
ARO1282		6km/h (4mph)	○
ARO1285		10 km/h (6mph) (includes light system and indicators - please select option ARM1367 or ARM1482)	◆
ARO1224		12 km/k (8mph) (4 pole) (includes LED light system and indicators)	€527
<b>BATTERY*</b>			
ARO1260		73.5 Ah Battery	◆
<b>PERFORMANCE ACCESSORIES</b>			
ARO1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	◆
<b>CHARGER*</b>			
ARO1343		OptiCharge Charger 8A (200 - 250V)	◆
ARO1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€174
<b>LIGHTS*</b>			
ARO1437		Light system with indicators & LED Headlights (Only available with REM24SD or G90A joystick controls)	(NC 2203) €463
ARO1254		Without Light system	◆
<b>TRANSPORTATION KIT*</b>			
ARO1450		Preparation for DAHL Docking Station: Docking station for car implementation, please order docking station separately. (Adequate headrest is compulsory.	€596
ARO1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	(NC 2053) ◆
ARO1401		6-Point-tie--down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19-2001) (Adequate headrest is compulsory)	○
<b>FIXATION BELT*</b>			
ARO2113		Lap belt with auto lock (flexible mount in side rails)	(NC 2701, 2702, 2703, 2704, 2705, 2706) ◆

<b>BODYPOINT PELVIC BELTS (includes mounting kit)</b>			
ARO2701		Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2025, 2027) €117
ARO2702		Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2028, 2029, 2030, 2031, 2032, 2033) €117
ARO2703		2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2025, 2027) €108
ARO2704		2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2028, 2029, 2030, 2031, 2032, 2033) €108
ARO2705		2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2025, 2027) €108
ARO2706		2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)	(NC 2035, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113, 2028, 2029, 2030, 2031, 2032, 2033) €108
<b>REAR VIEW MIRROR</b>			
ARO1415		Rear view Mirror RH - Must be selected with 10 km/h Motor	(NC 1207, 1707, 2075) ○
ARO1339		Rear view Mirror LH - Must be selected with 10 km/h Motor	(NC 1208, 1708, 2076) ○
<b>TRAY</b>			
ARO1404		Tray Swing Away RH	(NC 1207, 1707, 2033, 2014, 2018, 2209, 2207, 2213, 2214, 2241, 2242, 2275, 2276) €355
ARO1405		Tray Swing Away LH	(NC 1208, 1708, 2033, 2014, 2018, 2210, 2208, 2213, 2214, 2241, 2242, 2275, 2276) €355

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

ARO2053	300mm seat lifter with centre of gravity tilt (0 to 45° )
ARO0339	Seat-to Floor Height 480mm
ARO2055	10° Fixed Anterior Assist Mount
ARO2056	5° Pre-Tilt
ARO2021	405mm - 505mm (increments of 25mm) (16" - 20")
ARX2003	480mm - 555mm (Increments of 25mm) 19" - 22")
ARX2022	505mm - 610mm (increments of 25mm) (20" - 24")
ARO2024	480mm - 580mm (increments of 25mm) (19" - 23")
ARO2025	405mm (16") SEAT WIDTH
ARO2026	430mm (17") SEAT WIDTH
ARO2027	455mm (18") SEAT WIDTH
ARO2028	480mm (19") SEAT WIDTH
ARO2029	505mm (20") SEAT WIDTH
ARO2030	530mm (21") SEAT WIDTH
ARO2031	555mm (22") SEAT WIDTH
ARO2032	580mm (23") SEAT WIDTH
ARO2033	640mm (24") SEAT WIDTH
ARO2035	405mm (16") SEAT WIDTH
ARO2036	430mm (17") SEAT WIDTH
ARO2037	455mm (18") SEAT WIDTH
ARO2040	530mm (21") SEAT WIDTH
ARO2041	555mm (22") SEAT WIDTH
ARO2042	580mm (23") SEAT WIDTH
ARO2146	Matrx Flo-tech Contour Visco 43cm x 43cm
ARO2147	Matrx Libra Cushion 90mm thick
ARO2050	High Back, height 410mm (Mesh fabric)
ARO2051	High Back height 460mm (Mesh fabric)
ARO2320	High Back, height 505mm / 20" (Mesh fabric)
ARO2550	Matrx Elite Back, height 355mm / 14"
ARO2551	Matrx Elite Back, height 405mm / 15"
ARO2552	Matrx Elite Back, height 455mm / 18"
ARO2553	Matrx Elite Back, height 505mm / 20"
ARO2554	Matrx Elite HD Back, height 405mm / 16"
ARO2555	Matrx Elite HD Back, height 505mm / 20"
ARO2556	Matrx Elite TR Back, height 405mm / 16" (Mesh fabric)
ARO2557	Matrx Elite TR Back, height 455mm / 18"(Mesh fabric)
ARO2177	Matrx Elite TR Back height 505mm (Mesh fabric)
ARO2563	Matrx Elite TR HD Back, height 505mm / 20" (Mesh fabric)
ARO2558	Matrx Elite Deep Back, height 405mm / 16" (Mesh fabric)
ARO2559	Matrx Elite Deep Back, height 455mm / 18" (Mesh fabric)
ARO2178	Matrx Elite Deep Back, height 500mm (Mesh fabric)
ARO2174	Matrx Posture Back, height 405mm (Mesh fabric max dimension 455mm wide)
ARO2175	Matrx Posture Back height 505mm (Mesh fabric max dimension 455mm wide)

ARO2560	Matrx PB Deep Back, height 405mm / 16"
ARO2561	Matrx PB Deep Back, height 505mm / 20"
ARO2562	Matrx Elite Deep HD Back, height 505mm / 20"
ARO2070	560mm (22") measured from seat pan to top of back shell
ARO2071	585mm (23") measured from seat pan to top of back shell
ARO2073	635mm (25") measured from seat pan to top of back shell
ARO2074	660mm (26") measured from seat pan to top of back shell
ARO2014	Ergonomic arm though 30mm thick
ARO2018	Ergonomic arm though 30mm thick
ARO2347	Elbow Block XS (110mm wide x 83mm high)
ARO2348	Elbow Block S (135mm wide x 90mm high)
ARO2349	Elbow Block M (145mm wide x 95mm high)
ARO2350	Elbow Block L (165mm wide x 110mm high)
ARO2351	Elbow Block XS (110mm wide x 83mm high)
ARO2352	Elbow Block S (135mm wide x 90mm high)
ARO2353	Elbow Block M (145mm wide x 95mm high)
ARO2354	Elbow Block L (165mm wide x 110mm high)
ARO2184	Hip Support LH pad size 140 x 100 mm
ARO2186	Hip Support LH, pad size 230 x 100 mm
ARO2188	Hip Support LH, pad size 330 x 100 mm
ARO2185	Hip Support RH, pad size 140 x 100 mm
ARO2187	Hip Support RH, , pad size 230 x 100 mm
ARO2189	Hip Support RH, pad size 330 x 100 mm
ARO2092	Onyx Auto Style headrest
ARO2324	Elan Standard Pad 255mm / 10" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW)
ARO2325	Elan Occipital Pad 230mm x 180mm / 9" x 7" w/standard cover (EOP23) complete with Elan Headrest Hardware (MEHW)
ARO1495	Head Control or equivalent Headrest, provided independently from this order (ACCORDING TO THE REQUIREMENTS OF ISO 7176-19, A HEADREST IS COMPULSORY FOR USING THE WHEELCHAIR AS VEHICLE SEAT)
ARO0592	LNx Power centre mount legrest with telescoping foot platform (lowers to the ground) inc calf pads: Length 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting
ARO2365	Angle & width adjustable footplates Small (127mm width x 190mm depth / 5" x 7.5")
ARO0559	Angle & width adjustable footplate (240mm depth)
ARO0516	Angle & width adjustable footplates (290mm depth)
ARO0517	Angle adjustable foot platform (240mm depth, 305mm width)
ARO0537	Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARO0539	ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm
ARO0535	Without Hanger and footplates
ARO0547	Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARO0549	ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm
ARO0545	Without Hanger and footplates
ARO0505	Fixed Footplate
ARO0509	Angle Adjustable Footplate
ARO0511	Angle and depth adjustable Footplate
ARO0502	Fixed Footplate
ARO0506	Angle Adjustable Footplate



ARO0508	Angle and depth adjustable Footplate
ARO2203	REM211 Drive & Actuators
ARO2204	REM216 Drive, Actuators & Lights
ARO1207	Standard Bracket RH
ARO1208	Standard Bracket LH
ARO1707	Swing Away + height adjustable Bracket RH High end
ARO1708	Swing Away + height adjustable Bracket LH High end
ARO2075	Swing Away + height / width adjust mount on RH for user, Quad link design
ARO2076	Swing Away + height / width adjust mount on LH for user, Quad link design
ARO2227	Centre Mount swing-away Bracket to RH (Nucleus)
ARO2228	Centre Mount swing-away Bracket to LH (Nucleus)
ARO2209	REM500 Swing away RH multi-adjustable
ARO2210	REM500 Swing away LH multi-adjustable
ARO2213	REM500 Centre mount swing away to RH Nucleus
ARO2214	REM500 Centre mount swing away to LH Nucleus
ARO2075	Swing Away + height / width adjustable Bracket Mount on RH, Quad link design
ARO2076	Swing Away + height / width adjustable Bracket Mount on RH Quad link design
ARO2229	Proportional Chin Joystick (ASL138) (Requires Primary Remote REM400 or REM500) (Please select chin control below)
ARO1287	Chin Control Swing Away Bracket RH Powered
ARO1286	Chin Control Swing Away Bracket LH Powered
ARO2219	Compact Dual Control (DLX-ACU200) (Please select Dual control bracket below)
ARO1207	Dual Control Standard Bracket RH
ARO1208	Dual Control Standard Bracket LH
ARO2245	IDC Dual control (Please select AQE1221 OR AQE1222 for bracket option below)
ARO1221	IDC Dual control Detachable Bracket RH
ARO1222	IDC Dual control Detachable Bracket LH
ARO2231	Proportional Compact Remote DLX-CR400 (Please select Proportional Bracket)
ARO2232	Proportional Compact Remote Low Force DLX-CR400LF (Please select Proportional Bracketry)
ARO2207	Proportional Compact RH Mount Bracket
ARO2208	Proportional Compact LH Mount Bracket
ARO2213	Proportional Compact Centre Mount swing away to RH (Nucleus)
ARO2214	Proportional Compact Centre Mount swing away to LH (Nucleus)
ARX2215	Individual mounting option per request (Special) Requires additional information / specification
ARO2230	Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M4
ARO2247	Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings
ARO2233	MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry)
ARO2213	Centre Mount, swing away to RH, Nucleus
ARO2214	Centre Mount, swing away to LH, Nucleus
ARO2241	Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARO2242	Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARO2215	Individual mounting option per request (Special) Requires additional information / specification
ARO2240	Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LiNX Heavy Duty Bracketry)
ARO2213	Centre Mount, swing away to RH, Nucleus
ARO2214	Centre Mount, swing away to LH, Nucleus

ARO2241	Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARO2242	Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARX2215	Individual mounting option per request (Special) Requires additional information / specification
ARO2227	10-Way Switch Bracket RH, Mounted for Attendant
ARO2228	10-Way Switch Bracket LH, Mounted for Attendant
ARO2257	Four-Switch Proximity Array (ASL106) NEW!
ARO2275	Centre Mount, swing-away to RH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm)
ARO2276	Centre Mount, swing-away to LH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm)
ARO2215	Individual mounting option per request (Special) - Requires additional information / specification
ARO2243	Proportional Paediatric Compact (ASK132), Includes input Module NEW!
ARO2213	Centre Mount, swing-away to RH, Nucleus
ARO2214	Centre Mount, swing-away to LH, Nucleus
ARO2241	Centre Mount, swing-away to RH, Nucleus (incl Tray (ASL603, 100mm x 150mm)
ARO2242	Centre Mount, swing-away to RH, Nucleus (incl Tray (ASL603, 100mm x 150mm)
ARO2215	Individual mounting option per request (Special) - Requires additional information / specification
ARO2255	Sip'n Puff Head Control (ASL109) NEW! includes Input Module and Lip Switch (ASL 301 LIP) (Combines simple sip and puff controls with head movements)
ARO2256	Sip'n Puff Control NEW! (including Sip'n Puff kit, Lip Switch ASL301LIP and Input module)
ARO2271	Input Module DLX-IN500-A (sub D, Sip'n Puff) NEW!
ARO2253	USB Charger Port DLX-USB02-A
ARO2207	USB Charger Port RH Mount
ARO2208	USB Charger Port LH Mount
ARO2261	Stereo Toggle - direct control for Tilt only
ARO2262	Stereo Button - direct control for Tilt only
ARO2263	Four-Way Toggle - direct control for Tilt & Recline
ARO2264	Four-Way Button - direct control for Tilt & Recline
ARO2207	Toggle / Button RH Mount
ARO2208	Toggle / Button LH Mount
ARO1437	Light system with indicators & LED Headlights (Only available with REM24SD or G90A joystick controls)
ARO2113	Lap belt with auto lock (flexible mount in side rails)
ARO2701	Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ARO2702	Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ARO2703	2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ARO2704	2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ARO2705	2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)
ARO2706	2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ARO1415	Rear view Mirror RH - Must be selected with 10 km/h Motor
ARO1339	Rear view Mirror LH - Must be selected with 10 km/h Motor
ARO1404	Tray Swing Away RH
ARO1405	Tray Swing Away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):