



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

02 DECEMBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm 4 MAX

Standard Features

- Max User Weight 220KG (34.5 Stone)
- Max Seating
- Seat widths from 490mm - 610mm
- Seat depths from 560mm - 630mm
- Black Multi Stretch PU Coated Wipe Down Fabric TR26
- Seat heights from 480mm - 510mm
- Manual Settable Tilt with Manual Settable Recline
- Reinforced Armrests
- Reinforced Legrests with angle adjustable flip-back footboard
- Headrest, sideways and depth adjustable
- LiNX DLX-REM110 with standard bracket
- Choice of 3 colours
- Puncture Proof Castors
- Puncture Proof Drive Wheels
- 6 km/h (4mph) Motors
- 8A (200 - 250V) Charger
- 4 Point tie down system
- 73.5 Ah (Gel) batteries
- Bolted on anti tippers

Standard Base price: €10,817



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
AWA0000	Storm 4 MAX	€10,817
AWZ0000	Special Storm 4 MAX (Extended Lead Time). Please provide details in comments section.	POA
SEATING TYPE*		
AWA1145	MAX SEAT UNIT	◆
BACK CUSHION*		
AWA1106	Tension adjustable 630mm - 620mm high with SW55 / 610mm 570mm High with SW490mm	◆

SEAT WIDTH*			
AWA0378	<input type="checkbox"/>	Seat Width 550mm	○
AWA0380	<input type="checkbox"/>	Seat Width 610mm	○
SEAT CUSHION*			
AWA1459	<input type="checkbox"/>	Seat cushion Max (550mm x 635mm)	○
AWA1460	<input type="checkbox"/>	Seat cushion Max (610mm x 635mm)	○
SEAT DEPTH*			
AWA0396	<input type="checkbox"/>	Seat Depth 560mm	○
AWA0397	<input type="checkbox"/>	Seat Depth 595mm	○
ADAPTER TO REDUCE SEAT DEPTH			
AWA0399	<input type="checkbox"/>	Adapter to shorten the seat depth by 130mm	€59
UPHOLSTERY COLOUR*			
AWA0278	<input type="checkbox"/>	Black Multi Stretch PU Coated Wipe Down Fabric TR26	◆
ARMRESTS INCLUDING REINFORCED DEPTH ADJUSTABLE BRACKET*			
AWA0635	<input type="checkbox"/>	Fixed armrest 250mm - 340mm height adjustable removable	○
ARMPADS*			
AWA0647	<input type="checkbox"/>	Standard Armpad 325mm	○
AWA0649	<input type="checkbox"/>	Armpad extended plus 50mm	○
PUSH BAR*			
AWA1307	<input type="checkbox"/>	Push bar - height adjustable	○
HEADREST / NECKREST (Adequate headrest option must be selected)*			
AWA1432	<input type="checkbox"/>	Headrest sideways and depth adjustable	◆
AWA1427	<input type="checkbox"/>	Headrest with cheek support sideways adjustable	€242
AWA1495	<input type="checkbox"/>	Equivalent non Invacare Headrest supporting the requirements of ISO7176-19:2001. Provided independently from this order	○
SEAT HEIGHT*			
AWA1025	<input type="checkbox"/>	Fixed seat height 480mm	○
AWA1006	<input type="checkbox"/>	Fixed seat height 510mm	○
MANUAL / POWERED SEATING OPTIONS*			
AWA1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	○
AWA1014	<input type="checkbox"/>	Electrical Tilt 0° to 15°	€632
BACKREST RECLINER*			
AWA0958	<input type="checkbox"/>	Settable recliner by spindle (0° - 25°)	○
AWA0947	<input type="checkbox"/>	Electrical recliner (with high pivot, shearing compensation)	€557
LEG REST TYPE*			
AWA0564	<input type="checkbox"/>	Side Mounted: Reinforced hanger fixed angle (not detachable without tools) (NC 0522)	◆
CALF PADS			
AWA0583	<input type="checkbox"/>	Calf Pads adjustable	○
FOOTBOARD*			
AWA0599	<input type="checkbox"/>	Angle adjustable flip back heavy duty footboard	◆
AWA0522	<input type="checkbox"/>	Angle adjustable footboard (for max user weight 200kg) (NC 0564)	○
LINX CONTROLS*			
AWA2201	<input type="checkbox"/>	DLX-REM110 (Drive Only) (NC 1437)	◆
AWA2203	<input type="checkbox"/>	DLX-REM211 (Drive & Actuators) (NC 1437)	○
AWA2204	<input type="checkbox"/>	DLX-REM216 (Drive, Actuators & Lights)	€154

CONTROL BRACKETRY*			
AWA1207	<input type="checkbox"/>	Standard Bracket RH	○
AWA1208	<input type="checkbox"/>	Standard Bracket LH	○
AWA1707	<input type="checkbox"/>	Swing Away RH Bracket	€151
AWA1708	<input type="checkbox"/>	Swing Away LH Bracket	€151
AWA1709	<input type="checkbox"/>	Matrx Resolve Swing in/out Bracket RH	€281
AWA1710	<input type="checkbox"/>	Matrx Resolve Swing in/out Bracket LH	€281
AWA2227	<input type="checkbox"/>	Attendant Mounted RH (NC 2219, 2220)	€68
AWA2228	<input type="checkbox"/>	Attendant Mounted LH (NC 2219, 2220)	€68
DUAL CONTROLS			
AWA2220	<input type="checkbox"/>	LiNX DLX-REM050 (Dual Control used as secondary remote) (NC 2227, 2228, 2245, 1221, 1222)	€466
AWA2219	<input type="checkbox"/>	LiNX DLX-ACU200 (Compact Dual Control) (NC 2227, 2228, 2245, 2201, 1221, 1222)	€444
AWA2245	<input type="checkbox"/>	LiNX DLX-IDC (Intuitive Dual Control) (NC 2227, 2228, 2220, 2219, 2201, 1207, 1208)	€712
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)			
AWA1207	<input type="checkbox"/>	Dual Control Standard Bracket RH (NC 2227, 2228, 2245)	○
AWA1208	<input type="checkbox"/>	Dual Control Standard Bracket LH	○
AWA1221	<input type="checkbox"/>	IDC Dual Control Detachable Bracket RH (NC 2227, 2228, 2220, 2219)	○
AWA1222	<input type="checkbox"/>	IDC Dual Control Detachable Bracket LH (NC 2227, 2228, 2220, 2219)	○
LINX 10-WAY SWITCH			
AWA2251	<input type="checkbox"/>	10-Way Switch DLX-SW10-A for direct Powered Seating Access	€513
LINX 10-WAY SWITCH BRACKET			
AWA2217	<input type="checkbox"/>	Bracket RH, mounted for user	○
AWA2218	<input type="checkbox"/>	Bracket LH, mounted for user	○
AWA2227	<input type="checkbox"/>	Bracket RH, mounted for attendant	€30
AWA2228	<input type="checkbox"/>	Bracket LH, mounted for attendant	€30
LINX ACCESSORY CABLE			
AWA2252	<input type="checkbox"/>	Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	€105
REMOTE STOP SWITCH			
AWZ2258	<input type="checkbox"/>	Remote Stop Switch (ASL504)	€1,019
FURTHER LINX COMPONENTS / ACCESSORIES			
AWA2302	<input type="checkbox"/>	Not Necessary (NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716)	◆
AWA2300	<input type="checkbox"/>	To be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	○
BUDDY BUTTONS			
AWA1685	<input type="checkbox"/>	Red Buddy Button	€86
AWA1686	<input type="checkbox"/>	Yellow Buddy Button	€86
AWA1687	<input type="checkbox"/>	Green Buddy Button	€86
AWA1690	<input type="checkbox"/>	Blue Buddy Button	€86
EGG SWITCH			
AWA2254	<input type="checkbox"/>	Egg Switch, black	€86

PIKO BUTTONS			
AWA1711		Red Piko Button	€153
AWA1712		Yellow Piko Button	€153
AWA1713		Green Piko Button	€153
AWA1714		Orange Piko Button	€153
AWA1715		Blue Piko Button	€153
AWA1716		Black Piko Button	€153
TOGGLES / BUTTONS			
AWZ2261		Stereo Toggle - direct control for Tilt only	€250
AWZ2262		Stereo Button - direct control for Tilt only	€250
AWZ2263		Four-Way Toggle - direct control for Tilt & Recline	€326
AWZ2264		Four-Way Button - direct control for Tilt & Recline	€348
AWZ2207		Toggle / Button RH Mount	○
AWZ2208		Toggle / Button LH Mount	○
AWZ2225		Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
AWA2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	€71
AWA2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	€71
AWA2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	€71
AWA2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	€71
SHROUD COLOUR*			
AWA0202		Blue	○
AWA0228		Black	○
AWA0204		Transparent Cherry Red	○
AWZ0299		Special Shroud Colour (RAL Code must be supplied)	POA
CASTORS*			
AWA0840		10" Grey Pneumatic	○
AWA0858		10" Grey Puncture Proof Sealed	○
AWA0842		10" Grey Puncture Proof Solid Inlay	◆
DRIVE WHEELS*			
AWA0771		14" Grey Pneumatic	○
AWA0773		14" Grey Puncture Proof Sealed	○
AWA0750		14" Grey Puncture Proof Solid Inlay	◆
MOTOR TYPE*			
AWA1277		7 km/h	◆
BATTERY*			
AWA1364		73.5 AH Gel	◆
SUSPENSION*			
AWA1203		Heavy Duty Suspension	◆
G TRAC			
AWA1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	€865

CHARGER*		
AWA1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V) ◆
AWA1313	<input type="checkbox"/>	OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster €166
LIGHTS*		
AWA1254	<input type="checkbox"/>	No Lights ◆
AWA1437	<input type="checkbox"/>	LED Light system with indicators (30 LUX) €479
TRANSPORTATION KIT*		
AWA1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. ◆
FIXATION BELT*		
AWA1455	<input type="checkbox"/>	Fixation Belt Extra long with Auto Lock ◆
CANE HOLDER		
AWA1312	<input type="checkbox"/>	Cane holder RH €41
AWA1363	<input type="checkbox"/>	Cane holder LH €41
WHEEL FENDERS		
AWA1417	<input type="checkbox"/>	Fender for Drive Wheels €128

TOTAL PRICE	
--------------------	--

Non Conforming Options:

AWA0564	Side Mounted: Reinforced hanger fixed angle (not detachable without tools)
AWA0522	Angle adjustable footboard (for max user weight 200kg)
AWA2201	DLX-REM110 (Drive Only)
AWA1207	Standard Bracket RH
AWA1208	Standard Bracket LH
AWA2227	Attendant Mounted RH
AWA2228	Attendant Mounted LH
AWA2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
AWA2219	LiNX DLX-ACU200 (Compact Dual Control)
AWA2245	LiNX DLX-IDC (Intuitive Dual Control)
AWA1207	Dual Control Standard Bracket RH
AWA1208	Dual Control Standard Bracket LH
AWA1221	IDC Dual Control Detachable Bracket RH
AWA1222	IDC Dual Control Detachable Bracket LH
AWA2227	Bracket RH, mounted for attendant
AWA2228	Bracket LH, mounted for attendant
AWA1686	Yellow Buddy Button
AWA1687	Green Buddy Button
AWA1690	Blue Buddy Button
AWA2254	Egg Switch, black
AWA1711	Red Piko Button
AWA1712	Yellow Piko Button
AWA1713	Green Piko Button
AWA1714	Orange Piko Button
AWA1715	Blue Piko Button
AWA1716	Black Piko Button
AWA1437	LED Light system with indicators (30 LUX)



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):