



Tel: 01 810 7084 Fax: 01 810 7085
Online Spares available at: www.invacare.ie
Email: ordersireland@invacare.com
Jan-20

Customer Ref:	Account No:
Contact No:	Order No:
Delivery Address:	Order Qty:
	Order Date:

Invacare® Topaz
RETAIL
Max user weight 318Kg (50 stone)

Standard Features



- Choice of 6 Seat Widths 500mm - 760mm (20" - 30")
- Choice of 2 Seat Depths 460mm or 500mm (18" or 20")
- Choice of 3 Seat to Ground Heights 395mm - 495mm (15.5" - 19.5")
- Full length, fixed height armrest
- Adjustable height, 10° backrest 405mm (16") 430mm (17") or 460mm (18")
- Wet Black Frame Colour
- 200mm x 40mm Semi Pneumatic Front Castors
- 610mm HD Spoked Solid Rear Wheels with Chrome Handrim
- Swing away, detachable legrests with extra large footplate
- Heel loops
- Pull to lock brakes



Standard Base Price : €1,998

* Configurations may vary to the image shown

Code	Description (NB: All sizes are approximate)	NC's (Non Conformance)	Retail
<p>X Please make your selection with a cross in the box ♦ Standard feature ○ Free of charge (NC:9999) Non-compatibilities</p> <p>Customer weight = Kg</p>			
HTZ0001	Topaz		€1,998
Topaz Special Basic Model - PLEASE COMPLETE DETAILS AT THE END OF THE FORM			
SEAT			
SEAT WIDTH			
HTZ0002	500mm (20")		○
HTZ0003	550mm (22")		○
HTZ0004	610mm (24")		○
HTZ0005	660mm (26")		○
HTZ0037	710mm (28") (Extended lead time)		€362
HTZ0038	760mm (30") (Extended lead time)		€362
SEAT DEPTH			
HTZ0006	460mm (18")		○
HTZ0007	500mm (20")		○
SEAT HEIGHT			
HTZ0008	19½" Adult		○
HTZ0009	17½" Hemi		○
ARMREST TYPE			
HTZ0011	Full length, fixed height, conventional		♦
HTZ0012	Full length, adjustable height, conventional		€177
BACKREST TYPE			
HTZ0013	Adjustable height backrest 16", 17" or 18" with 10° Bend		♦
COLOURS			
HTZ0014	Wet Black		♦

FRONT WHEELS			
HTZ0016		8" x 1.75" Semi Pneumatic (black)	◆
REAR WHEELS			
HTZ0036		24" HD Spoked wheel solid tyre	◆
HANDRIM OPTION			
HTZ0019		Aluminium	○
UPHOLSTERY TYPE			
HTZ0020	X	Black Vinyl	◆
HANGERS/FOOTPLATES			
HTZ0022		One Piece Footboard (NC: 0026)	○
HTZ0023		Swing away footplates (with extra large footplates)	◆
HTZ0024		Detachable Swingaway Elevating (with extra large footplates)	€129
LEGREST ACCESSORY			
HTZ0025		Calf strap	◆
HTZ0026	X	Heel loops	◆
BACKREST			
HTZ0028		Push bar	€210
SEAT ACCESSORY			
HTZ0029		Solid seat insert	€63
ACCESSORIES			
HTZ0030		Oxygen holder package assembly	€37
HTZ0031		Rear anti tippers	€30
HTZ0032		Pull to lock brake	◆
HTZ0033		Crutch and cane holder	€28
HTZ0034		Brake extensions	€10
HTZ0035		Non folding device	€22
			Total Price:

Seat-to-Floor Matrix 8" Casters			
Seat to Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17.5"	Top	24"	Top
19.5"	Bottom	24"	Bottom



Please refer to our Invacare Unique special order form, for more complex custom requirements. Available at www.invacare.co.uk

Special Instructions - Please supply the following information:

Name:

Email:

Phone:

Fax:

Description of special request (please include relevant option code or part number):