



Yes, you can:

# Powerchairs Ireland Prescription Form

NOVEMBER 2023  
Customer Ref:

Please refer to [www.invacare.ie](http://www.invacare.ie) for Local Distributor

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## AVIVA RX40 HD with Modulite HD seating

### Standard Features

- Max user weight 160kg (25 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite HD adjustable seating system
- HD Telescopic Seat Frame and Plate with Side Rails width 530 to 550 mm with an additional 50mm, (2") by arms
- Seat depths 510mm - 580mm stepless adjustment
- Tension adjustable backrest
- Standard T-Armrests, detachable with short armpad
- Manual settable seat Tilt (0° - 20°)
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM110 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h (4mph) 4 pole Durawatt motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: €7,720**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =  KG

Code	Description	
<b>BASE MODEL*</b>		
JBH0002	AVIVA RX 40 HD Base Model Max user weight up to 160 kg HD Version, includes RX 40 Performance Package Wide Base only	€7,720
JBX1300	Special Extended lead time: + 10 days	POA
<b>CHASSIS WIDTH*</b>		
JBH2830	Chassis Wide (635 mm/25")	<input checked="" type="checkbox"/>

<b>SEATING TYPE*</b>			
JBX2101	<input type="checkbox"/>	Invacare® Modulite HD Seat Unit On RX 40 HD only Tension Adjustable Upholstry only JBX2171	<input type="radio"/>
<b>SEAT FRAME*</b>			
JBX2102	<input type="checkbox"/>	HD Telescopic Seat Frame and Plate with Side Rails width 530 to 550 mm	<input type="radio"/>
<b>SEAT WIDTH*</b>			
JBX2124	<input type="checkbox"/>	530 mm Seat Plate (NC 2702, 2704, 2706)	<input type="radio"/>
JBX2104	<input type="checkbox"/>	550 mm Seat Plate	<input type="radio"/>
<b>SEAT DEPTH*</b>			
JBX2136	<input type="checkbox"/>	Seat Depth 510 mm	<input type="radio"/>
JBX2137	<input type="checkbox"/>	580 mm	<input type="radio"/>
<b>SEAT CUSHION*</b>			
JBX2148	<input type="checkbox"/>	Seat Cushion (90 mm high)	<input type="radio"/>
JBX2140	<input type="checkbox"/>	Without seat cushion	<input type="radio"/>
<b>BACKREST SUPPORT*</b>			
JBX2162	<input type="checkbox"/>	Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrix Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	<input checked="" type="radio"/>
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
JBX2171	<input type="checkbox"/>	Tension Adjustable Upholstery (NC 2230, 2247) nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm MUST BE SELECTED WITH HD SEATING JBX2102	<input checked="" type="radio"/>
<b>UPHOLSTERY TYPE*</b>			
JBX2156	<input type="checkbox"/>	Black Nylon	<input checked="" type="radio"/>
<b>ARMREST (Pair)*</b>			
JBX2191	<input type="checkbox"/>	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	<input type="radio"/>
JBX2181	<input type="checkbox"/>	T-Arm detachable, Extended height Adjustment range from seat plate: 295 - 365 mm	€50
<b>ARMPAD (Pair)*</b>			
JBX2195	<input type="checkbox"/>	Short Armpad (325 mm)	<input checked="" type="radio"/>
JBX2196	<input type="checkbox"/>	Long Armpad (375 mm)	<input type="radio"/>
JBX2197	<input type="checkbox"/>	Long Ergonomic Armpad (397 mm) (NC 2230, 2247)	<input type="radio"/>
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
JBX1350	<input type="checkbox"/>	Neckrest depth adjustable (NC 2230, 2247)	€103
JBX1432	<input type="checkbox"/>	Headrest sideways adjustable (NC 2230, 2247)	€285
JBX1427	<input type="checkbox"/>	Headrest with cheek support (NC 2230, 2247)	€285
JBX2324	<input type="checkbox"/>	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware (NC 2230, 2247)	€454
JBX2325	<input type="checkbox"/>	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware (NC 2230, 2247)	€454
JBX2327	<input type="checkbox"/>	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware (NC 2230, 2247)	€615
JBX1495	<input type="checkbox"/>	Non Invacare Headrest - no mounting hardware With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	<input type="radio"/>

SEAT LIFTER / SEAT TILT*			
JBX1091		Manually settable Seat Tilt 0 - 15° (NC 1401)	○
JBX1092		Powered Tilt (0° - 20°) (NC 2201)	€507
JBX1093		Powered Tilt (0° - 30°) with CoG shift (NC 2201)	€1,348
JBX1094		Seat lifter with powered tilt (300 mm, 30°) with CoG shift, 6 point Tie downs Must be selected (NC 2201, 2220)	€2,377
SEAT TO GROUND HEIGHT*			
JBX2857		Seat-to-Floor Height 435 mm	○
JBX2858		Seat-to-Floor Height 460 mm	○
JBX2859		Seat-to-Floor Height 485 mm	○
BACKREST RECLINER*			
JBX0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
JBX0977		Powered Recliner 30° (NC 2201)	€634
CENTRE MOUNTED LEGREST*			
JBH0556		Centre mount legrest, angle settable Length: 280 - 385 mm (NC 0517, 1244, 1416)	€468
JBH0589		LNx Power centre mount legrest, length compensating Length: 340 - 410 mm (NC 2201, 1244, 1416)	€2,060
JBH2855		Without Legrests	○
CENTRE MOUNT FOOTPLATES / FOOT PLATFORM*			
JBH2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	€130
JBH0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBH0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")	€170
JBH0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	€339
JBX3001		Please select If Centre mount footplate is not required. (NC 0589)	○
FOOTREST HANGER TYPE*			
JBH0503		Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools	◆
JBH0533		Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad	€292
JBH0531		Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad (NC 2201)	€841
JBH0501		Without Hanger pair (NC 0504, 0513)	○
FOOTPLATE TYPE*			
JBH0512		Fixed Footplate pair (NC 0533, 0531)	◆
JBH0504		Angle Adjustable Footplate pair	€40
JBH0513		Angle and depth adjustable Footplate pair	€149
JBX3003		No Side Mounted Footplate required (NC 0503, 0533, 0531)	○
CONTROLS*			
JBH2201		Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module (NC 2837, 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2258, 2219, 2245, 2270, 1908)	◆
JBH2203		Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators (NC 1437, 2231, 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2270, 2258, 1908, 1909)	○
JBH2204		Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights (NC 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2270, 2258)	€191
JBH2205		Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets (NC 2220)	€428
JBH2206		Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233 (NC 2220, 2270)	€974

<b>REM400 TOGGLE OPTION</b>			
JBH2270		Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)	(NC 2203, 2204, 2206) €397
<b>CONTROL BRACKETRY*</b>			
JBH1207		Standard Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) ○
JBH1208		Standard Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) ○
JBH1707		Swing Away Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) €167
JBH1708		Swing Away Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) €167
JBH2213		Centre Mount, swing-away to RH, Nucleus	(NC 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 1405) €626
JBH2214		Centre Mount, swing-away to LH, Nucleus	(NC 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1404, 1405) €626
JBH2227		Bracket RH, mounted for Attendant	(NC 1287, 2220, 2128, 2219, 1221, 1312) €81
JBH2228		Bracket LH, mounted for Attendant	(NC 1286, 2220, 2128, 2219, 1222, 1363) €81
JBX2239		Individual mounting option per request (Special) Requires additional information/specification	POA
<b>DUAL CONTROL</b>			
JBH2220		Dual Control LiNX (DLX-REM050 as secondary remote)	(NC 2273, 2205, 2206, 1221, 1222, 2227, 2228) €497
JBH2219		Compact Dual Control LiNX (DLX-ACU200)	(NC 1221, 1222, 2227, 2228) €497
<b>INVACARE LINX IDC DUAL CONTROL</b>			
JBH2245		Invacare LiNX IDC Dual control	(NC 2259, 2260) €943
<b>DUAL CONTROL BRACKETRY</b>			
JBH2259		Bracket RH, Attendant	(NC 2291, 2245, 1312) ○
JBH2260		Bracket LH, Attendant	(NC 2292, 2245, 1363) ○
<b>INVACARE LINX IDC DUAL CONTROL BRACKETRY</b>			
JBH1221		Bracket RH, detachable	(NC 2291, 2220, 2219, 2227, 1312) ○
JBH1222		Bracket LH, detachable	(NC 2292, 2220, 2219, 2228, 1363) ○
<b>LINX COMPACT REMOTE</b>			
JBH2231		Proportional Compact Remote DLX-CR400	(NC 2203, 2204) €903
JBH2232		Proportional Compact Remote Low Force DLX-CR-400LF	(NC 2203, 2204) €1,129
<b>LINX COMPACT REMOTE BRACKET</b>			
JBH2207		RH Mount	(NC 1707, 2213, 2217, 1404) ○
JBH2208		LH Mount	(NC 1708, 2214, 2218, 1405) ○
JBH2265		Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1707, 2213, 2217, 1404, 1405, 1415) €626
JBH2266		Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1708, 2214, 2218, 1404, 1405, 1339) €626
JBX2215		Individual mounting option per request (Special) Requires additional information/specification	POA
<b>CHIN CONTROL</b>			
JBH2229		Proportional Chin Joystick (ASL138) SELECTED HEADREST ALWAYS WITH STD MOUNTING HARDWARE	(NC 2201, 2203, 2204, 2230, 2247, 2233, 2257, 2240, 2128) €2,178
<b>CHIN CONTROL BRACKETRY MANUAL</b>			
JBH1289		Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427) ○
JBH1288		Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427) ○
<b>CHIN CONTROL BRACKETRY POWERED</b>			
JBH1287		Swing Away Bracket RH, powered	(NC 2259, 1221, 2291, 1312, 1415) €2,743
JBH1286		Swing Away Bracket LH, powered	(NC 2260, 1222, 2292, 1363, 1339) €2,743

LINX HEAD CONTROL			
JBH2230		Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3	(NC 2203, 2204, 2229, 2233, 2324, 2325, 2327, 2257, 2240, 2171, 2197, 1350, 1432, 1427) €4,673
JBH2247		Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings	(NC 2203, 2204, 2229, 2233, 2324, 2325, 2327, 2257, 2240, 2171, 2197, 1350, 1432, 1427) €5,472
LINX FOUR -SWITCH PROXIMITY ARRAY			
JBH2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2229, 2230, 2247, 2233, 2240, 2128) €7,236
JBH2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1707, 2213, 2217, 1404, 1405, 1415) €626
JBH2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1708, 2214, 2218, 1404, 1405, 1339) €626
JBX2269		Individual mounting option per request (Special) Requires additional information/specification	POA
LINX WIRELESS ACCESSORIES FOR HEAD CONTROL AND PROXIMITY CONTROL			
JBH2249		Wireless Triple Switch Receiver (ASL557-3)	€428
JBH2250		Wireless Mouse Emulator (ASL558)	€1,027
LINX MICRO EXTREMITY CONTROL BRACKET			
JBH2233		MEC Proportional Joystick (Micro Extremity Control ASL130)	(NC 2203, 2204, 2229, 2230, 2247, 2257, 2128) €4,537
LINX HEAVY DUTY REMOTE			
JBH2240		Proportional Heavy Duty Remote Joystick w mode switch (ASL133)	(NC 2203, 2204, 2229, 2230, 2247, 2257, 2128) €2,529
LINX MICRO EXTREMITY AND HEAVY DUTY CONTROL BRACKET			
JBH2285		Tray Mount RH (100x150 mm)	(NC 1707, 2213, 2217, 1404, 1405, 1415) ○
JBH2286		Tray Mount LH (100x150 mm)	(NC 1708, 2214, 2218, 1404, 1405, 1339) ○
JBH2283		Centre Mount, swing-away to RH, Nucleus	(NC 1707, 2213, 2217, 1404, 1405, 1415) €626
JBH2284		Centre Mount, swing-away to LH, Nucleus	(NC 1708, 2214, 2218, 1404, 1405, 1339) €626
JBH2287		Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1707, 2213, 2217, 1404, 1405, 1415) €626
JBH2288		Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1708, 2214, 2218, 1404, 1405, 1339) €626
JBX2289		Individual mounting option per request (Special) Requires additional information/specification	POA
LINX 10-WAY SWITCH			
JBH2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	€634
LINX 10-WAY SWITCH BRACKET			
JBH2217		Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2265, 1404) ○
JBH2218		Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2266, 1405) ○
JBH2291		Bracket RH, mounted for attendant	(NC 2259, 1221, 1287) €81
JBH2292		Bracket LH, mounted for attendant	(NC 2260, 1222, 1288) €81
LINX USB CHARGER PORT			
JBH2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706) €207
LINX USB CHARGER PORT BRACKET			
JBH2281		RH Mount	○
JBH2282		LH Mount	○
REMOTE STOP SWITCH			
JBH2258		Remote Stop Switch (ASL504)	(NC 2203, 2204) €1,319
LINX IN200 MODULE			
JBH2273		IN200 Module (Jack port for external power button + control inputs via DB9)	€303
FURTHER LINX COMPONENTS / ACCESSORIES			
JBH2300		To be enclosed only	○
JBX2301		To be connected and programmed (Special)	POA

<b>BUDDY BUTTONS</b>			
JBH1685		Red Buddy Button	€110
JBH1686		Yellow Buddy Button	€110
JBH1687		Green Buddy Button	€110
JBH1690		Blue Buddy Button	€110
<b>EGG SWITCH</b>			
JBH2254		Egg Switch, black	€110
<b>PIKO BUTTONS</b>			
JBH1711		Red Piko Button	€198
JBH1712		Yellow Piko Button	€198
JBH1713		Green Piko Button	€198
JBH1714		Orange Piko Button	€198
JBH1715		Blue Piko Button	€198
JBH1716		Black Piko Button	€198
<b>BODYPOINT JOYSTICK KNOBS</b>			
JBH2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	€79
JBH2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	€79
JBH2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	€79
JBH2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	€79
<b>5:1 SWITCH ADAPTOR</b>			
JBH1681		5:1 Switch Adapter ( 9-pin sub D to 5 jack ports)	€152
<b>SHROUD COLOUR*</b>			
JBH2871		Mid range - Black/Anthracite Shroud	○
JBH2872		Mid range - Silver Shroud	○
JBH2873		Mid range - Bengal Red Shroud	○
JBH2874		Mid range - Mercury Blue Shroud	○
JBH2875		Mid range - Midnight Purple Shroud	○
JBX0299		Premium - Individual Shroud Colour (includes coloured spreader bars on Modulite back) RAL number required Extended lead lime: + 15 days	€269
JBH2865		Premium - Racing Green Shroud	€131
JBH2866		Premium - Deep Ocean Blue Satin Shroud	€131
JBH2867		Premium - Frosty Blue Shroud	€131
JBH2868		Premium - Ultimate Red Shroud	€131
JBH2869		Premium - Falcon White Shroud	€131

<b>RIM COLOUR*</b>			
JBH2882		Mid range - Black Rim Inserts	○
JBH2883		Mid range - Silver Rim Inserts	○
JBH2884		Mid range - Bengal Red Rim Inserts	○
JBH2885		Mid range - Mercury Blue Rim Inserts	○
JBH2886		Mid range - Midnight Purple Rim Inserts	○
JBX2887		Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days	○
JBH2888		Premium - Racing Green Rim Inserts	○
JBH2889		Premium - Deep Ocean Blue Satin Rim Inserts	○
JBH2890		Premium - Frosty Blue Rim Inserts	○
JBH2891		Premium - Ultimate Red Rim Inserts	○
JBH2892		Premium - Falcon White Rim Inserts	○
<b>FORKS*</b>			
JBH2631		Single sided fork	€285
JBH2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
<b>FRONT WHEEL*</b>			
JBH0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	○
JBH0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆
<b>CASTOR LOCKS</b>			
JBH1490		Castor Lock	€248
<b>DRIVE WHEEL*</b>			
JBH0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
JBH0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
<b>C.T.C SUSPENSION*</b>			
JBH2870		CTC Suspension Control - Traction - Comfort	◆
<b>SUSPENSION SETTING*</b>			
JBH1916		Suspension up to 160 kg User Weight For RX 40 HD only	◆
<b>MOTOR TYPE*</b>			
JBH1294		Indoor setting, reduced speed approx. 3 km/h	○
JBH1282		6 km/h	◆
JBH1908		10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2203, 1437, 1254, 2837) €683
JBH1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2201, 2203, 1437, 1254, 2837) €820
<b>BATTERY TYPE*</b>			
JBH1259		60 Ah Gel	◆
JBH1260		73,5 Ah Gel	(NC 1450) €116
<b>G TRAC</b>			
JBH1215		G-Trac Module	€406
<b>CHARGER*</b>			
JBH1343		Opti Charge Charger 8A (IP44)	◆
JBH1313		Opti Charge Charger 10A (IP44)	€191
<b>ANTI-TIPPERS*</b>			
JBH1205		Anti-tipper	◆

LIGHTS				
JBH1437		LED Light system with Indicators 30 Lux headlights	(NC 2203, 1908, 1909)	€602
JBH2837		LED Light system extra bright with indicators 60 Lux headlights	(NC 2203)	€180
JBH1254		Without Light system	(NC 1909)	○
KERB RIDER				
JBH1244		Kerb Rider With park position feature	(NC 0556, 0589, 1294)	€266
TRANSPORTATION KIT*				
JBH1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering		€736
JBH1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory, Not compatible with Seat lifter		◆
JBH1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory	(NC 1091, 1092, 1093)	○
LAP BELT*				
JBH1329		Lap belt with Auto Lock Standard mount on backrest		○
JBH2113		Lap belt with Auto Lock Flexible mount in side rails		○
POSITIONING BELT				
JBH2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2253)	€97
JBH2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2124, 2253)	€97
JBH2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2253)	€100
JBH2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2124, 2253)	€100
JBH2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2253)	€100
JBH2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2124, 2253)	€100
REAR VIEW MIRROR				
JBH1415		Rear view mirror RH	(NC 2217, 2285, 2213, 2265, 2275, 2283, 2287)	€53
JBH1339		Rear view mirror LH	(NC 2218, 2286, 2214, 2266, 2276, 2284, 2288)	€53
OTHER OPTIONS				
JBH2128		Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight	(NC 2227, 2228, 2229, 2233, 2240, 2257, 2291, 2292)	€97
JBH1312		Cane holder RH	(NC 2227, 2259, 1221)	€53
JBH1363		Cane holder LH	(NC 2228, 2260, 1222)	€53
JBH1416		Fender for Castors		€12
JBH1417		Fender for Drive wheels		€213
JBH1404		Tray, Swing Away RH	(NC 1707, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2207, 2213, 2214, 2217)	€210
JBH1405		Tray, Swing Away LH	(NC 1708, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2208, 2213, 2214, 2218)	€210

<b>TOTAL PRICE</b>	
--------------------	--



**Non Conforming Options:**

JBX2124	530 mm Seat Plate
JBX2171	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm MUST BE SELECTED WITH HD SEATING JBX2102
JBX2197	Long Ergonomic Armpad (397 mm)
JBX1350	Neckrest depth adjustable
JBX1432	Headrest sideways adjustable
JBX1427	Headrest with cheek support
JBX2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware
JBX2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware
JBX2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware
JBX1091	Manually settable Seat Tilt 0 - 15°
JBX1092	Powered Tilt (0° - 20°)
JBX1093	Powered Tilt (0° - 30°) with CoG shift
JBX1094	Seat lifter with powered tilt (300 mm, 30°) with CoG shift, 6 point Tie downs Must be selected
JBX0977	Powered Recliner 30°
JBH0556	Centre mount legrest, angle settable Length: 280 - 385 mm
JBH0589	LNx Power centre mount legrest, length compensating Length: 340 - 410 mm
JBH0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBH0503	Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools
JBH0533	Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad
JBH0531	Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad
JBH0504	Angle Adjustable Footplate pair
JBH0513	Angle and depth adjustable Footplate pair
JBH2201	Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module
JBH2203	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators
JBH2204	Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights
JBH2205	Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets
JBH2206	Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233
JBH2270	Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)
JBH1707	Swing Away Bracket RH
JBH1708	Swing Away Bracket LH
JBH2213	Centre Mount, swing-away to RH, Nucleus
JBH2214	Centre Mount, swing-away to LH, Nucleus
JBH2227	Bracket RH, mounted for Attendant
JBH2228	Bracket LH, mounted for Attendant
JBH2220	Dual Control LiNX (DLX-REM050 as secondary remote)
JBH2219	Compact Dual Control LiNX (DLX-ACU200)
JBH2245	Invacare LiNX IDC Dual control
JBH2259	Bracket RH, Attendant

JBH2260	Bracket LH, Attendant
JBH1221	Bracket RH, detachable
JBH1222	Bracket LH, detachable
JBH2231	Proportional Compact Remote DLX-CR400
JBH2232	Proportional Compact Remote Low Force DLX-CR-400LF
JBH2207	RH Mount
JBH2208	LH Mount
JBH2265	Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBH2266	Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBH2229	Proportional Chin Joystick (ASL138) SELECTED HEADREST ALWAYS WITH STD MOUNTING HARDWARE
JBH1288	Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest
JBH1287	Swing Away Bracket RH, powered
JBH1286	Swing Away Bracket LH, powered
JBH2230	Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3
JBH2247	Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings
JBH2257	Four-Switch Proximity Array (ASL106)
JBH2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2233	MEC Proportional Joystick (Micro Extremity Control ASL130)
JBH2240	Proportional Heavy Duty Remote Joystick w mode switch (ASL133)
JBH2285	Tray Mount RH (100x150 mm)
JBH2286	Tray Mount LH (100x150 mm)
JBH2283	Centre Mount, swing-away to RH, Nucleus
JBH2284	Centre Mount, swing-away to LH, Nucleus
JBH2287	Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm
JBH2288	Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm
JBH2217	Bracket RH, mounted for user
JBH2218	Bracket LH, mounted for user
JBH2291	Bracket RH, mounted for attendant
JBH2292	Bracket LH, mounted for attendant
JBH2253	USB Charger Port DLX-USB02-A
JBH2258	Remote Stop Switch (ASL504)
JBH2273	IN200 Module (Jack port for external power button + control inputs via DB9)
JBH1294	Indoor setting, reduced speed approx. 3 km/h
JBH1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1437	LED Light system with Indicators 30 Lux headlights
JBH2837	LED Light system extra bright with indicators 60 Lux headlights
JBH1254	Without Light system
JBH1244	Kerb Rider With park position feature
JBH1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering

JBH1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory
JBH2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JBH2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JBH2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JBH2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JBH2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JBH1415	Rear view mirror RH
JBH1339	Rear view mirror LH
JBH2128	Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight
JBH1312	Cane holder RH
JBH1363	Cane holder LH
JBH1416	Fender for Castors
JBH1404	Tray, Swing Away RH
JBH1405	Tray, Swing Away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):