



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Kuschall The KSL 2.0

Standard Features

- Max User Weight 100kg (16 stone)
- Active, Fully Welded, Lightweight, Rigid Wheelchair with Height Adjustable Rigid Footrest
- Fully Hydroformed Aluminium Front Frame 75° or 90° Front Angle
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims & Pneumatic Tyres
- Choice of 6 Seat Widths 340mm to 440mm (13" - 17")
- Choice of 5 Seat Depths 385mm to 485mm (15" - 19")
- Range of Front Seat Heights from 450mm to 530mm (17" - 21")
- Range of Rear Seat Heights from 400mm to 480mm (15" - 19")
- Lockdown Folding, Angle & Height Adjustable Backrest
- Black Nylon Upholstery
- Without Push Handles, Carbon Clothesguards
- Carbon Fibre Axle Tube with 3° Camber
- Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit

Standard Base price: €4,057



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
DSA0085	<input type="checkbox"/> KÜSCHALL THE KSL 2.0	€4,057
UNIQUE SPECIALS (For special modifications only) Additional lead time of 4 weeks		
DSA0000	<input type="checkbox"/> UNIQUE Special Request - For use when requesting any special options or builds	POA
UNIQUE CUSTOMISATION AND PERSONALISATION OPTIONS Additional 4 weeks lead time for Unique options		
USA0001	<input type="checkbox"/> UNIQUE Personalisation - For use when selecting USA customisation and personalisation options	<input type="radio"/>

FRAME*			
DSA1000		Aluminium Hydroformed V Style Frame Angle 75° to the ground	◆
DSA1010		Aluminium Hydroformed Tapered/Adducted Frame Angle 75° to the ground	○
DSA1020		Aluminium Hydroformed V Style Frame Angle 90° to the ground	○
DSA1030		Aluminium Hydroformed Tapered/Adducted Frame Angle 90° to the ground	○
		(NC 5300, 5350, 5400, 5450, 5500)	
WELDED UNDERSEAT RIGID BAR*			
SWO3333		Under seat rigidiser bar (Welded)	€128
SWO5555		Without under seat rigidiser bar	◆
		(NC 4444, 9999)	
SEAT WIDTH*			
DSA1315		340 mm	○
		(NC 2710, 2720, 2750, 2748)	
DSA1320		360 mm	○
		(NC 2748)	
DSA1325		380 mm	○
DSA1330		400 mm	◆
DSA1335		420 mm	○
DSA1340		440 mm	○
SEAT DEPTH*			
DSA1430		385 mm	○
		(NC 2720, 7810, 7820, 7830)	
DSA1440		410 mm	◆
DSA1455		435 mm	○
DSA1470		460 mm	○
DSA1485		485 mm	○
		(NC 2710, 2748)	
BACK HEIGHT*			
DSA1605		270 mm	○
		(NC 2240)	
DSA1610		285 mm	○
		(NC 2240)	
DSA1615		300 mm	◆
		(NC 2240)	
DSA1620		315 mm	○
		(NC 2240)	
DSA1625		330 mm	○
		(NC 2240)	
DSA1630		345 mm	○
		(NC 2240)	
DSA1635		360 mm	○
		(NC 2240)	
DSA1640		375 mm	○
		(NC 2240)	
DSA1645		390 mm	○
		(NC 2240)	
DSA1650		405 mm	○
		(NC 2240)	
DSA1655		420 mm	○
		(NC 2240)	
FRONT SEAT HEIGHT. Lower leg must be front height -100mm To -30mm on 75° frames and -80mm To -40mm on 90° frames*			
DSA1740		450 mm	○
		(NC 1840, 1845, 1850, 1855, 1860, 6420)	
DSA1745		460 mm	○
		(NC 1845, 1850, 1855, 1860)	
DSA1750		470 mm	○
		(NC 1850, 1855, 1860)	
DSA1755		480 mm	○
		(NC 1855, 1860)	
DSA1760		490 mm	○
		(NC 1860)	
DSA1765		500 mm	◆
DSA1770		510 mm	○
		(NC 1820)	
DSA1775		520 mm	○
		(NC 1820, 1825)	
DSA1780		530 mm	○
		(NC 1820, 1825, 1830, 6410)	

REAR SEAT HEIGHT. Rear height = "Front height -100 MM" TO "Front height -20 MM"*			
DSA1820		400 mm	(NC 6420) ○
DSA1825		410 mm	○
DSA1830		420 mm	○
DSA1835		430 mm	◆
DSA1840		440 mm	○
DSA1845		450 mm	○
DSA1850		460 mm	○
DSA1855		470 mm	○
DSA1860		480 mm	○
LOWER LEG LENGTH Seat canvas to footplate measurement, cushion height must be deducted.*			
DSA1930		350 mm	○
DSA1932		360 mm	○
DSA1934		370 mm	○
DSA1936		380 mm	○
DSA1938		390 mm	○
DSA1940		400 mm	○
DSA1942		410 mm	○
DSA1944		420 mm	○
DSA1946		430 mm	◆
DSA1948		440 mm	○
DSA1950		450 mm	○
DSA1952		460 mm	○
DSA1954		470 mm	○
DSA1956		480 mm	○
DSA1958		490 mm	○
DSA1960		500 mm	○
BACKREST TYPE*			
DSA2110		Fixed Aluminium Hydroformed backrest in frame colour	(NC 6666, 7777) ◆
SWO4444		K-Series style lockdown folding backrest, supplied with carbon fibre flip up, fold in clothesguards and tension adjustable upholstery	(NC 2200, 2220, 2240, 3410, 3430) €491
SWO9999		K-Series style lockdown folding backrest, supplied with carbon fibre removable mudguards and tension adjustable upholstery	(NC 2200, 2220, 2240, 3410, 3430) €491
BACKREST UPHOLSTERY*			
DSA2200		Standard Light Black Nylon upholstery (No tension adjustment)	(NC 0002, 0003, 0004, 2241, 2242, 2243, 2244, 2245) ◆
DSA2220		Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener	(NC 0002, 0003, 0004, 2241, 2242, 2243, 2244, 2245) €160
DSA2230		Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener	€253
DSA2240		Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required	(NC 0002, 0003, 0004, 2241, 2242, 2243, 2244, 2245) ○
SWO7777		Nylon tension adjustable upholstery for options 4444 & 9999	(NC 0002, 0003, 0004, 2241, 2242, 2243, 2244, 2245) ○

BACK UPHOLSTERY STITCHING COLOUR			
USA2241		Neon Yellow	€96
USA2242		Neon Pink	€96
USA2243		Blue	€96
USA2244		White	€96
USA2245		Red	€96
USA0002		Your personal colour stitching. Add RAL colour of your choice in special comments box, for a personalised colour stitching on back upholstery.	€96
BACK UPHOLSTERY EMBROIDERY			
USA0003		Your personal design embroidery. Supply your design with your order, subject to copyright. Design size up to 250mm x 250mm (NC 0004)	€256
BACK UPHOLSTERY FLOCKING			
USA0004		Your personal flock design. Supply your design with your order, subject to copyright. Design size up to 280mm x 400mm	€127
PUSH HANDLES*			
DSA2400		Without Push Handles	◆
DSA2410		Push Handles Fixed Short/Mini Hand Grip	€149
DSA2420		Push Handles Fixed Long Hand Grip	€149
DSA2440		Push Handles Height Adjustable, Integrated in backrest	€225
DSA2450		Push Handles Foldable	€287
BACKREST ANGLE*			
DSA2500		74° angle from the seat	○
DSA2520		78° angle from the seat	○
DSA2540		82° angle from the seat	◆
DSA2560		86° angle from the seat	○
DSA2580		90° angle from the seat	○
SEAT UPHOLSTERY*			
DSA2600		Nylon Seat Upholstery with one 5cm central band (Loop)	○
DSA2601		Standard Seat Upholstery with 1 central loop and 2 side hook	◆
DKA2610		Tension Adjustable Seat Upholstery with one 5cm central band (Loop)	€54
DSA2630		Light Seat Upholstery, mesh material to reduce weight	€80
SEAT CUSHION			
DSA2700		Standard Seat Cushion, 60mm high Double Foam with Visco Top Layer - Comfort/Low Risk (NC 2810, 2820, 2830, 2840, 2850)	€74
DSA2710		Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk (NC 2810, 2820, 2830, 2840, 2850)	€168
DSA2720		Matrx Libra Cushion, Black Smooth Cover 70mm High - Very High Risk	€645
DSA2748		Contour Visco NG Cushion 100mm High - High Risk (NC 2810, 2820, 2830, 2840, 2850)	€236
DSA2750		Vicair Active O2 lightweight cushion 60mm High - High risk (NC 2810, 2820, 2830, 2840, 2850, 2860)	€731
MATRX LIBRA CUSHION ACCESSORIES			
DSA2810		Fluid supplement sac (Pair)	€123
DSA2820		Tapered Wedge Pair	€28
DSA2830		Abductor	€58
DSA2840		Lateral Side Wedge	€46
DSA2850		Anterior Front Wedge	€43

SEAT RIGIDISER			
DSA2860		Seat Rigidiser	€79
FOOTREST*			
DSA3200		Footrest Tube, Rigid - Aluminium Hydroformed profile. Weight approx. 260g	◆
DSA3210		Footrest Tube, Rigid - Titanium Rounded profile tube. (NC 6962) Weight approx. 210g	€301
DSA3220		Carbon Footplate, Angle Adjustable (NC 3290, 6962)	€225
FOOTREST COVER			
DSA3290		Footrest Cover	€41
CALF STRAP*			
DSA3300		Calf strap Mounted on the frame with hook and loop	◆
DSA3305		Calf strap Mounted on the supporter with bracket	○
CLOTHESGUARDS/MUDGUARDS*			
DSA3410		Carbon Clothesguards, fixed between the back rest and seat rails	◆
DSA3430		Carbon Mudguards, fixed between the backrest and seat rails (NC 7230)	€327
SWO6666		Carbon fibre flip up clothesguards or removable carbon mudguards, for use with option SWO4444 or SWO9999 K-Series style folding backrest	○
FRAME COLOUR*			
DSA4000		Space Black (RAL 9005) (NC 4066)	○
DSA4010		Anthracite Grey (RAL 7016) (NC 4066)	○
DSA4020		Pure White (RAL 9016) (NC 4066, 4220)	○
DSA4030		Carmine Red (RAL 3002) (NC 4066)	○
DSA4040		Azure Blue (RAL 5017) (NC 4066)	○
DSA4060		Flashy Yellow (RAL 1516) (Additional 15 day lead time) (NC 4066)	€39
DSA4100		Black Glossy (RAL 9005) (Additional 15 day lead time) (NC 4066)	€39
DSA4130		Silver Grey (RAL 9006) (Additional 15 day lead time) (NC 4066)	€39
DSA4145		Pearl Night Blue (RAL 5026) (Additional 15 day lead time) (NC 4066)	€39
DSA4150		Ruby Red (RAL 1503) (Additional 15 day lead time) (NC 4066)	€39
DSA4170		Dynamite Orange (Additional 15 day lead time) (NC 4066)	€39
DSA4180		Forest Green (RAL 6029) (Additional 15 day lead time) (NC 4066)	€39
USA4101		Diamond Silver/Mirror Polish with standard varnish finish (Premium ceramic glaze finish available, for this option select USA4066)	€513
USA4102		Pastel Green with standard varnish finish (Premium ceramic glaze finish available, for this option select USA4066)	€413
USA4103		Sky Blue with standard varnish finish (Premium ceramic glaze finish available, for this option select USA4066)	€413
USA4104		Plum Purple with standard varnish finish (Premium ceramic glaze finish available, for this option select USA4066)	€413
USA4105		Salmon with standard varnish finish (Premium ceramic glaze finish available, for this option select USA4066)	€413
USA0008		Your personal frame colour. Add RAL colour of your choice in special comments box, for a personalised colour frame	€384
CERAMIC GLAZE FRAME FINISH			
USA4066		Ceramic glaze finish on frame, for added durability	€432
FRAME DECAL TYPE*			
DSA4190		Bold Decal (Full coloured transfer on frame) (NC 4300, 4310, 1111)	€39
DSA4195		Tag (Small Kuschall brand transfer on front of frame) (NC 4200, 4210, 4220, 4230, 4240, 4250, 2222)	◆
DSA4197		Without Tag (Kuschall branding only) (NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310)	○

BOLD DECAL COLOUR*			
DSA4200	<input type="checkbox"/>	Black	○
DSA4210	<input type="checkbox"/>	Dark Grey	○
DSA4220	<input type="checkbox"/>	Light Grey	○
DSA4230	<input type="checkbox"/>	Red	○
DSA4240	<input type="checkbox"/>	Blue	○
DSA4250	<input type="checkbox"/>	Gold	○
DSA1111	<input type="checkbox"/>	No bold colour, select only with Option DSA4195 Tag or DSA4197 Without Tag (NC 4190)	○
TAG COLOUR*			
DSA4300	<input type="checkbox"/>	Black	○
DSA4310	<input type="checkbox"/>	Red	○
DSA2222	<input type="checkbox"/>	No tag colour, select only with Option DSA4190 Bold Decal or DSA4197 Without Tag (NC 4195)	○
FRONT CASTORS*			
DSA5000	<input type="checkbox"/>	75mm x 25mm (3" x 1") Supersport (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DSA5050	<input type="checkbox"/>	100mm x 30mm (4" x 1.2") Everyday Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DSA5100	<input type="checkbox"/>	100mm x 30mm (4" x 1.2") Comfort Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DSA5150	<input type="checkbox"/>	100mm x 25mm (4" x 1") Starec Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	€102
DSA5200	<input type="checkbox"/>	100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options	€114
DSA5250	<input type="checkbox"/>	100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) (NC 5502, 5504, 5506, 5507, 5508, 5509)	€114
DSA5300	<input type="checkbox"/>	125mm x 30mm (5" x 1.2") Everyday Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DSA5350	<input type="checkbox"/>	125mm x 30mm (5" x 1.2") Comfort Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	◆
DSA5400	<input type="checkbox"/>	125mm x 25mm (5" x 1") Starec Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	€97
DSA5450	<input type="checkbox"/>	125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options	€114
DSA5500	<input type="checkbox"/>	125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs) (NC 5502, 5504, 5506, 5507, 5508, 5509)	€114
ANODISED COLOUR SOFTROLL CASTORS			
USA5502	<input type="checkbox"/>	Pink	€64
USA5504	<input type="checkbox"/>	Orange	€64
USA5506	<input type="checkbox"/>	Blue	€64
USA5507	<input type="checkbox"/>	Purple	€64
USA5508	<input type="checkbox"/>	Red	€64
USA5509	<input type="checkbox"/>	Gold	€64
CERAMIC CASTOR BEARINGS			
USA5510	<input type="checkbox"/>	Ceramic castor bearings	€64
ADJUSTABLE REAR WHEEL POSITION (SELECT POSITION FROM THE 5 OPTIONS BELOW)			
USA6050	<input type="checkbox"/>	Adjustable centre of gravity (adjustable +/- 1 position) Select from positions 1-5 below (15mm increments adjustment)	€572
SEAT POSITION (Measured from centre of backpost to centre of axle)*			
DSA6100	<input type="checkbox"/>	Position 1 = 140 mm MOST ACTIVE POSITION	○
DSA6110	<input type="checkbox"/>	Position 2 = 125 mm	○
DSA6120	<input type="checkbox"/>	Position 3 = 110 mm	○
DSA6130	<input type="checkbox"/>	Position 4 = 95 mm	◆
DSA6140	<input type="checkbox"/>	Position 5 = 80 mm LEAST ACTIVE POSITION	○

CAMBER OF CARBON AXLE*			
DSA6200		0° Camber	○
DSA6210		3° Camber	◆
DSA6220		6° Camber	○
QUAD QUICK RELEASE AXLE			
DSA6300		Quad Quick Release/Tetra Clip	(NC 6550, 6551, 6552, 6553, 6554, 0907) €96
REAR WHEEL SIZE*			
DSA6410		24"	◆
DSA6420		25"	(NC 6510, 6600, 6630, 6700, 6710, 6720, 6760, 6770, 6775, 6900, 6910) ○
REAR WHEEL TYPE*			
DSA6510		24" KÜSCHALL Starec Wheel, Black Cross Spokes, Rim & Hub.	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579) ◆
DSA6520		24", 25" High Performance Wheel Black hub and rim	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579) €311
DSA6550		24", 25" Spinergy L XK, Black hub, rim and 18 spokes in Black	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6910) €746
USA6551		24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6910) €1,009
USA6552		24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6910) €1,084
USA6553		24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6910) €1,147
USA6554		24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard)	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6910) €1,393
SWO0907		24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKÜSCHALL	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6571, 6572, 6573, 6575) €543
ANODISED COLOUR SPINERGY RIM (DESIGN LAB)			
USA6557		Blue	€230
USA6558		Gold	€230
USA6559		Green	€230
USA6560		Orange	€230
USA6561		Red	€230
USA6562		Purple	€230
ANODISED COLOUR SPINERGY HUB (DESIGN LAB)			
USA6563		Blue	€230
USA6564		Gold	€230
USA6565		Silver	€230
USA6566		Green	€230
USA6567		Orange	€230
USA6568		Red	€230
USA6569		Purple	€230

SPINERGY COLOURED SPOKES			
USA6570		Blue	€108
USA6571		Green	€108
USA6572		Orange	€108
USA6573		Pink	€108
USA6574		Red	€108
USA6575		Purple	€108
USA6576		White	€108
USA6577		Yellow	€108
CERAMIC BEARING REAR WHEEL			
USA6579		Ceramic Bearing for Spinergy rear wheels	€218
REAR WHEEL TYRE*			
DSA6600		24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection	◆
DSA6625		24", 25" x 1" Maxxis Alpha Protect. Pneumatic and puncture resistant.	€122
DSA6630		24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)	€107
UNIQUE CUSTOMISATION COLOUR PACKAGE*			
USA6960		UNIQUE KIT 1 - Includes custom colour forks and bearing cover. Select colour from 10 choices below	€167
USA6961		UNIQUE KIT 2 - Includes custom colour forks, bearing cover, backrest bracket with K-Series backrest, brake clamps. Select colour from the 10 choices below	€366
USA6962		UNIQUE KIT 3 - Includes custom colour forks, bearing cover, backrest bracket with K-Series backrest, brake clamps, rear wheel axle clamps, calf strap bracket. Select colour from the 10 choices below	€439
USA1234		Without colour package	○
UNIQUE CUSTOMISATION COLOURS			
USA6963		Red	○
USA6964		Purple	○
USA6965		Copper	○
USA6966		Green	○
USA6967		Orange	○
USA6968		Gold	○
USA6969		Champagne	○
USA6970		Pink	○
USA6971		Blue	○
USA6972		Mirror Finish	○

HANDRIM*			
DSA6700		24" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance	€334
DSA6710		24" Steel, mounted at 3 cm distance	€68
DSA6720		24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance	€200
DSA6740		24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, profile of 32mm diameter.	€280
DSA6750		24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, smaller profile of 24mm diameter for smaller hands.	€280
DSA6760		24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape, all metal black coated, low-friction surface for complete hand function for low-effort propulsion	€250
DSA6770		24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics without tricep function. Features a unique angle shaped profile and anti-slip coating providing a larger contact surface designed to be propelled using the heel of the hand or thumb pad.	€290
DSA6775		24" Carbolife Curve L Tetra handrim, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics with tricep function, but limited hand strength. Ergonomic shape with a flat, non-slip coated pushing surface for the thumb and a curved side for the fingers to wrap around.	€200
DSA6780		24", 25" Light Aluminium Grey, long tab with 2 holes for 3 cm or 4cm mountings, mounted at 3cm distance	◆
DSA6790		25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting	€51
SPOKEGUARD			
DSA6900		24" Transparent Spokeguards (When ordered with Carbolife handrims, handrim distance is set at 4cm)	€155
DSA6910		24" Black Spokeguard (When ordered with Carbolife handrims, handrim distance is set at 4cm)	€155
BRAKES*			
DSA7210		Standard Push to Lock Brakes, Lever Bent	◆
DSA7230		Performance Brakes, Push to Lock. Brake barrel flips away from tyre	€94
DSA7240		Active/Scissor Brakes	€94
DSA7250		Light/Compact Active, Aluminium Brakes	€248
SAFETY OPTIONS			
DSA7500		Anti-Tipper Right	€149
DSA7510		Anti-Tipper Left	€149
DSA7520		Anti-Tipper Both sides	€299
DSA7810		Bodypoint Posture Belt Evoflex 2 points, Small, Flat-Mount with Push Button. User Weight: 25-50KG	€171
DSA7820		Bodypoint Posture Belt Evoflex 2 points, Medium, Flat-Mount with Rehab Latch. User Weight: 50-75KG	€171
DSA7830		Bodypoint Posture Belt Evoflex, 2 points, Large, Flat-Mount with Rehab Latch. User Weight: 75-100KG	€171

OTHER OPTIONS			
DSA8010		Front Frame Protectors; 2 Piece. Removable with hook and loop	€43
DSA8500		Air Pump	€17
DSA8550		Tool Kit	◆
DSA8560		Kuschall branded Wheel bag (Fits all wheel sizes)	€58
DSA8680		24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included	€335
BODYPOINT ACCESSORIES			
DSA8800		Bodypoint Calf Strap	€80
DSA8850		Bodypoint Mobility Bag	€80

TOTAL PRICE	
--------------------	--

Non Conforming Options:

DSA1820	400 mm
DSA1825	410 mm
DSA1830	420 mm
DSA1840	440 mm
DSA1845	450 mm
DSA1850	460 mm
DSA1855	470 mm
DSA1860	480 mm
SWO4444	K-Series style lockdown folding backrest, supplied with carbon fibre flip up, fold in clothesguards and tension adjustable upholstery
SWO9999	K-Series style lockdown folding backrest, supplied with carbon fibre removable mudguards and tension adjustable upholstery
DSA2200	Standard Light Black Nylon upholstery (No tension adjustment)
DSA2220	Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener
DSA2240	Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required
SWO7777	Nylon tension adjustable upholstery for options 4444 & 9999
USA2241	Neon Yellow
USA2242	Neon Pink
USA2243	Blue
USA2244	White
USA2245	Red
USA0002	Your personal colour stitching. Add RAL colour of your choice in special comments box, for a personalised colour stitching on back upholstery.
USA0003	Your personal design embroidery. Supply your design with your order, subject to copyright. Design size up to 250mm x 250mm
USA0004	Your personal flock design. Supply your design with your order, subject to copyright. Design size up to 280mm x 400mm
DSA2710	Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk
DSA2720	Matrix Libra Cushion, Black Smooth Cover 70mm High - Very High Risk
DSA2748	Contour Visco NG Cushion 100mm High - High Risk
DSA2750	Vicair Active O2 lightweight cushion 60mm High - High risk
DSA2810	Fluid supplement sac (Pair)
DSA2820	Tapered Wedge Pair
DSA2830	Abductor
DSA2840	Lateral Side Wedge
DSA2850	Anterior Front Wedge
DSA2860	Seat Rigidiser
DSA3290	Footrest Cover
DSA3410	Carbon Clothesguards, fixed between the back rest and seat rails
DSA3430	Carbon Mudguards, fixed between the backrest and seat rails
SWO6666	Carbon fibre flip up clothesguards or removable carbon mudguards, for use with option SWO4444 or SWO9999 K-Series style folding backrest
USA4066	Ceramic glaze finish on frame, for added durability
DSA4190	Bold Decal (Full coloured transfer on frame)
DSA4195	Tag (Small Kuschall brand transfer on front of frame)
DSA4200	Black
DSA4210	Dark Grey
DSA4220	Light Grey

DSA4230	Red
DSA4240	Blue
DSA4250	Gold
DSA1111	No bold colour, select only with Option DSA4195 Tag or DSA4197 Without Tag
DSA4300	Black
DSA4310	Red
DSA2222	No tag colour, select only with Option DSA4190 Bold Decal or DSA4197 Without Tag
DSA5300	125mm x 30mm (5" x 1.2") Everyday Black hub and tyre
DSA5350	125mm x 30mm (5" x 1.2") Comfort Black hub and tyre
DSA5400	125mm x 25mm (5" x 1") Starec Black hub and tyre
DSA5450	125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options
DSA5500	125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs)
USA5502	Pink
USA5504	Orange
USA5506	Blue
USA5507	Purple
USA5508	Red
USA5509	Gold
USA5510	Ceramic castor bearings
DSA6410	24"
DSA6420	25"
DSA6510	24" KÜSCHALL Starec Wheel, Black Cross Spokes, Rim & Hub.
DSA6550	24", 25" Spinergy L XK, Black hub, rim and 18 spokes in Black
USA6551	24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
USA6552	24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
USA6553	24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
USA6554	24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard)
SWO0907	24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKÜSCHALL
USA6557	Blue
USA6558	Gold
USA6559	Green
USA6560	Orange
USA6561	Red
USA6562	Purple
USA6563	Blue
USA6564	Gold
USA6565	Silver
USA6566	Green
USA6567	Orange
USA6568	Red
USA6569	Purple
USA6570	Blue

USA6571	Green
USA6572	Orange
USA6573	Pink
USA6574	Red
USA6575	Purple
USA6576	White
USA6577	Yellow
USA6579	Ceramic Bearing for Spinergy rear wheels
DSA6600	24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection
DSA6630	24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)
USA6962	UNIQUE KIT 3 - Includes custom colour forks, bearing cover, backrest bracket with K-Series backrest, brake clamps, rear wheel axle clamps, calf strap bracket. Select colour from the 10 choices below
DSA6700	24" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance
DSA6710	24" Steel, mounted at 3 cm distance
DSA6720	24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance
DSA6760	24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape, all metal black coated, low-friction surface for complete hand function for low-effort propulsion
DSA6770	24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics without tricep function. Features a unique angle shaped profile and anti-slip coating providing a larger contact surface designed to be propelled using the heel of the hand or thumb pad.
DSA6775	24" Carbolife Curve L Tetra handrim, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics with tricep function, but limited hand strength. Ergonomic shape with a flat, non-slip coated pushing surface for the thumb and a curved side for the fingers to wrap around.
DSA6900	24" Transparent Spokeguards (When ordered with Carbolife handrims, handrim distance is set at 4cm)
DSA6910	24" Black Spokeguard (When ordered with Carbolife handrims, handrim distance is set at 4cm)
DSA7230	Performance Brakes, Push to Lock. Brake barrel flips away from tyre
DSA7810	Bodypoint Posture Belt Evoflex 2 points, Small, Flat-Mount with Push Button. User Weight: 25-50KG
DSA7820	Bodypoint Posture Belt Evoflex 2 points, Medium, Flat-Mount with Rehab Latch. User Weight: 50-75KG
DSA7830	Bodypoint Posture Belt Evoflex, 2 points, Large, Flat-Mount with Rehab Latch. User Weight: 75-100KG



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):