



Yes, you can:

Please refer to www.invacare.ie for a local distributor

NOVEMBER 2023
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action Ampla

Standard Features

- Max user weight: 300kg (47 stone)
- Crash tested ISO 7176-Part 19
- 24" quick release rear wheel solid tyre
- 6" front castor solid tyre
- 6 available seat widths: 20" or 22" standard (24" - 30" seat widths chargeable)
- 2 Seat depths 18" - 20"
- Folding backrest with cushion
- Swing in/out legrest hanger
- Flip up, angle and depth adjustable footplates
- Transfer aid/Armrest
- Choice of 2 colours
- Push to lock brakes
- Transportation kit (successful crash test achieved without posture belt)
- Lap belt
- Anti-tippers



Standard Base price: €4,130

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
SMB0001	ACTION AMPLA	<input type="radio"/>
SMBSWO	SPECIAL	POA
PROPULSION*		
SMB0500	SELF-PROPEL	€4,130

SEAT WIDTH*			
SMB1000	<input type="checkbox"/>	SW 20" (505mm)	◆
SMB1010	<input type="checkbox"/>	SW 22" (555mm)	○
SMB1020	<input type="checkbox"/>	SW 24" (605mm) TRANSPORTATION BEARING SYSTEM RECOMMENDED	€56
SMB1030	<input type="checkbox"/>	SW 26" (655mm) TRANSPORTATION BEARING SYSTEM RECOMMENDED	€56
SMB1040	<input type="checkbox"/>	SW 28" (705mm) TRANSPORTATION BEARING SYSTEM RECOMMENDED	€62
SMB1050	<input type="checkbox"/>	SW 30" (755mm) TRANSPORTATION BEARING SYSTEM RECOMMENDED	€62
SEAT DEPTH*			
SMB1100	<input type="checkbox"/>	SD 18" FROM FRONT OF BACKREST CUSHION TO FRONT OF CANVAS / 19" FROM FRONT OF BACKPOST TO FRONT OF CANVAS	◆
SMB1110	<input type="checkbox"/>	SD 20" FROM FRONT OF BACKREST CUSHION TO FRONT OF CANVAS / 21" FROM FRONT OF BACKPOST TO FRONT OF CANVAS	○
SEAT HEIGHT*			
SMB1200	<input type="checkbox"/>	SH 15"	○
SMB1210	<input type="checkbox"/>	SH 16"	○
SMB1220	<input type="checkbox"/>	SH 17"	◆
SMB1230	<input type="checkbox"/>	SH 18"	○
SEAT UPHOLSTERY*			
SMB1300	<input type="checkbox"/>	SEAT UPHOLSTERY	◆
BACKREST*			
SMB2000	<input type="checkbox"/>	FOLDING 17"	◆
SMB2010	<input type="checkbox"/>	FOLDING 19"	○
BACKREST SUPPORT*			
SMB2100	<input type="checkbox"/>	BACKREST CUSHION	◆
LEGREST HANGER RH*			
SMB3110	<input type="checkbox"/>	LEGREST RH	◆
LEGREST HANGER LH*			
SMB3210	<input type="checkbox"/>	LEGREST LH	◆
FOOTREST RH*			
SMB3400	<input type="checkbox"/>	FLIP UP RH LL 365-445MM	(NC 3550) ○
SMB3450	<input type="checkbox"/>	FLIP UP RH LL 445-525MM	(NC 3500) ○
FOOTREST LH*			
SMB3500	<input type="checkbox"/>	FLIP UP LH LL 365-445MM	○
SMB3550	<input type="checkbox"/>	FLIP UP LH LL 445-525MM	○
TRANSFER AID RH*			
SMB4110	<input type="checkbox"/>	TRANSFER AID RH	◆
TRANSFER AID LH*			
SMB4210	<input type="checkbox"/>	TRANSFER AID LH	◆
SEAT CUSHION*			
SWO8888	<input type="checkbox"/>	WITHOUT SEAT CUSHION	○
SMB5000	<input type="checkbox"/>	STANDARD	€37
FRAME COLOUR*			
SMB6000	<input type="checkbox"/>	TECH WHITE	○
SMB6010	<input type="checkbox"/>	TECH GREY BLACK	○
SMB6099	<input type="checkbox"/>	SPECIAL	POA

DECAL*			
SMB6500	<input type="checkbox"/>	BLUE	<input type="radio"/>
SMB6510	<input type="checkbox"/>	PINK	<input type="radio"/>
SMB6520	<input type="checkbox"/>	GREEN	<input type="radio"/>
FRONT CASTORS*			
SMB7010	<input type="checkbox"/>	6" X 1.3/8"	<input type="checkbox"/>
REAR WHEELS*			
SMB7100	<input type="checkbox"/>	24" X 1"	<input type="checkbox"/>
REAR WHEEL AXLE*			
SMB7200	<input type="checkbox"/>	QUICK RELEASE	<input type="checkbox"/>
BRAKES*			
SMB8000	<input type="checkbox"/>	PUSH TO LOCK	<input type="checkbox"/>
SMB8010	<input type="checkbox"/>	ATTENDANT BRAKES	(NC 9050) €265
TRANSPORTATION KIT*			
SMB8500	<input type="checkbox"/>	TRANSPORTATION KIT	<input type="checkbox"/>
POSTURE BELT*			
SMB8510	<input type="checkbox"/>	LAP BELT	<input type="checkbox"/>
SAFETY OPTIONS			
SMB8520	<input type="checkbox"/>	ANTI-TIPPER PAIR - CANNOT BE USED IN ADDITION TO TRANSPORTATION BEARING SYSTEM	(NC 9060) <input type="checkbox"/>
OTHER OPTIONS			
SMB9050	<input type="checkbox"/>	HT ADJ PUSHBAR	<input type="radio"/>
SMB9060	<input type="checkbox"/>	TRANSPORTATION BEARING SYSTEM RECOMMENDED FOR SEAT WIDTH 24" - 30" - CANNOT BE USED WITH STANDARD ANTI-TIPPERS	€208
SWO9999	<input type="checkbox"/>	HEADREST - EXTENDED LEAD TIME - DE259646	€166

TOTAL PRICE	
--------------------	--

Non Conforming Options:

SMB3500 FLIP UP LH LL 365-445MM

SMB3550 FLIP UP LH LL 445-525MM

SMB9050 HT ADJ PUSHBAR

SMB9060 TRANSPORTATION BEARING SYSTEM RECOMMENDED FOR SEAT WIDTH 24" - 30" - CANNOT BE USED WITH STANDARD ANTI-TIPPERS



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):