

# Scorpius A

AUTO-FOLD | COMPACT | STYLISH



The *Invacare Scorpius A* has been specifically designed for individuals who wish to maintain an active, independent lifestyle without compromising on space. Ideal for daily use, and fully flight safe, these folding mobility scooters are also the perfect travel companion for those who wish to take their scooter abroad.

# Easily folded and transportable

Compact enough to fit into the rear of a normal car the versatile **Scorpius A** is perfect for those individuals who are constantly on the move.

**Scorpius A** can be compactly folded in under 10 seconds with the single push of a button, or by using the remote control provided.



Compact when folded, the scooter can be easily placed inside a boot of any car. Weighing less than 30kg with the batteries removed, it can be effortlessly lifted when necessary. The impressively small footprint when folded means that it can be conveniently stored in small spaces.

#### **Airline Compliant**

For independence abroad, the standard 11.6Ah lithiumion battery weighs just 2.8kgs and is IATA and MSDS compliant for overseas travel on aeroplanes\*.

\*It is recommended to check with your airline before travelling.











### **Powered Auto Folding**

- ▶ Effortlessly folded for travel and storage
- ► Scorpius A folded dimensions = 775 (I) x 512 (w) x 540 (h) mm



#### **Touch LED panel**

- Intuitive, touch screen technology
- ▶ Able to fold the **Scorpius A** in less than 8 seconds!
- Equipped with an anti-theft device

### **Airline compliant battery**

- ► Standard 11.6Ah lithium-ion battery can be taken on aeroplanes
- ▶ IATA and MSDS compliant

▶ With standard 8" PU wheels, the Scorpius is suitable for various types of terrain

#### **Additional features**

Storage bags

lights

Optional front, rear and under seat storage bags are available as well as a convenient storage cover

▶ 17.6Ah battery

A larger lithium-ion battery available for those individuals who wish to travel longer distances (14 Miles)

\*not airline compliant



**Durability** 

#### Technical data •

For more information about this product, including the product's user manual, please visit www.invacare.co.uk











width 🔻

Seat depth •

height 🔻

 $(L \times W \times H)$ 

width 🔻

Scorpius A

340mm

419 mm

560mm

775 x 512 x 540 mm

512mm

953 mm



Total height v



Total weight\* ▼



Max. user weight v



Battery capacity •



Turning diameter v



Max. climbable obstacle height\*\* ▼

Scorpius A

953 mm

28 kg

115kg / 18 st

1 x 24V/11.6Ah

1140mm

50mm / 2 inch\*\*



Ground clearance 🔻



Max. safe slope



Max. driving



Max. driving speed ▼



Motor capacity ▼

Scorpius A

43mm /1.7in

6°


distance ▼ 9 miles

6kmp/h

250W, 4600 r.p.m

\*Weight without batteries \*\* Measured with run-up

### Colour of shrouds •



Invacare Limited Pencoed Technology Park - Pencoed Bridgend CF35 5AQ - United Kingdom Tel: +44 1 656 776 222 Fax: +44 1656 776 220 E-mail: uk@invacare.com Sales Order E-mail: ordersuk@invacare.com

## Upholstery •



**Black mesh** fabric

Invacare Ireland Ltd Unit 5 - Seatown Business Campus Seatown Road - SWORDS County Dublin - Ireland Tel. +353 18107084 Fax +353 1 8107085 Email: Ireland@invacare.com







