

# INVACARE® MATRIX® E2™ DEEP BACK

## HARDWARE & ACCESSORY REFERENCE GUIDE

		WIDTH								HARDWARE OPTIONS			
		14"	15"	16"	17"	18"	19"	20"					
HEIGHT	6"	MiniSet Hardware*	●	●	●	●	●	●	●	  	MiniSet	Fixed MiniSet	Heavy-Duty MiniSet
		Fixed MiniSet Hardware**	●	●	●	●	●	●	●				
		Heavy Duty MiniSet Hardware*	●	●	●	●	●	●	●				
	8"	MiniSet Hardware*	●	●	●	●	●	●	●				
		Fixed MiniSet Hardware**	●	●	●	●	●	●	●				
		Heavy Duty MiniSet Hardware*	●	●	●	●	●	●	●				
	10"	MiniSet Hardware*	●	●	●	●	●	●	●				
		Fixed MiniSet Hardware**	●	●	●	●	●	●	●				
		Heavy Duty MiniSet Hardware*	●	●	●	●	●	●	●				
	12"	MiniSet Hardware*	●	●	●	●	●	●	●		EasySet	Fixed EasySet	
		Fixed MiniSet Hardware**	●	●	●	●	●	●	●				
		Heavy Duty MiniSet Hardware*	●	●	●	●	●	●	●				
		EasySet Hardware**	●	●	●	●	●	●	●				
	14"	Heav-Duty EasySet Hardware**	●	●	●	●	●	●	●				
		MiniSet Hardware*	●	●	●	●	●	●	●				
		Fixed MiniSet Hardware**	●	●	●	●	●	●	●				
		Heavy Duty MiniSet Hardware*	●	●	●	●	●	●	●				
	16"	EasySet Hardware**	●	●	●	●	●	●	●		Heavy-Duty EasySet	Split	
		Heav-Duty EasySet Hardware**	●	●	●	●	●	●	●				
		MiniSet Hardware*	●	●	●	●	●	●	●				
		Fixed MiniSet Hardware**	●	●	●	●	●	●	●				
	18"	Heavy Duty MiniSet Hardware*	●	●	●	●	●	●	●				
		EasySet Hardware**	●	●	●	●	●	●	●				
		Heav-Duty EasySet Hardware**	●	●	●	●	●	●	●				
MiniSet Hardware*		●	●	●	●	●	●	●					
20"	Fixed MiniSet Hardware**	●	●	●	●	●	●	●					
	Heavy Duty MiniSet Hardware*	●	●	●	●	●	●	●					
	EasySet Hardware**	●	●	●	●	●	●	●					
	Heav-Duty EasySet Hardware**	●	●	●	●	●	●	●					

● Standard Hardware ● Optional Hardware

NOTE: \* Hardware Fits Chairs 1" narrower or 1" wider than back model width  
 \*\* Hardware Fits Chairs 1" narrower or 2" wider than back model width

		WIDTH							
		14"	15"	16"	17"	18"	19"	20"	
6"	Back Bone PSIS/Sacrolumbar Support	●	●	●	●	●	●	●	
	Privacy Flap	●	●	●	●	●	●	●	
	Posture Flex Abdominal Support/Quad Release	●	●	●	●	●	●	●	
8"	Back Bone PSIS/Sacrolumbar Support	●	●	●	●	●	●	●	
	Privacy Flap	●	●	●	●	●	●	●	
	Posture Flex Abdominal Support/Quad Release	●	●	●	●	●	●	●	
10"	EasyFix Laterals*	●	●	●	●	●	●	●	
	Standard Swing-Away Laterals*	●	●	●	●	●	●	●	
	Back Bone PSIS/Sacrolumbar Support*	●	●	●	●	●	●	●	
	Privacy Flap	●	●	●	●	●	●	●	
12"	Posture Flex Abdominal Support/Quad Release	●	●	●	●	●	●	●	
	EasyFix Laterals*	●	●	●	●	●	●	●	
	Standard Swing-Away Laterals*	●	●	●	●	●	●	●	
	Back Bone PSIS/Sacrolumbar Support*	●	●	●	●	●	●	●	
14"	Back Bone PSIS/Sacrolumbar Support	●	●	●	●	●	●	●	
	Privacy Flap	●	●	●	●	●	●	●	
	Posture Flex Abdominal Support/Quad Release	●	●	●	●	●	●	●	
	Headrest Interface Plate	●	●	●	●	●	●	●	
	Heavy-Duty Headrest Interface Plate	●	●	●	●	●	●	●	
	Chest Harness Interface	●	●	●	●	●	●	●	
	EasyFix Laterals*	●	●	●	●	●	●	●	
	Standard Swing-Away Laterals*	●	●	●	●	●	●	●	
	Offset Swing-Away Laterals*	●	●	●	●	●	●	●	
16"	Back Bone PSIS/Sacrolumbar Support	●	●	●	●	●	●	●	
	Privacy Flap	●	●	●	●	●	●	●	
	Posture Flex Abdominal Support/Quad Release	●	●	●	●	●	●	●	
	Headrest Interface Plate	●	●	●	●	●	●	●	
	Heavy-Duty Headrest Interface Plate	●	●	●	●	●	●	●	
	Chest Harness Interface	●	●	●	●	●	●	●	
	EasyFix Laterals*	●	●	●	●	●	●	●	
	Standard Swing-Away Laterals*	●	●	●	●	●	●	●	
	Offset Swing-Away Laterals*	●	●	●	●	●	●	●	
18"	Back Bone PSIS/Sacrolumbar Support	●	●	●	●	●	●	●	
	Privacy Flap	●	●	●	●	●	●	●	
	Posture Flex Abdominal Support/Quad Release	●	●	●	●	●	●	●	
	Headrest Interface Plate	●	●	●	●	●	●	●	
	Heavy-Duty Headrest Interface Plate	●	●	●	●	●	●	●	
	Chest Harness Interface	●	●	●	●	●	●	●	
	EasyFix Laterals*	●	●	●	●	●	●	●	
	Standard Swing-Away Laterals*	●	●	●	●	●	●	●	
	Offset Swing-Away Laterals*	●	●	●	●	●	●	●	
20"	Back Bone PSIS/Sacrolumbar Support	●	●	●	●	●	●	●	
	Privacy Flap	●	●	●	●	●	●	●	
	Posture Flex Abdominal Support/Quad Release	●	●	●	●	●	●	●	
	Headrest Interface Plate	●	●	●	●	●	●	●	
	Heavy-Duty Headrest Interface Plate	●	●	●	●	●	●	●	
	Chest Harness Interface	●	●	●	●	●	●	●	
	EasyFix Laterals*	●	●	●	●	●	●	●	
	Standard Swing-Away Laterals*	●	●	●	●	●	●	●	
	Offset Swing-Away Laterals*	●	●	●	●	●	●	●	

## ACCESSORIES



Headrest Interface Plate  
Standard Heavy-Duty



Chest Harness Interface



EasyFix



Standard Swing-Away



Offset Swing-Away



Back Bone



Privacy Flap



Posture Flex Abdominal Support Quad Release Tab

● Standard ● Optional **NOTE:** \* Lateral options are NOT available on Matrix E2 Deep Backs

Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1656 776 222  
Fax: +44 1656 776 220  
E-mail: uk@invacare.com  
Sales Order E-mail: ordersuk@invacare.com  
www.invacare.co.uk

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: Ireland@invacare.com  
www.invacare.ie

**matrix**  
SEATING SERIES

**INVACARE**

Yes, you can.®