



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

07 JANUARY 2026  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

**Action 5 Folding**

**Standard Features**



- Max user weight 100kg - 130kg (Pending seat width)
- Crash tested ISO 7176 - Part 19
- Self Propulsion with 22" or 24" quick release rear wheels
- Choice of 7 seat widths 15" - 21" (380mm to 530mm)
- Choice of 7 seat depths 14" - 20" (350mm to 500mm)
- Choice of 8 front seat to floor heights 14" - 21" (360mm to 535mm)
- Choice of 7 rear seat to floor heights 14" - 20" (360mm to 510mm)
- Nylon seat upholstery
- Fixed backrest (90° to the seat rail)
- Fixed backrest, nylon and mesh padded cover
- Fixed front end or 80° swing in/swing out legrest hangers
- Flip up footplates
- Sideguards / skirtguards
- Choice of colours and 2 component accents
- 2° Camber
- Push to lock brakes
- Transportation kit with lap belt
- Adjustable heel straps
- User tool kit

**Standard Base price: €1,695**

Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
SJC0004	ACTION 5 FOLDING	<input type="radio"/>
SJCSPEC	SPECIAL	POA
<b>FRAME*</b>		
SJC0200	SWING AWAY FRAME <small>(NC 3000, 3100, 3230, 3440, 3450)</small>	<input type="radio"/>
SJC0210	FIXED FRONT END 80 DEGREE ANGLE FRAME <small>(NC 1111, 3010, 3020, 3030, 3040, 3050, 3110, 3120, 3130, 3140, 3150, 3333, 3400, 3410, 3420, 3430)</small>	<input type="radio"/>
<b>PROPULSION*</b>		
SJC0300	SELF PROPULSION	€1,695

<b>SEAT WIDTH*</b>			
SJC1030		15" SW	○
SJC1040		16" SW	○
SJC1050		17" SW	○
SJC1060		18" SW	○
SJC1070		19" SW	○
SJC1075		20" SW	(NC 1110, 1170) ○
SJC1078		21" SW	(NC 1110, 1170) ○
<b>SEAT DEPTH*</b>			
SJC1100		14" SD	(NC 1270, 3040, 3140, 7080, 7100, 7110, 7600, 7610) ○
SJC1110		15" SD	(NC 1270, 7600) ○
SJC1120		16" SD	(NC 1270) ○
SJC1130		17" SD	○
SJC1140		18" SD	○
SJC1150		19" SD	○
SJC1160		20" SD	○
SJC1170		14" - 16" ADJ SD	(NC 1510, 7600) ○
SJC1180		16" - 18" ADJ SD	(NC 1510) ○
SJC1190		18" - 20" ADJ SD	(NC 1510) ○
<b>FRONT SEAT HEIGHT*</b>			
SJC1200		14" FSH	(NC 1310, 1320, 1330, 1340, 1350, 1360, 3420, 3430, 3450, 7040, 7055, 7070, 7080, 7100, 7110) ○
SJC1210		15" FSH	(NC 1320, 1330, 1340, 1350, 1360, 3430, 7000, 7040, 7070, 7080, 7100, 7110) ○
SJC1220		16" FSH	(NC 1330, 1340, 1350, 1360, 7080, 7100, 7110) ○
SJC1230		17" FSH	(NC 1300, 1340, 1350, 1360, 7110) ○
SJC1240		18" FSH	(NC 1300, 1310, 1350, 1360) ○
SJC1250		19" FSH	(NC 1300, 1310, 1320, 1360) ○
SJC1260		20" FSH	(NC 1300, 1310, 1320, 1330, 7000, 7040) ○
SJC1270		21" FSH	(NC 1300, 1310, 1320, 1330, 1340, 3030, 3130, 7000, 7040, 7055, 7070, 7080, 7100) ○
<b>REAR SEAT HEIGHT - MUST BE EQUAL TO OR LOWER THAN FRONT SEAT HEIGHT*</b>			
SJC1300		14" RSH	(NC 7250, 7260) ○
SJC1310		15" RSH	○
SJC1320		16" RSH	○
SJC1330		17" RSH	○
SJC1340		18" RSH	○
SJC1350		19" RSH	(NC 7200) ○
SJC1360		20" RSH	(NC 7200, 7230, 7240) ○
<b>SEAT UPHOLSTERY*</b>			
SJC1500		DEPTH ADJUSTABLE	€26
SJC1510		NYLON FIXED	◆
SJC1520		TENSION AND DEPTH ADJUSTABLE	€70
<b>BACKREST ATTACHMENT*</b>			
SJC2000		FIXED 90°	◆
SJC2020		ANGLE ADJUSTABLE (-15° TO +15°) AND FOLDING	€83

<b>BACKREST FIXED HEIGHT - INCLUDING FIXED PUSH HANDLES*</b>				
SJC2143		14" BACK HT	(NC 2300, 2310, 2400, 2410, 2420, 2505, 8160, 8200, 8888)	○
SJC2150		15" BACK HT	(NC 2300, 2310, 2400, 2410, 2420, 8200, 8888)	○
SJC2153		16" BACK HT	(NC 2300, 2310, 2400, 2410, 2420, 8888)	○
SJC2160		17" BACK HT	(NC 2300, 2310, 2400, 2410, 2420, 8888)	○
SJC2163		18" BACK HT	(NC 2300, 2310, 2400, 2410, 2420, 8888)	○
SJC2170		19" BACK HT	(NC 2300, 2310, 2400, 2410, 2420, 8888)	○
SJC2173		20" BACK HT	(NC 2300, 2310, 2400, 2410, 2420, 8888)	○
SJC7777		No Fixed backrest, use only when selecting height adjustable backrest, option 2300 oe 2310	(NC 8888)	○
<b>BACKREST HEIGHT ADJUSTABLE - PUSH HANDLES MUST ALSO BE SELECTED IN THE NEXT SECTION BELOW*</b>				
SJC2300		HT ADJ 14.5" - 16"	(NC 2505)	€44
SJC2310		HT ADJ 16" - 19.5"		€44
SJC8888		No height adjustable backrest, use only when selecting Fixed backrest, options 2143-2173		○
<b>PUSH HANDLES - WHEN SELECTED WITH HEIGHT ADJUSTALE BACKREST</b>				
SJC2400		W/O PUSH HANDLES	(NC 8160, 9100, 9200)	○
SJC2410		REAR SET	(NC 2505)	€93
SJC2420		INTEGRATED		€148
<b>BACKREST UPHOLSTERY*</b>				
SJC2500		W/O BACK UPH - DELIVERED WITH SIMPLE SLEEVE		○
SJC2505		FIXED NYLON AND MESH	(NC 5500)	◆
SJC2510		CONTOUR ADJUSTABLE		€89
<b>LEGREST HANGER RH*</b>				
SJC3000		W/O HANGER RH	(NC 3230, 5400)	○
SJC3010		70° RH	(NC 1111, 3440, 3450, 22222)	○
SJC3020		80° RH	(NC 1111, 3440, 3450, 22222)	○
SJC3030		90° RH	(NC 1111, 3200, 3440, 3450, 7080, 7100, 7110, 22222)	○
SJC3040		ELEVATING LEGREST RH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GOUND IF UL IS MORE THAN 400MM	(NC 3230, 3400, 3440, 3450, 5400, 22222, 8120)	€135
SJC3050		AMPUTEE LEGREST RH	(NC 1111, 3200, 3210, 3230, 5400, 8100, 8120)	€156
<b>LEGREST HANGER LH*</b>				
SJC3100		W/O HANGER LH	(NC 5400)	○
SJC3110		70° LH	(NC 3333, 3440, 3450)	○
SJC3120		80° LH	(NC 3333, 3440, 3450)	○
SJC3130		90° LH	(NC 3300, 3333, 3440, 3450, 7080, 7100, 7110)	○
SJC3140		ELEVATING LEGREST LH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GOUND IF UL IS MORE THAN 400MM	(NC 3400, 3440, 3450, 44444, 5400, 8120)	€135
SJC3150		AMPUTEE LEGREST LH	(NC 3333, 3300, 3310, 5400, 8100, 8120)	€156

FOOTREST RH*			
SJC3200		FLIP UP RH (NC 3410)	◆
SJC3210		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS.	€51
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)	€51
SJC3230		ONE PIECE (NC 3300, 3310, 3333, 3400, 3410, 3420, 3430)	€216
SJC22222		NO FOOTPLATE REQUIRED	○
FOOTREST LH*			
SJC3300		FLIP UP LH (NC 3410)	◆
SJC3310		FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS.	€51
SWO3333		FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)	€51
SJC44444		NO FOOTPLATE REQUIRED	○
FOOTREST TUBE*			
SJC3400		SWING AWAY FRAME HIGH MOUNTED TUBE	○
SJC3410		SWING AWAY FRAME SHORT TUBE	○
SJC3420		SWING AWAY FRAME STANDARD TUBE	○
SJC3430		SWING AWAY FRAME LONG TUBE	○
SJC3440		FIXED FRAME STANDARD TUBE	○
SJC3450		FIXED FRAME LONG TUBE	○
MUDGUARDS/SIDEGUARDS*			
SJC3700		W/O SIDEGUARD/MUDGUARD	○
SJC3750		SIDEGUARDS (NC 4010, 4020, 4030, 4040, 4110, 4120, 4130, 4140)	◆
SJC3800		MUDGUARDS (NC 4010, 4020, 4030, 4040, 4110, 4120, 4130, 4140, 8160)	€88
ARMREST RH*			
SJC4000		W/O ARMREST	○
SJC4010		T-ARM DESK ARMPAD RH (NC 7600, 7610)	€88
SJC4020		T-ARM FULL ARMPAD RH	€88
SJC4030		SWING AWAY DESK ARMPAD RH (NC 5700)	€120
SJC4040		SWING AWAY FULL ARMPAD RH (NC 5700)	€120
ARMREST LH*			
SJC4100		W/O ARMREST	○
SJC4110		T-ARM DESK ARMPAD LH (NC 7600, 7610)	€88
SJC4120		T-ARM FULL ARMPAD LH	€88
SJC4130		SWING AWAY DESK ARMPAD LH (NC 5700)	€120
SJC4140		SWING AWAY FULL ARMPAD LH (NC 5700)	€120

<b>SEAT CUSHION*</b>			
SJC55555	<input type="checkbox"/>	WITHOUT SEAT CUSHION	○
SJC5000	<input type="checkbox"/>	BASIC FOAM CUSHION	€24
<b>SEATING OPTIONS</b>			
SJC5100	<input type="checkbox"/>	CLIP TO KEEP THE CHAIR FOLDED	€19
SJC5400	<input type="checkbox"/>	CALF STRAP	○
SJC5500	<input type="checkbox"/>	FOLDING BACKBRACE	€24
SJC5700	<input type="checkbox"/>	W/O AUTOLOCK SAFETY DEVICE	○
<b>FRAME COLOUR*</b>			
SJC6000	<input type="checkbox"/>	PEARL WHITE	○
SJC6010	<input type="checkbox"/>	ANTHRACITE BLACK	○
SJC6015	<input type="checkbox"/>	BLUE	○
SJC6040	<input type="checkbox"/>	ELECTRIC GREEN	○
SJC6045	<input type="checkbox"/>	DARK PURPLE	○
SJC6060	<input type="checkbox"/>	PEARL GREY	○
SJC6070	<input type="checkbox"/>	MATT BLACK	○
SJC6099	<input type="checkbox"/>	SPECIAL COLOUR	€164
SWO10002	<input type="checkbox"/>	SLATE GREY - SPECIAL COLOUR, 3 WEEK LEAD TIME	€164
SWO10003	<input type="checkbox"/>	RED MATT - SPECIAL COLOUR, 3 WEEK LEAD TIME	€164
SWO10004	<input type="checkbox"/>	POLAR BLUE - SPECIAL COLOUR, 3 WEEK LEAD TIME	€164
SWO10005	<input type="checkbox"/>	CANDY PINK - SPECIAL COLOUR, 3 WEEK LEAD TIME	€164
<b>ACCENT COLOUR*</b>			
SJC6100	<input type="checkbox"/>	WITHOUT ACCENT COLOUR SILVER FORK AND REARWHEEL SUPPORT	○
SJC6110	<input type="checkbox"/>	BLUE ACCENT	○
SJC6200	<input type="checkbox"/>	PURPLE ACCENT	○
<b>FRONT CASTORS*</b>			
SJC7000	<input type="checkbox"/>	4" X 1.1/4" SOFTROLL	€46
SJC7040	<input type="checkbox"/>	4" X 1.3/8" ALU SOFTROLL	€46
SJC7055	<input type="checkbox"/>	5" X 1" SOLID	◆
SJC7070	<input type="checkbox"/>	5" X 1.3/8" ALU SOFTROLL	€46
SJC7080	<input type="checkbox"/>	6" X 1" SOLID	○
SJC7100	<input type="checkbox"/>	6" X 1.3/8" ALU SOFTROLL	€46
SJC7110	<input type="checkbox"/>	7" X 1.1/4" SOLID	○

<b>REAR WHEELS*</b>			
SJC7200		22" (560 MM) STANDARD REINFORCED, 6 PTS HANDRIM DISTANCE WHEEL HANDRIM 11MM OR 16MM (SET AT 11MM) CROSS SPOKED WHEEL	(NC 7320, 7410) ○
SJC7230		24" (610 MM) STANDARD REINFORCED, 6 PTS HANDRIM DISTANCE WHEEL HANDRIM 11MM OR 16MM (SET AT 11MM) CROSS SPOKED WHEEL	(NC 7320) ◆
SJC7240		24" (610 MM) LIGHTWEIGHT, RADIAL SPOKED, BLACK RIM AND HUB - 6 POINTS HANDRIM DISTANCE WHEEL HANDRIM 11MM OR 16MM (SET AT 11MM) CROSS SPOKED WHEEL	(NC 7410, 8000, 8160) €70
SJC7250		25" (635 MM) LIGHTWEIGHT, RADIAL SPOKED, BLACK RIM AND HUB - 6 POINTS HANDRIM DISTANCE WHEEL HANDRIM 11MM OR 16MM (SET AT 11MM) CROSS SPOKED WHEEL	(NC 7302, 7304, 7315, 7410, 8000, 8160) €70
SJC7260		25" (635 MM) SPINERGY WHEELS 6 PTS - BLACK DISTANCE WHEEL HANDRIM 16MM OR 22MM (SET AT 16MM) CROSS SPOKED WHEEL	(NC 7302, 7304, 7315, 7410, 7420, 8000, 8160) €688
<b>REAR WHEEL TYRE*</b>			
SJC7302		SOLID BLACK 1.3/8"	◆
SJC7304		PNEUMATIC TYRE BLACK 1.3/8"	○
SJC7315		PNEUMATIC LOW PROFILE BLACK 1"	○
SJC7320		PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM	€48
<b>HANDRIM*</b>			
SJC7400		ALUMINIUM HANDRIM	◆
SJC7410		STAINLESS STEEL	○
SJC7420		HANDRIM COATED IN BLACK ANTI-SLIP	€62
<b>REAR WHEEL AXLE*</b>			
SJC7500		Q/R AXLE	◆
<b>REAR WHEEL POSITION*</b>			
SJC7600		VERY ACTIVE	○
SJC7610		ACTIVE	○
SJC7620		MEDIUM ACTIVE	○
SJC7630		STANDARD	◆
<b>CAMBER ANGLE*</b>			
SJC7700		0° CAMBER	○
SJC7710		1° CAMBER	○
SJC7720		2° CAMBER	◆
SJC7730		3° CAMBER	○
SJC7740		4° CAMBER	○
<b>CHASSIS OPTION</b>			
SJC8000		SPOKE PROTECTORS	€69
<b>PUMP</b>			
SJC8030		PUMP - HIGH PRESSURE	€48
<b>BRAKES*</b>			
SJC8100		PUSH TO LOCK BRAKES	◆
SJC8110		PULL TO LOCK BRAKES	○
SJC8120		PUSH TO LOCK BRAKES WITH EXTENTION	€16
SJC8130		PULL TO LOCK BRAKES WITH EXTENTION	€16
SJC8160		ATTENDANT HUB BRAKES INCLUDES STANDARD USER PUSH TO LOCK BRAKE (OR PULL VERSION IF IT APPEARS AS THE ONLY CHOICE ABOVE)	€186

TRANSPORTATION KIT			
SJC8200	<input type="checkbox"/>	TRANSPORTATION BRACKETS	◆
SAFETY OPTIONS*			
SJC8300	<input type="checkbox"/>	LAP BELT	◆
SJC8500	<input type="checkbox"/>	REAR ANTI-TIPPER ACTIVE RH	(NC 9200) €76
SJC8550	<input type="checkbox"/>	REAR ANTI-TIPPER ACTIVE LH	(NC 9100) €76
SJC8600	<input type="checkbox"/>	REFLECTORS KIT (WHEELS)	€17
ATTENDANT OPTION			
SJC9000	<input type="checkbox"/>	CANE HOLDER	€20
SJC9100	<input type="checkbox"/>	TIPPING AID LH	€62
SJC9200	<input type="checkbox"/>	TIPPING AID RH	€62
OTHER OPTIONS			
SWO10001	<input type="checkbox"/>	ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	€105

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

SJC1110	15" SD
SJC1170	14" - 16" ADJ SD
SJC1270	21" FSH
SJC1300	14" RSH
SJC1310	15" RSH
SJC1320	16" RSH
SJC1330	17" RSH
SJC1340	18" RSH
SJC1350	19" RSH
SJC1360	20" RSH
SJC1510	NYLON FIXED
SJC2300	HT ADJ 14.5" - 16"
SJC2310	HT ADJ 16" - 19.5"
SJC8888	No height adjustable backrest, use only when selecting Fixed backrest, options 2143-2173
SJC2400	W/O PUSH HANDLES
SJC2410	REAR SET
SJC2420	INTEGRATED
SJC2505	FIXED NYLON AND MESH
SJC3000	W/O HANGER RH
SJC3010	70° RH
SJC3020	80° RH
SJC3030	90° RH
SJC3040	ELEVATING LEGREST RH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GOUND IF UL IS MORE THAN 400MM
SJC3050	AMPUTEE LEGREST RH
SJC3100	W/O HANGER LH
SJC3110	70° LH
SJC3120	80° LH
SJC3130	90° LH
SJC3140	ELEVATING LEGREST LH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GOUND IF UL IS MORE THAN 400MM
SJC3150	AMPUTEE LEGREST LH
SJC3200	FLIP UP RH
SJC3210	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
SJC3230	ONE PIECE
SJC22222	NO FOOTPLATE REQUIRED
SJC3300	FLIP UP LH
SJC3310	FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO3333	FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
SJC44444	NO FOOTPLATE REQUIRED

SJC3400	SWING AWAY FRAME HIGH MOUNTED TUBE
SJC3410	SWING AWAY FRAME SHORT TUBE
SJC3420	SWING AWAY FRAME STANDARD TUBE
SJC3430	SWING AWAY FRAME LONG TUBE
SJC3440	FIXED FRAME STANDARD TUBE
SJC3450	FIXED FRAME LONG TUBE
SJC4010	T-ARM DESK ARMPAD RH
SJC4020	T-ARM FULL ARMPAD RH
SJC4030	SWING AWAY DESK ARMPAD RH
SJC4040	SWING AWAY FULL ARMPAD RH
SJC4110	T-ARM DESK ARMPAD LH
SJC4120	T-ARM FULL ARMPAD LH
SJC4130	SWING AWAY DESK ARMPAD LH
SJC4140	SWING AWAY FULL ARMPAD LH
SJC5400	CALF STRAP
SJC5500	FOLDING BACKBRACE
SJC5700	W/O AUTOLOCK SAFETY DEVICE
SJC7000	4" X 1.1/4 SOFTROLL
SJC7040	4" X 1.3/8" ALU SOFTROLL
SJC7055	5" X 1" SOLID
SJC7070	5" X 1.3/8" ALU SOFTROLL
SJC7080	6" X 1" SOLID
SJC7100	6" X 1.3/8" ALU SOFTROLL
SJC7110	7" X 1.1/4" SOLID
SJC7200	22" (560 MM) STANDARD REINFORCED, 6 PTS HANDRIM DISTANCE WHEEL HANDRIM 11MM OR 16MM (SET AT 11MM) CROSS SPOKED WHEEL
SJC7230	24" (610 MM) STANDARD REINFORCED, 6 PTS HANDRIM DISTANCE WHEEL HANDRIM 11MM OR 16MM (SET AT 11MM) CROSS SPOKED WHEEL
SJC7240	24" (610 MM) LIGHTWEIGHT, RADIAL SPOKED, BLACK RIM AND HUB - 6 POINTS HANDRIM DISTANCE WHEEL HANDRIM 11MM OR 16MM (SET AT 11MM) CROSS SPOKED WHEEL
SJC7250	25" (635 MM) LIGHTWEIGHT, RADIAL SPOKED, BLACK RIM AND HUB - 6 POINTS HANDRIM DISTANCE WHEEL HANDRIM 11MM OR 16MM (SET AT 11MM) CROSS SPOKED WHEEL
SJC7260	25" (635 MM) SPINERGY WHEELS 6 PTS - BLACK DISTANCE WHEEL HANDRIM 16MM OR 22MM (SET AT 16MM) CROSS SPOKED WHEEL
SJC7302	SOLID BLACK 1.3/8"
SJC7304	PNEUMATIC TYRE BLACK 1.3/8"
SJC7315	PNEUMATIC LOW PROFILE BLACK 1"
SJC7320	PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM
SJC7410	STAINLESS STEEL
SJC7420	HANDRIM COATED IN BLACK ANTI-SLIP
SJC7600	VERY ACTIVE
SJC7610	ACTIVE
SJC8000	SPOKE PROTECTORS
SJC8100	PUSH TO LOCK BRAKES
SJC8120	PUSH TO LOCK BRAKES WITH EXTENTION

SJC8160	ATTENDANT HUB BRAKES INCLUDES STANDARD USER PUSH TO LOCK BRAKE (OR PULL VERSION IF IT APPEARS AS THE ONLY CHOICE ABOVE)
SJC8200	TRANSPORTATION BRACKETS
SJC9100	TIPPING AID LH
SJC9200	TIPPING AID RH



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):