



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JULY 2022
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Bora Modulite Stock

Standard Features

- Max User Weight 130Kg (20.5 stone)
- Crash tested ISO7176-Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable Seat set at width 430mm (17") depth 460mm (18") with an additional 50mm, (2") by arms
- Seat Depths 410mm - 510mm Stepless Adjustment (16" - 20")
- Tension Adjustable backrest
- Without seat cushion
- Standard T-Armrests, detachable with Short armpad
- Manual settable seat Tilt (0° - 20°)
- Manual settable backrest Recline
- 80° hangers with fixed Footplate
- LiNX Control REM211 with swing away bracket RH
- Black Body Shroud, Black Frame
- 200mm x 45mm (8" x 1.175") Puncture proof solid PU-foam, black castors
- 318mm x 57mm (12.5" x 2.25") Puncture proof solid inlay Black drive wheels
- Silver Rims
- Standard Rear suspension
- 6 km/h (4mph) 2 Pole Motors and 50Ah Batteries
- 8A (200 - 250V) Charger
- Bolt on or quick release anti tippers
- 4 Point tie down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock (adjustable both sides)
- Kerbclimber with park position - (NORMALLY A COST OPTION OF £74)

Standard Base price: €2,450



Q'STRAIT
QLK-150
Compatible

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code **Description**

Bora Modulite Stock	
1661626 <input type="checkbox"/>	Bora Modulite Stock €2,450

TOTAL PRICE

Non Conforming Options:



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):