

Spare Parts Catalogue

Version: 2013/01/09a EN





Yes, you can.®

General Information

For any order of spare parts the following information is required:

- Product type
- Product name & model
- Serial Number
- Part-Number
- Part-Description

Exploded view drawings and code numbers are not contractual and subject to modification, in relation to Invacare® continuous product improvement policy to meet end users needs.

The reproduction or use of this work, or of any substantial part of this work, in any form or by any electronic, mechanical or other means, including photocopying, scanning and recording, and in any information storage and retrieval system, including electronic transmission, is forbidden without the prior written permission of Invacare® International Sarl.



Yes, you can.®

TABLE OF CONTENTS

Drive unit complete viamobil V25

Drive unit complete viamobil V25	5
Drive unit complete viamobil V25 spare parts	6

Control unit viamobil V25 complete

Control unit viamobil V25 complete	7
Control unit viamobil V25 complete spare parts	8

Charger for viamobil V25 PS4515

Charger for viamobil V25 PS4515	9
---	---


Battery Pack for viamobil V25 complete

Battery Pack for viamobil V25 complete	10
Battery Pack for viamobil V25 complete spare parts	11


Others options & accessories

Others options & accessories	12
Additional push handles complete viamobil V25 pair comprising right and left spare parts	13
Anti-tippers V25 spare parts	14

Drive unit complete viamobil V25




Pos.	Part number	Description	Valid from → until	Qty
		Drive unit complete viamobil V25 spare parts		
	1521005	Drive unit complete viamobil V25		1




Pos.	Part number	Description	Valid from → until	Qty
	1456191	Cylinder head screw DIN 912 - M 6 x 20		1
	1456550	Set collar viamobil		1
	1459449	Threaded pin M6		1
	1520662	Hoist Motor V25		1
	1520666	Mounting profile V25 upper part		1
	1520667	Mounting profile V25 lower part		1
	1520698	Counterpart Profile V25 lower part		1
	1520725	Counterpart Profile V25 upper part		1
	1520765	Motor wire viamobil V25		1
	1520859	Cable harness with contact plate viamobil V25		1
	1521028	Bolt set viamobil V25 upper part		1
	1521029	Bolt set viamobil V25 lower part		1
	1521033	Cover without opening for push handles V25		1
	1521035	Drive Unit Cover V25		1
	1520575	Wheel incl. Tyre viamobil V25		1

Control unit viamobil V25 complete




Pos.	Part number	Description	Valid from → until	Qty
		Control unit viamobil V25 complete spare parts		
	1520689	Control unit viamobil V25 complete		1

Control unit viamobil V25 complete spare parts




Pos.	Part number	Description	Valid from → until	Qty
	1460702	Inlet control unit and handle for viamobil V15 chrome-plated		1
	1461578	Plug		1
	1463520	Screw Ejot PT K 30 x 14 WN		1
	1467795	Housing handle upper part V25/E25/E26		1
	1467803	Throttle Control Unit V25/E25/E26		1
	1467964	Housing Control Unit handle lower part E25/E26		1
	1520830	Cover Control Unit V25 bottom side		1
	1520833	Housing Control Unit V25 bottom part		1
	1520840	Switch ON/OFF Control Unit viamobil V25		1
	1520860	Cable for control unit V25 / E35 / E36		1
	1521034	Housing Control Unit V25 upper part		1
	1521042	Roll Support V25 complete		1

Charger for viamobil V25 PS4515




Pos.	Part number	Description	Valid from → until	Qty
	1520861	Charger for viamobil V25 PS4515		1
	1520182	Mains cable for charger 1.0A EU		1

Battery Pack for viamobil V25 complete




Pos.	Part number	Description	Valid from → until	Qty
	1564967	Battery Pack for viamobil V25 complete spare parts Battery Pack for viamobil V25 complete		1

Battery Pack for viamobil V25 complete spare parts




Pos.	Part number	Description	Valid from → until	Qty
	1521043	Handle of Battery Pack viamobil V25 complete		1

Others options & accessories




Pos.	Part number	Description	Valid from → until	Qty
		Additional push handles complete viamobil V25 pair comprising right and left spare parts		
	1521030	Additional push handles complete viamobil V25 pair comprising right and left		1
		Anti-tippers V25 spare parts		
	1521032	Anti-tippers V25		1

Additional push handles complete viamobil V25 pair comprising right and left spare parts



Pos.	Part number	Description	Valid from → until	Qty
	1459382	Cylinder head screw DIN 912 - M 6 x 80 - 8.8		1
	1460704	Ballcup viamobil chromium-plated		1
	1460719	Bolt		1
	1467795	Housing handle upper part V25/E25/E26		1
	1467802	Handle lower part V25/E25/E26		1
	1520852	Cover for handle viamobil V25		1
	1520979	Cover with opening for push handles V25		1
	1520982	Connecting piece left viamobil V25		1
	1520983	Wedge		1
	1520988	Connecting piece right viamobil V25		1
	1521033	Cover without opening for push handles V25		1
	1521047	Extra handle viamobil V25 left		1
	1521048	Extra handle viamobil V25 right		1
	1521049	Mounting for clamping		1
	1521125	Fixation set für V25 additional push handles		1
	1521212	Push handle for viamobil V25, complete		1

Anti-tippers V25 spare parts



Pos.	Part number	Description	Valid from → until	Qty
	1459382	Cylinder head screw DIN 912 - M 6 x 80 - 8.8		1
	1459944	Washer DIN 988 - 10 x 16 x 1		1
	1467845	Support angle right for anti-tippers		1
	1467846	Support angle left for anti-tippers		1
	1489365	Wheel for anti-tippers		1
	1519603	Securing disk DIN 6799-Gr.8 vz		1
	1521049	Mounting for clamping		1
	1521069	Anti-tipper viamobil V25 left		1
	1521070	Anti-tipper viamobil V25 right		1
	1521071	Fixation set for two anti-tippers V25		1
	1521081	Locking mechanism for anti-tipper V25 left		1
	1521088	Wedge		1
	1521089	Locking mechanism for anti-tipper V25 right		1
	1521113	Sleeve anti-tipper		1
	1521309	Toothed spline for V25 antitipper, assembled		1



Yes, you can.®



Belgium & Luxemburg:

Invacare nv, Autobaan 22, B-8210 Loppem
Tel: (32) (0)50 83 10 10, Fax: (32) (0)50 83 10 11
belgium@invacare.com / www.invacare.be

Denmark:

Invacare A/S, Sdr. Ringvej 37, DK-2605 Brøndby
Tel: (45) (0)36 90 00 00, Fax: (45) (0)36 90 00 01
denmark@invacare.com / www.invacare.dk

Deutschland:

Invacare Aquatec GmbH, Alemannenstraße 10, D-88316 Isny
Tel: (49) (0)75 62 7 00 0, Fax: (49) (0)75 62 7 00 66
info@invacare-aquatec.com / www.invacare-aquatec.de

European Distributor Organisation:

Invacare, Kleiststraße 49, D-32457 Porta Westfalica
Tel: (49) (0)57 31 754 540, Fax: (49) (0)57 31 754 541
edo@invacare.com / www.invacare.eu.com

España:

Invacare SA, c/Areny s/n, Polígon Industrial de Celrà, E-17460 Celrà (Girona)
Tel: (34) (0)972 49 32 00, Fax: (34) (0)972 49 32 20
contactsp@invacare.com / www.invacare.es

France:

Invacare Poirier SAS, Route de St Roch, F-37230 Fondettes
Tel: (33) (0)2 47 62 64 66, Fax: (33) (0)2 47 42 12 24
contactfr@invacare.com / www.invacare.fr

Ireland:

Invacare Ireland Ltd, Unit 5 Seatown Business Campus, Seatown Road, Swords, County Dublin - Ireland
Tel: (353) 1 810 7084, Fax: (353) 1 810 7085
ireland@invacare.com / www.invacare.ie

Italia:

Invacare Mecc San s.r.l., Via dei Pini 62, I-36016 Thiene (VI)
Tel: (39) 0445 38 00 59, Fax: (39) 0445 38 00 34
italia@invacare.com / www.invacare.it

Nederland:

Invacare BV, Celsiusstraat 46, NL-6716 BZ Ede
Tel: (31) (0)318 695 757, Fax: (31) (0)318 695 758
nederland@invacare.com, csede@invacare.com / www.invacare.nl

Norge:

Invacare AS, Grensesvingen 9, Postboks 6230, Etterstad, N-0603 Oslo
Tel: (47) (0)22 57 95 00, Fax: (47) (0)22 57 95 01
norway@invacare.com, island@invacare.com / www.invacare.no

Österreich:

Invacare Austria GmbH, Herzog Odilostrasse 101, A-5310 Mondsee
Tel: (43) 6232 5535 0, Fax: (43) 6232 5535 4
info-austria@invacare.com / www.invacare.at

Portugal:

Invacare Lda, Rua Estrada Velha, 949, P-4465-784 Leça do Balio
Tel: (351) (0)225 1059 46/47, Fax: (351) (0)225 1057 39
portugal@invacare.com / www.invacare.pt

Sverige & Suomi:

Invacare AB, Fagerstagatan 9, S-163 91 Spånga
Tel: (46) (0)8 761 70 90, Fax: (46) (0)8 761 81 08
sweden@invacare.com, finland@invacare.com / www.invacare.se

Switzerland:

Invacare AG, Benkenstrasse 260, CH-4108 Witterswil
Tel: (41) (0)61 487 70 80, Fax: (41) (0)61 487 70 81
switzerland@invacare.com / www.invacare.ch

United Kingdom:

Invacare Limited, Pencoed Technology Park, Pencoed, Bridgend CF35 5AQ
Switchboard Tel: (44) (0)1656 776200, Fax: (44) (0)1656 776201
Customer services Tel: (44) (0)1656 776222, Fax: (44) (0)1656 776220
uk@invacare.com / www.invacare.co.uk